



Aug 04, 2021

# Integration of a control brick

In 1 collection

Carolyn N Bayer<sup>1</sup>, Maja Rennig<sup>1</sup>, Anja Ehrmann<sup>1</sup>, Morten Norholm<sup>1</sup><sup>1</sup>Technical University of Denmark

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bvman42e](https://dx.doi.org/10.17504/protocols.io.bvman42e)

SEGA

Carolyn Bayer  
Technical University of Denmark

## ABSTRACT

SEGA, the Standardized Genome Engineering Architecture, is a comprehensive strain collection that enables genome engineering by combining only two reagents: a DNA fragment that can be ordered from a commercial vendor and a stock solution of bacterial cells followed by incubation on agar plates. Recombinant genomes are identified by visual inspection using green-white colony screening akin to classical blue-white screening for recombinant plasmids. The modular nature of SEGA allows precise multi-level control of transcriptional, translational, and post-translational regulation. The SEGA architecture simultaneously supports increased standardization of genetic designs and a broad application range by utilizing well-characterized parts optimized for robust performance in the context of the bacterial genome

This protocol describes the process of integrating a SEGA control brick. A control brick is integrated using *galk* counterselection.

DOI

[dx.doi.org/10.17504/protocols.io.bvman42e](https://dx.doi.org/10.17504/protocols.io.bvman42e)

## PROTOCOL CITATION

Carolyn N Bayer, Maja Rennig, Anja Ehrmann, Morten Norholm 2021. Integration of a control brick.

**protocols.io**<https://dx.doi.org/10.17504/protocols.io.bvman42e>

## COLLECTIONS ⓘ

 **SEGA protocol collection**

## KEYWORDS

SEGA, genome engineering, recombineering, E. coli, synthetic biology

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 08, 2021

## LAST MODIFIED

Aug 04, 2021

PROTOCOL INTEGER ID

50562

PARENT PROTOCOLS

Part of collection

[SEGA protocol collection](#)

BEFORE STARTING

Transform a SEGA strain, harbouring the upstream gadget *galK*, with pSIM19 (Spectinomycin resistance). From now on cultures have to be kept at **30 °C** to retain the plasmid (temperature-sensitive *ori*).

#### preculture and DNA fragment- Day 1

1d










1d

- 1 Prepare a PCR product of the control elements that need to be integrated and purify it from an agarose gel.
- 2 Setup a preculture of the strain harbouring pSIM19 in LB medium supplemented with Spectinomycin **0.05 mg/ml**. Incubate overnight at **250 rpm, 30°C**

#### Recombineering- Day 2







1d

- 3 Prepare:
  - 3.1 Cold sterile water
  - 3.2 Cold Glycerol **15 % volume**
  - 3.3 Pre-chilled centrifuge and tabletop centrifuge at **4 °C**
  - 3.4 M63 agar plates supplemented with **0.2 Mass / % volume** 2-deoxy-galactose, **0.2 % volume** glycerol and **5 Milimolar (mM)** L-rhamnose
- 4 Inoculate **50 mL** LB-Medium supplemented with Spectinomycin (**0.05 mg/ml**) with **500 µl** of the preculture from step 3
  - 4.1 Incubate at **250 rpm, 30°C** until cultures reached an OD<sub>600</sub> of 0.5
- 5 Induce expression by transferring the culture to a shaking water bath at **150 rpm, 42°C, 00:20:00**

- 6 Transfer culture to prechilled  **50 mL** falcon tubes and put on ice for  **00:15:00** 15m
- 7 Spin the culture down at  **4000 x g, 4°C, 00:05:00** and discard the supernatant 5m
- 8 Add  **1 mL** of ice cold water, resuspend and transfer to a 1.5 ml tube
- 9 30s  
Spin at  **11000 x g, 4°C, 00:00:30** in a tabletop centrifuge
- 10 Wash pellet twice with  **1 mL** ice cold water  **go to step #9**
- 11 Resuspend the pellet in  **600 µl** cold glycerol (  **15 % volume** )



11.1 Unused cells can be stored at  **-80 °C**

This is not possible for *E. coli* Nissle

- 12 Electroporate  **50 µl** of cells with  **200 ng** of purified PCR product from step 2 or  **2 µl** of a  **100 Micromolar (µM)** single-stranded oligonucleotide
- 13 Transfer cells into  **50 mL** LB medium in a 250 ml baffled conical shake flask and recover overnight at  **250 rpm, 30°C**

Alternatively, cells can be recovered in  **10 mL** medium for  **04:30:00**

Plating- Day 3 1d

- 14 Wash  **1 mL** of the recovered cells twice with 1X M9 salts. Centrifuge at  **11000 rpm, 20°C, 00:00:30** 30s

- 15 Make a dilution series and plate **100 µl** of the 1:100 - 1:1000 dilution on M63 agar supplemented with **0.2 Mass / % volume** 2-deoxy-galactose, **0.2 % volume** glycerol and **5 Milimolar (mM)** L-rhamnose.

In case the cells were recovered for **04:30:00** dilution 1:10 and 1:100 need to be plated

- 16 incubate the plates at **30 °C** for **48:00:00** to **72:00:00**

5d

Cells can be incubated at **37 °C** if it is the last step of recombination and pSIM19 is no longer needed in the cells

Screening- Day 4-5

1d

- 17 Screen for positive colonies by "green-white screening" on a blue-light table and perform colony PCR on the colorless colonies to identify the correct recombinants