



3 ▼

Jan 05, 2022

Cloning shRNA Oligos into pLKO.1 V.3

Addgene The Nonprofit Plasmid Repository¹

¹Addgene

1



dx.doi.org/10.17504/protocols.io.b3hxqj7n



Addgene The Nonprofit Plasmid Repository
Addgene

This is the protocol accompanying the "pLKO.1 – TRC Cloning Vector". For information about the pLKO.1-TRC cloning vector and tips on designing shRNA oligos for pLKO.1 see Addgene's website: <http://www.addgene.org/tools/protocols/plko/>

DOI

dx.doi.org/10.17504/protocols.io.b3hxqj7n

<http://www.addgene.org/tools/protocols/plko/>

Addgene The Nonprofit Plasmid Repository 2022. Cloning shRNA Oligos into pLKO.1. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b3hxqj7n>
Addgene The Nonprofit Plasmid Repository



protocol ,

Jan 05, 2022

Jan 05, 2022

56599

STEP MATERIALS

 [NEBuffer 3 - 5.0 ml](#) **New England**

Biolabs Catalog #B7003S

 [NEBuffer 1 - 5.0 ml](#) **New England**

Biolabs Catalog #B7001S Step 10

 [Agel - 300 units](#) **New England**

Biolabs Catalog #R0552S Step 11

 [EcoRI - 10,000 units](#) **New England**

Biolabs Catalog #R0101S In 2 steps

 [EcoRI - 10,000 units](#) **New England**

Biolabs Catalog #R0101S In 2 steps





 [T4 DNA Ligase Reaction Buffer - 6.0 ml](#) **New England**

Biolabs Catalog #B0202S Step 26

 [T4 DNA Ligase - 20,000 units](#) **New England**

Biolabs Catalog #M0202S Step 27

Annealing Oligos

- 1 Resuspend oligos in ddH₂O to a concentration of 20 µM.
- 2 Add 5ul Forward oligo
 **5 µL**
- 3 Add 5ul Reverse oligo
 **5 µL**
- 4 Add 5ul 10X NEB Buffer 2
 **5 µL**
- 5 Add 35 µL ddH₂O
 **35 µL**

- 6 Incubate for 4 minutes at 95°C in a PCR machine or in a beaker of boiling water.

🕒 00:04:00

- 7 Incubate the sample at 70°C for 10 minutes in a PCR machine.

🕒 00:10:00

- 8 Slowly cool to room temperature over the period of several hours.

🕒 03:00:00

This will take a few hours, but it is important for the cooling to occur slowly for the oligos to anneal.

If using a beaker of water, remove the beaker from the flame, and allow the water to cool to room temperature.

Digesting pLKO.1 TRC Cloning Vector

- 9 Mix: 6 µg pLKO.1 TRC-cloning vector (maxiprep or miniprep DNA)

🧴 6 µg

- 10 with 5 µL 10x NEB buffer 1

🧴 5 µL

🔗 NEBuffer 1 - 5.0 ml New England

Biolabs Catalog #B7001S

- 11 with 1 µL AgeI

🧴 1 µL





🔗 AgeI - 300 units New England

Biolabs Catalog #R0552S

- 12 Bring up to 50 µl with ddH₂O

- 13 Incubate at 37°C for 2 hours.

🕒 02:00:00

- 14 Purify with Qiaquick gel extraction kit, eluting in 30 μ L of ddH₂O.
- 15 Digest eluate with EcoRI by mixing: 30 μ L pLKO.1 TRC-cloning vector digested with AgeI
- 16 with 5 μ L 10x NEB buffer for EcoRI
 **5 μ L**
[EcoRI - 10,000 units New England](#)
Biolabs Catalog #R0101S
- 17 with 1 μ L EcoRI
 **1 μ L**
[EcoRI - 10,000 units New England](#)
Biolabs Catalog #R0101S
- 18 and 14 μ L ddH₂O
 **14 μ L**
- 19 Incubate at 37°C for 2 hours.
 **02:00:00**
- 20 Run digested DNA on 0.8% low melting point agarose gel until you can distinctly see 2 bands, one 7kb and one 1.9kb.
- When visualizing DNA fragments to be used for ligation, use only long-wavelength UV light. Short wavelength UV light will increase the chance of damaging the DNA.
- 21 Cut out the 7kb band and place in a sterile microcentrifuge tube.
- 22 Purify the DNA using a Qiaquick gel extraction kit. Elute in 30 μ L of ddH₂O.

23 Measure the DNA concentration.

Ligating and Transforming into Bacteria

24 Use your ligation method of choice. For a standard T4 ligation, mix: 2 μ L annealed oligo from "Annealing Oligos" section above.

 2 μ L

25 With 20 ng digested pLKO.1 TRC-cloning vector from the "Digesting pLKO.1 TRC Cloning Vector" section above.

 20 ng

If you were unable to measure the DNA concentration, use 1 μ L

26 With 2 μ L 10x NEB T4 DNA ligase buffer

 2 μ L

 [T4 DNA Ligase Reaction Buffer - 6.0 ml New England](#)

Biolabs Catalog #B0202S

27 With 1 μ L NEB T4 DNA ligase

 1 μ L

 [T4 DNA Ligase - 20,000 units New England](#)

Biolabs Catalog #M0202S

28 Bring up to 20ul with ddH2O

29 Incubate at 16°C for 4-20 hours.

 04:00:00

30 Transform 2 μ L of ligation mix into 25 μ L competent cells, following manufacturer's protocol.

Due to the long terminal repeats found in lentiviral plasmids, we recommend using a

strain that reduces the frequency of homologous recombination of unstable regions, such as Invitrogen Stbl3™ or NEB Stable cells. This will ensure that the repeats will be maintained and often results in a greater yield of DNA.

- 31 Plate on LB agar plates containing 100 µg/mL ampicillin or carbenicillin (an ampicillin analog).