



JAN 18, 2024

Single-Molecule Antibody Slides For Fluorescence Microscopy

Rebecca
Andrews¹

¹Univeristy of Cambridge



Rebecca Andrews
Univeristy of Cambridge

ABSTRACT

This protocol describes how to create single-molecule antibody slides for fluorescence microscopy.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1

Protocol Citation: Rebecca Andrews 2024. Single-Molecule Antibody Slides For Fluorescence Microscopy .

protocols.io

<https://dx.doi.org/10.17504/protocols.io.ewov1qbbygr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 18, 2024


Last Modified: Jan 18, 2024

PROTOCOL integer ID:
93745

Keywords: ASAPCRN



Coverslip Cleaning

1h 0m 10s

- 1 Coverslips (24 x 50 mm, #1, VWR,Catalogue Number 48404-453) were argon plasma cleaned (Ar plasma cleaner, PDC-002, Harrick Plasma) for  00:30:00 . 30m

Surface preparation


1h 0m 10s

- 2 A trimmed gasket was placed on top of the slide (CultureWell™ Reusable Gasket, 6mm diameter, Grace Bio-Labs, SKU: 103280).
- 3 Poly-L-Lysine  30 µL (0.01 % w/v PLL, P4707, Sigma-Aldrich) was placed in the wells for  00:30:00 . 30m
- 4 The PLL was removed, the wells washed three times with PBS (pH 7.4, 1x Gibco, Thermo Fisher Scientific, Catalogue number 10010023)

Depositing Antibodies

1h 0m 10s

- 5 A secondary antibody of choice was added at a final concentration of 0.0002 mg/ml in PBS.
- 6 The antibodies were left in the wells for 5-10 seconds for sufficient surface density before the wells were washed three times with PBS. 10s

7 PBS  30 µL was left in the wells for imaging