



Aug 29, 2024

DFP-induced mitophagy experiments

DOI

dx.doi.org/10.17504/protocols.io.e6nvw11m9lmk/v1

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.e6nvw11m9lmk/v1

Protocol Citation: Elias Adriaenssens 2024. DFP-induced mitophagy experiments. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvw11m9lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 101114

Keywords: ASAPCRN

Funders Acknowledgement:

**Aligning Science Across
Parkinson's (ASAP)**

Grant ID: ASAP-000350

**Marie Skłodowska-Curie
MSCA Postdoctoral
fellowship**

Grant ID: 101062916

Abstract

This protocol details the DFP-induced mitophagy experiments.

Materials

Deferiprone (DFP) (379409, Sigma Aldrich)

 3-Hydroxy-1,2-dimethyl-4(1H)-pyridone **Merck MilliporeSigma (Sigma-Aldrich) Catalog #379409**







DFP-induced mitophagy

2d

1 To induce mitophagy,

1d

- Treat the cells for  24:00:00 with  1 millimolar (mM) Deferiprone (DFP) (379409, Sigma Aldrich), an iron chelator that mimics hypoxic conditions through stabilization of the transcription factor HIF1 α and subsequent upregulation of NIX and BNIP3.

2 Freshly dissolve the compound in sterile dH₂O on the day itself and boil for  00:01:00 at  95 °C to dissolve the compound.

1m

