



Jan 28, 2022

DNA extraction from colonial tunicates

Marta Wawrzyniak¹, Simon Blanchoud¹¹University of Fribourg

1

dx.doi.org/10.17504/protocols.io.b33sqqne

Blanchoud lab, UNIFR

Marta Wawrzyniak
University of Fribourg

This protocol has been successfully used with *Botrylloides diegensis* and was adapted to our needs based on the HotPhenol DNA extraction protocol.

DOI

dx.doi.org/10.17504/protocols.io.b33sqqne

Marta Wawrzyniak, Simon Blanchoud 2022. DNA extraction from colonial tunicates. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.b33sqqne>

DNA extraction, colonial tunicates, ascidians

protocol ,

Jan 20, 2022

Jan 28, 2022

57170


Change gloves frequently, particularly as the protocol progresses from crude extracts to more purified materials. Use sterile tubes. Perform all steps on ice and use RNase-free and DNase-free water unless otherwise stated.

Heat bath setup at 70 °C (If working with fresh samples: glass beads 0.1mm and eppendorf thermal shaker)
phenol pH 8 (4 °C)
Ph-Ch-IA: phenol:chloroform:isoamyl alcohol (25:24:1, best prepared fresh)
Ch-IA: chloroform:isoamyl alcohol (24:1)
SDS-Lysis buffer : 10mL Lysis buffer, 4mL SDS 10%
Lysis buffer 50mL: 12.3g 3M Sodium acetate (pH 5.2), 7.3g 0.5M EDTA , Nuclease-free water
12.3g 3M sodium acetate in 50mL Nuclease-free water
80% Ethanol
ultra pure water

- 1 This protocol was developed to extract both RNA and DNA in parallel (See [RNA extraction from colonial tunicates](#) , steps 1-7) using the same samples. However it could be run directly on fresh samples (steps 1.1-1.4).

- 1.1 Clean the slide from which you will take the colony of your interest. See [Cleaning colonial ascidians](#).

- 1.2 Isolate a cleaned colony composed of approx. 20 zooids.

- 1.3 Transfer to a tube and spin at maximum speed for  **00:02:00** . 2m

- 1.4 Remove the excess water.

Cell Lysis 4m 30s

- 2 Add  **500 µL** of Ph-Ch-IA solution and  **350 µL** of SDS-Lysis buffer.

- 3 Heat the tube at  **70 °C** for  **00:02:00** . 2m

- 3.1 If working with fresh samples, add glass beads 0.1mm to the tube and shake on eppendorf thermal shaker at 900rpm.

4 Mix by vortexing at maximum speed for ⌚ 00:00:45 . 45s

5 Cool for ⌚ 00:01:00 on ice. 1m

6 Mix again by vortexing at maximum speed for ⌚ 00:00:45 . 45s

7 Heat the tubes at 🔥 70 °C for ⌚ 00:10:00 , mix regularly by inversion. 10m

7.1 If working with fresh samples, shake the tubes on eppendorf thermal shaker at 900rpm.

8 Mix again by vortexing at maximum speed for ⌚ 00:00:45 . 45s












9 Cool for ⌚ 00:01:00 on ice. 1m

10 Mix again by vortexing at maximum speed for ⌚ 00:00:45 . 45s

DNA extraction 3m

11 Centrifuge at 🔥 Room temperature at maximum speed for ⌚ 00:03:00 . 3m

12 Transfer 📄 400 µL of the upper aqueous phase to a new tube.

- 13 Add  **400 µL** of Ph-Cl-IA solution.
- 14 Shake the tube by inversion for  **00:00:30** . 30s
- 15 Centrifuge at maximum speed for  **00:03:00** . 3m
- 16 Transfer  **300 µL** of the upper aqueous phase to a new tube.
- 17 Add  **300 µL** of Ph-Cl-IA solution.
- 18 Shake the tube by inversion for  **00:00:30** . 30s
- 19 Centrifuge at maximum speed for  **00:03:00** . 3m
- 20 Transfer  **200 µL** of the upper aqueous phase to a new tube.
- 21 Add  **200 µL** of CHA solution.
- 22 Shake the tube by inversion for  **00:00:30** . 30s
- 23 Centrifuge at maximum speed for  **00:03:00** . 3m

24 Transfer aqueous phase to a new tube.

DNA precipitation

3h 30m

25 Add 2 volumes of [M]100 % volume Ethanol (typically ▢300-400 µL).

26 Add 0.1 volume of [M]3 Molarity (M) sodium acetate (typically ▢15-20 µL).

27 Mix by inversion.

28 Incubate at ⚡ -20 °C for ⌚ 03:00:00 .

3h

29 Centrifuge at maximum speed for ⌚ 00:20:00 at ⚡ 4 °C .

20m


30 Discard the supernatant.

31 Add ▢450 µL of cold [M]80 % volume ethanol.

32 Centrifuge at maximum speed for ⌚ 00:05:00 at ⚡ Room temperature .

5m


33 Discard the supernatant.

34 Add  **200 µL** of cold  **80 % volume** ethanol.

35 Centrifuge at maximum speed for  **00:05:00** at  **Room temperature** .

5m

36 Discard the supernatant.

37 Resuspend the pellet in ultra pure water (typically  **20-100 µL**).

38 Measure the DNA concentration using the NanoDrop.

39 Store at  **-20 °C** for short storage or at  **-80 °C** for long storage.