



Aug 05, 2024

Isolation and Transfection of *Nicotiana benthamiana* Mesophyll Protoplasts for Fluorescent Protein Visualization

This protocol is a draft, published without a DOI.

Namitha Nayak¹, Sandhya Mehrotra¹, Rajesh Mehrotra¹

¹BITS Pilani



Namitha Nayak

BITS Pilani

OPEN  ACCESS



Protocol Citation: Namitha Nayak, Sandhya Mehrotra, Rajesh Mehrotra 2024. Isolation and Transfection of *Nicotiana benthamiana* Mesophyll Protoplasts for Fluorescent Protein Visualization. **protocols.io** <https://protocols.io/view/isolation-and-transfection-of-nicotiana-benthamian-dihi4b4e>

Manuscript citation:

Nayak N, Mehrotra R, Mehrotra S (2024) The N-region sequence context impacts the chloroplast import efficiency of multi-TMD protein. *Plant molecular biology*, 114(4), 88. <https://doi.org/10.1007/s11103-024-01485-2>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 05, 2024

Last Modified: August 05, 2024

Protocol Integer ID: 104714

Keywords: Confocal microscopy, Fluorescent protein, *Nicotiana benthamiana*, Protoplast, Transfection

Funders Acknowledgement:

SERB

Grant ID: EMR/2016/002470



Abstract

Protoplast transfection is routinely used to study the sub-cellular localization of fluorescent-tagged proteins. Using *N. benthamiana* protoplasts is advantageous as it is a hardy plant that grows well, especially in tropical climates. Isolating protoplasts from the leaves of young 3-4-week-old *N. benthamiana* plants compensates for the benefits of using the ephemeral *Arabidopsis*. Moreover, the larger protoplast size of *N. benthamiana* offers better visualization of fluorescent proteins at lower magnifications. The protocol described here is an easy method for *N. benthamiana* protoplast isolation and transfection with simple and economical modifications to increase yield and transfection efficiency. The protocol is optimized for easy performance with minimal laboratory equipment.

Attachments



Protoplast isolation...

4MB

