

Apr 28, 2024 Version 1

MHV Tissue Titering Protocol V.1

This protocol is a draft, published without a DOI.

Siddharth Krishnamurthy¹

¹CU Anschutz



Siddharth

Siddharth Krishnamurthy

CU Anschutz, CU Anschutz

OPEN  ACCESS



Protocol Citation: Siddharth Krishnamurthy 2024. MHV Tissue Titering Protocol. [protocols.io](https://protocols.io/view/mhv-tissue-titering-protocol-dcuh2wt6) <https://protocols.io/view/mhv-tissue-titering-protocol-dcuh2wt6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 28, 2024

Last Modified: April 28, 2024

Protocol Integer ID: 98921

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Abstract

This protocol is for isolating total nucleic acid from soft tissues in mice, for subsequent analysis of Viral RNA levels

Protocol materials

⊗ RPMI 0% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R7755** In [2 steps](#)

⊗ 2.0 mm Zirconia beads **BioSpec Products Catalog #11079124ZX** Step 4

⊗ Sterile 24 Well Plate **VWR International Catalog #103348-844** Step 7

⊗ Deep 96 Well Plate **VWR International Catalog #10011-940** In [2 steps](#)

⊗ Ethyl alcohol, Pure **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023** Step 17

⊗ 1.1 mL Polypropylene Cluster Tubes, 12-Tube Strip Format, Nonsterile **VWR International Catalog #89005-574**

Step 2

⊗ 12-Well Cluster Tube Caps **VWR International Catalog #89005-728** Step 3

⊗ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632** Step 12

⊗ Buffer RLT Plus **Qiagen Catalog #1053393** Step 12

⊗ EconoSpin 96 Well DNA & RNA Binding Plate **Epoch Life Science Catalog #2020-001** Step 15

⊗ Buffer RW1 **Qiagen Catalog #1053394** Step 20



Day -1 (Or earlier): Design Well Layout for Tissue Collection tubes

- 1 Design the 96-well plate layout in which you will process (and eventually store) your samples

15m

Day -1: Prepare tissue collection tubes and plate

- 2 Label a sufficient number of 1.1 mL 12-well cluster tubes, and then place them in a new rack.

10m



1.1 mL Polypropylene Cluster Tubes, 12-Tube Strip Format, Nonsterile **VWR International Catalog #89005-574**

- 3 Retrieve enough 12-well cluster tube caps for your cluster tubes



12-Well Cluster Tube Caps **VWR International Catalog #89005-728**

- 4 Add 5-10 2.0 mm Zirconia beads to each tube that will have tissue in it

5m



2.0 mm Zirconia beads **BioSpec Products Catalog #11079124ZX**

- 5 Add  650 μ L RPMI 0% to each bead containing tube

5m



RPMI 0% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R7755**



- 6 Cover rack and tubes with resealable plate mat and leave in  4 °C until ready to use

5m



- 7 Add  6 mL (Approx) of RPMI 0% to a 24 well deep well plate



RPMI 0% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #R7755**



Sterile 24 Well Plate **VWR International Catalog #103348-844**



Day 0

1h 30m

- 8 Dissect the mouse and open the abdominal cavity without disturbing the adipose tissue




- 9 Collect mesenteric lymph nodes, Peyer's patches, and proximal colon

- 9.1 Orient mesenteric adipose tissue so mesenteric lymph nodes (mLN) are easily identifiable and place all mesenteric lymph nodes (with capsule) in the 1.1 cluster tube
- 9.2 Extract all the Peyer's patches and place them in the cluster tube
- 9.3 Cut  0.5 cm of proximal colon (the part of the colon that connects to the cecum)
- 10 Seal tube tightly with cap. Put whole weight on it, if necessary; it must be sealed by any means necessary 1m
- 11 Bead beat plate in SPEX MiniG Tissue Homogenizer for 5 min at 1500 rpm 5m
 1500 rpm, Room temperature , 00:05:00

Equipment

	NAME
SPEX	BRAND
SP 1600	SKU
MiniG 1600 Automated Tissue Homogenizer and Cell Lyser	SPECIFICATIONS



- 11.1 Seal the bead beater with the 96 well plate mat and a paper towel to examine if there is serious leakage
- 12 Thaw  2 Molarity (M) DTT to make a sufficient amount of RLT + 50 mM DTT

A	B	C
Number of Samples	Amount of RLT Plus (mL)	Amount of 2 M DTT (uL)
30	12.5	250

Spreadsheet to calculate how much RLT Plus and DTT will be needed

⌘ DTT **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D0632**

⌘ Buffer RLT Plus **Qiagen Catalog #1053393**

13 Aliquot 350 uL of RLT+DTT buffer to appropriate wells in a new 96 well deep well plate

⌘ Deep 96 Well Plate **VWR International Catalog #10011-940**

14 Once the plate is done beating: ⌘ 3000 x g, Room temperature, 00:02:00

2m

15 While the plate is spinning, set up the Multi-well Plate Manifold with a 96-well silica plate on top

⌘ EconoSpin 96 Well DNA & RNA Binding Plate **Epoch Life Science Catalog #2020-001**

Equipment

Multi-well plate manifold

NAME

Vacuum Manifold

TYPE

Pall

BRAND

5017

SKU

<https://www.cytivalifesciences.com/en/us/shop/lab-filtration/manifolds-and-accessories/microbiology-manifold/vacuum-manifold-and-accessories-p-36407>

LINK

16 Once the plate is done spinning, add 100 uL of tissue supernatant to the RLT-DTT plate



16.1 Use the pipet that you transfer the organ homogenate to pipet up and down to mix the RLT with the homogenate



16.2 Once the 100 uL has been transferred to the RLT-DTT plate, the RNA is stable, and you can transfer the 200 uL of the remaining homogenate to another deep 96-well plate for FFU titers



 Deep 96 Well Plate **VWR International Catalog #10011-940**

17 Add 350 uL of pure ethanol to the RLT-DTT well plate and do not mix here



 Ethyl alcohol, Pure **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7023**

18 Turn the vacuum on to prepare the silica plate

19 With a new pipet tip, mix the RLT+Ethanol solution and transfer to the silica plate



19.1 Transfer the solution with the same pipet you mix the ethanol with the pipet

19.2 The solution should take ~1 minute, usually less to suck through

20 Wash the silica plate with 350 uL RW1 Buffer



 Buffer RW1 **Qiagen Catalog #1053394**

21 Wash the silica plate with 800 uL RPE Buffer (10 mM Tris-Cl + 80% Ethanol)




22 Dry the plate  3000 x g, Room temperature, 00:02:00

2m



23 Place plate onto an elution skirted plate.

23.1 Elute by adding 75 uL of DEPC treated water to the wells and

 3000 x g, Room temperature, 00:02:00

2m

