

1 ▼

Oct 27, 2021

🌐 Artificial saliva V.1

Bjorn Bartholdy¹, a.g.henry¹

¹Leiden University



protocol .

BYOC



Bjorn Bartholdy

Creating an artificial saliva solution for oral biofilm growth.

This is a modified version of Sissons et al. 1991.

Bjorn Bartholdy, a.g.henry 2021. Artificial saliva. **protocols.io**
<https://protocols.io/view/artificial-saliva-bva9n2h6>



European Union's Horizon 2020 research and innovation program
 Grant ID: STG-677576 ("HARVEST")

protocol ,

May 26, 2021

Oct 27, 2021

50241

In steps of

[Amylase activity](#)

[Biofilm growth with starch treatment](#)

[Biofilm growth with starch treatment](#)

[Biofilm growth with starch treatment](#)

[Amylase activity](#)

[Biofilm growth with starch treatment](#)

This protocol explains how to mix  **1000 mL** of artificial saliva.




Chemicals

Mucin from porcine stomach (Type III)
Trypticase peptone
Proteose peptone
Yeast Extract
KCl
NaCl
CaCl₂
Na₂HPO₄
NaHCO₃
Hemin
Menadione
Urea
L-Arginine









Equipment

1000 mL beaker
2 x 1000 mL bottle
Magnetic stirrer (with heating element)
Autoclave
pH measure

Mix the solution under a fumehood. It can smell pretty bad.
Also have the incubator under a fumehood, if possible. I had no such luxury...













1 Add  **300 mL** distilled (or deionized) dH₂O to a  **1000 mL** beaker, with stirring and heat  **60 °C**.

2 Add:

-  **Mucin from porcine stomach (Type III) Sigma**
-  **2.5 g Aldrich Catalog #M1778**
-  **Trypticase™ Peptone Thermo**
-  **5 g Fisher Catalog #211921**
-  **Oxoid™ Proteose Peptone Thermo**
-  **10 g Fisher Catalog #LP0085B**
-  **5 g  Bacto Yeast Extract Becton-Dickinson**

Let the reagents completely dissolve before continuing to the next step

3 Add:

-  **2.5 g**  **KCl Contributed by users**
-  **0.35 g**  **NaCl Contributed by users**
-  **0.2 g**  **CaCl₂ Contributed by users**
 **Sodium phosphate dibasic Sigma**
-  **0.74 g** **Aldrich Catalog #7558-79-4**
-  **0.54 g**  **NaHCO₃ Contributed by users**
-  **2.5 mg**  **Hemin Contributed by users**

4 Add the remaining  **700 mL** distilled (or deionized) dH₂O

5 Adjust to  **pH 7** with  **NaOH Contributed by users** and stirring

6 Transfer to two 1000 ml bottles, so half of each bottle is filled.

15m







Autoclave at  **121 °C** ,  **1 Bar** for  **00:15:00** minutes



Do NOT screw bottle caps on tightly.

Loosely screw the caps on the bottles or cover the tops with foil

7 Once the solution has cooled, add:

-  **1 mg**  **Menadione Contributed by users**
-  **0.3 g**  **Urea Contributed by users**
 **L-Arginine Contributed by**
-  **0.17 g** **users Catalog #A5006**

8 Store in fridge at ca.  **4 °C**

Occasionally test the pH to ensure it stays around $\text{pH } 7$