



Oct 17, 2022

# Immunostaining of H&E Stained Paraffin Sections of Fly Heads

Mel Feany<sup>1,2</sup><sup>1</sup>Brigham and Women's Hospital; <sup>2</sup>Harvard Medical School

1 Works for me

 Share[dx.doi.org/10.17504/protocols.io.81wgby61nvpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgby61nvpk/v1)

Daniel's workspace

 Daniel El Kodsi

## ABSTRACT

This protocol describes how to perform immunostaining on H&E stained paraffin sections of fly heads.

## DOI

[dx.doi.org/10.17504/protocols.io.81wgby61nvpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgby61nvpk/v1)

## PROTOCOL CITATION

Mel Feany 2022. Immunostaining of H&E Stained Paraffin Sections of Fly Heads.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.81wgby61nvpk/v1>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 17, 2022

## LAST MODIFIED

Oct 17, 2022

## PROTOCOL INTEGER ID

71457

- 1 Embed the fly heads in paraffin and stain with H&E for neuronal counting, per the following protocol:

- 2 Microwave slides in [M]10 millimolar (mM) sodium citrate for ⌚00:15:00 . 35m  
Cool ⌚00:20:00 .
- Stock: [M]100 millimolar (mM) sodium citrate, pH 6.0  
Use at least 📏1 L of citrate solution in large glass box to avoid drying.
- 3 Block in PBST (PBS with 0.3% Triton) with 2% dry milk for ⌚00:30:00 to 1h 30m  
Stock: 10X PBS
- 4 Incubate with primary antibody ⌚Overnight at 🌡Room temperature .
- 5 Wash 3 x in PBST.
- 6 Incubate with appropriate biotinylated secondary (for immunohistochemistry) or fluorescent<sup>2h</sup> secondary (for immunofluorescence) antibody at 1:200 in PBST + milk for ⌚01:00:00 at 🌡Room temperature .  
For immunohistochemistry incubate in ABC reagent (Vector) for ⌚01:00:00 at 🌡Room temperature .
- 7 Rinse 3 x PBST.
- 8 Mount slides with antifading medium for immunofluorescence or dehydrate through ethanol series and xylenes and mount in Permount for immunohistochemistry.