



1 ▼

Jan 05, 2022

# SWAP mouse genotyping V.1

Peter W W Chomczynski<sup>1</sup>, Kianna M Vires<sup>2</sup>, Michal Rymascewski<sup>1</sup>,  
Judith A. Heiny<sup>2</sup>

<sup>1</sup>Molecular Research Center, Inc.; <sup>2</sup>University of Cincinnati

1



protocol .

Molecular Research Center, inc.



Peter W Chomczynski

The highly-conserved, cardiotonic steroid (CTS) binding site (also termed ouabain binding site) on the primary  $\alpha$  subunit of Na,K-ATPase (NKA) plays a receptor signaling role in a range of vital cell processes and is a therapeutic target for human disease. Mouse lines with altered affinity for CTS on the  $\alpha 1$  or  $\alpha 2$  subunit isoform of NKA, developed by the late Dr. Jerry Lingrel and collaborators, are a valuable tool for studying its physiological roles and drug actions. In one model, the normally ouabain resistant  $\alpha 1$  isoform was rendered sensitive to ouabain binding. In a second model, the normally sensitive  $\alpha 2$  isoform was rendered resistant to ouabain binding. Additional useful models are obtained by mating these mice. The Lingrel mice will soon become publicly available in a biorepository. To further advance their use, we developed a rapid real-time, direct PCR method that amplifies digested tail clip samples without the need to isolate DNA, and detects variant alleles using specific primers and fluorescent probes. PCR is performed in fast mode with up to 23 samples processed in 40 min. The method was validated using mice of known genotype and results were compared with a previous two-step method using PCR amplification followed by gel electrophoresis. We also performed Sanger sequencing to clarify inconsistencies in published sequences, update numbering to current reference sequences, and confirm the continued presence of the mutations in the colony. It is expected that the wider availability of these models and a more efficient genotyping protocol will advance studies of the NKA and its CTS receptor.

Peter W W Chomczynski, Kianna M Vires, Michal Rymascewski, Judith A. Heiny  
2022. SWAP mouse genotyping. **protocols.io**  
<https://protocols.io/view/swap-mouse-genotyping-b2kaqcse>



SWAP mouse, ATP1A1, ATP1A2, Ouabain binding site, Mouse models, NaK ATPase

protocol ,

Diagrams created by BioRender.com. Screenshots from ABI StepOne software.

Dec 03, 2021

Jan 05, 2022

55650

**Common reagents:**

1M Tris pH 8.0

0.5 M EDTA pH 8.0

5 M NaCl

SDS 10% w/v

NP40 detergent (10%)

Tween-20

proteinase K (20 mg/mL)

ddH<sub>2</sub>O

Bio-Rad iTaq Universal Probes Supermix (or compatible equivalent)

**Protocol-specific primers and probes:**

ATP1A1 FWD primer (100 mM): CAG CTC TTT GGA GGC TTT

ATP1A1 REV primer (100 mM): GCT ACC GTA ACT ACA CAA CTC

ATP1A1 WT probe (100 mM): /56-FAM/CA+T +CC+G +A+AG T+GC /3IABkFQ/

ATP1A1 mutant probe (100 mM): /56-FAM/TGG AAT +TC+A +G+AG T+GC /3IABkFQ/

ATP1A2 FWD primer (100 mM): TCC TCT GCT TCT TAG CCT ATG G

ATP1A2 REV primer (100 mM): CAG GGC TAT AAG CAG GTC CA

ATP1A2 WT probe (100 mM): /56-FAM/CAC ATT ATC /ZEN/GTT GGA TGG TTC GTC CTC  
C/3IABkFQ/

ATP1A2 mutant probe (100 mM): /56-FAM/CTC ACA TCA /ZEN/TCG TTC GAA GGC TCG  
TC/3IABkFQ/

Pre-experiment preparation

- 1 In the interest of time and consistency, it is recommended that certain stock solutions and buffers be prepared ahead.

- 1.1 Prepare a **10X stock solution of Tail Lysis Buffer** and store in freezer.<sup>15m</sup>

Reagent	Vol. to make 10 mL
1M Tris pH 8.0	1 mL
0.5 M EDTA pH 8.0	2 mL
5 M NaCL	2 mL
SDS 10% w/v	5 mL

#### Tail Lysis Buffer, 10X stock solution

## 1.2

20m

Prepare **assay mixes**.

Primers and probes should be pre-mixed and stored frozen for convenience. Be mindful of the reagent stock concentrations to ensure a successful outcome. Mixes can be scaled as necessary.

Reagent	For 200 rxn
ATP1A1 FWD primer (100 mM)	20 µL
ATP1A1 REV primer (100 mM)	20 µL
ATP1A1 WT probe (100 mM)	10 µL
ddH <sub>2</sub> O	150 µL

#### α1<sup>R</sup> assay mix

Reagent	For 200 rxn
ATP1A1 FWD primer (100 mM)	20 µL
ATP1A1 REV primer (100 mM)	20 µL
ATP1A1 mutant probe (100 mM)	10 µL
ddH <sub>2</sub> O	150 µL

#### α1<sup>S</sup> assay mix

Reagent	For 200 rxn
ATP1A2 FWD primer (100 mM)	20 $\mu$ L
ATP1A2 REV primer (100 mM)	20 $\mu$ L
ATP1A2 WT probe (100 mM)	10 $\mu$ L
ddH <sub>2</sub> O	150 $\mu$ L

**$\alpha$ 2<sup>S</sup> assay mix**

Reagent	For 200 rxn
ATP1A2 FWD primer (100 mM)	20 $\mu$ L
ATP1A2 REV primer (100 mM)	20 $\mu$ L
ATP1A2 mutant probe (100 mM)	10 $\mu$ L
ddH <sub>2</sub> O	150 $\mu$ L

**$\alpha$ 2<sup>R</sup> assay mix**


- 1.3 Prepare **Tail Digestion Buffer** fresh for each experiment, per the following<sup>15m</sup> table. Add proteinase K last from a frozen and thawed aliquot.

Reagent	Vol. to make 10 mL
10X Tail Lysis Buffer *	1 mL
NP40	45 $\mu$ L
Tween-20	45 $\mu$ L
proteinase K (20 mg/mL)	200 $\mu$ L
ddH <sub>2</sub> O	to 10 mL

**Tail Digestion Buffer**

Tail clip digestion 12h 45m

- 2 Clip 2-3 mm from the tail of each mouse to be genotyped. Place in a clean, labelled 1.5 mL<sup>30m</sup> microcentrifuge tube.




Steps 2-4 visual overview

- 3 Add **150 µL** of **Tail Digestion Buffer** containing proteinase K to each tube. 5m
- 4 Place the tail samples in a heat block or water bath. 3m
- 4.1 Incubate **Overnight** at **55 °C**. 15m


4.2 Raise the temperature to **100 °C** for **00:15:00** . 15m

4.3 Remove the tubes and allow to cool to **Room temperature** . 10m

5 Briefly vortex and centrifuge the tubes to pellet insoluble material. DNA will remain in the supernatant. 5m



6 Create a 1:100 dilution of each sample in a new, labelled 1.5 mL microcentrifuge tube. Add <sup>10m</sup>  
**198 µL** of TE buffer (or ddH<sub>2</sub>O) to each tube, then add **2.0 µL** of supernatant from the digested sample.



- 7 Store samples and dilutions at  $4^{\circ}\text{C}$  until ready to genotype. If storing samples longer than 14 days, keep in  $-20^{\circ}\text{C}$  freezer.

#### PCR machine setup

- 8 Create a new experiment on your PCR machine's software. Templates for ABI StepOne and StepOnePlus machines are attached.

[SWAP mouse genotype StepOne template.edt](#)

[SWAP mouse genotype StepOnePlus template.edt](#)

- 9 Enter your list of probes and samples into the machine software. Assign 8 wells per sample; 2 wells for  $a1^R$ , 2 wells for  $a1^S$ , 2 wells for  $a2^S$ , and 2 wells for  $a2^R$ .

A sample plate layout for genotyping 3 mice is shown (ABI StepOne software).

Define Targets and Samples **Assign Targets and Samples**

**Instructions:**  
 To set up standards: Click "Define and Set Up Standards."  
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as the task for each target assignment, then assign a sample.  
 To set up negative controls: Select wells, assign target(s), then select "N" (Negative Control) as the task for each target assignment.

**Assign target(s) to the selected wells.**

Assign	Target	Task	Quantity
<input checked="" type="checkbox"/>	a1R		
<input type="checkbox"/>	a1S		
<input type="checkbox"/>	a2S		
<input type="checkbox"/>	a2R		

☒ Mixed ☐ Unknown ☐ Standard ☐ Negative Control

**Assign sample(s) to the selected wells.**

Assign	Sample
<input checked="" type="checkbox"/>	Mouse 1
<input type="checkbox"/>	Mouse 2
<input type="checkbox"/>	Mouse 3
<input type="checkbox"/>	Neg. Ctrl.

**Assign sample(s) of selected well(s) to biological group.**

Assign:

**Select the dye to use as the passive reference.**

ROX: ☒

**View Plate Layout** **View Well Table**

Select Wells With:

Show in Wells:

	1	2	3	4	5	6	7	8
A	Mouse 1	Mouse 1	Mouse 1	Mouse 1	Mouse 1	Mouse 1	Mouse 1	Mouse 1
B	Mouse 2	Mouse 2	Mouse 2	Mouse 2	Mouse 2	Mouse 2	Mouse 2	Mouse 2
C	Mouse 3	Mouse 3	Mouse 3	Mouse 3	Mouse 3	Mouse 3	Mouse 3	Mouse 3
D	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.	Neg. Ctrl.
E								
F								

Wells: 24 Unknown 0 Standard 8 Negative Control

16 Empty

10 Set the thermocycling program as per the following table, then start the cycling.

Stage	Duration (m:ss)	Temperature	Cycles
Initial denaturation	3:00	95 C	1
Denaturation	0:03	95 C	35 Cycles
Annealing/Extension	0:30	95 C	

PCR reaction setup 2m





11 Prepare 4 reaction mixes as follows, 1 for each of the assay mixes. Create enough mix for the total number of reactions, plus 2 extra to account for potential pipetting loss.

15m

Reagent	1 rxn	50 rxn	100 rxn
iTaq Probe mix (2X)	10.0 µL	500 µL	1000 µL
Assay mix	1.0 µL	50 µL	100 µL
ddH2O	4.0 µL	200 µL	400 µL


PCR reaction mix



- 12 Pipet  **15  $\mu$ L** reaction mix into each well of a 96-well plate. 5m
- 13 Pipet  **5  $\mu$ L** of each diluted sample into its assigned wells. Pipet  **5  $\mu$ L** ddH<sub>2</sub>O into each negative control well. 7m
- 14 Cover the plate with optical sealing film. Ensure that the film is fully adhered to the plate. 1m
- 15 Briefly vortex the plate. 1m
- 16 Centrifuge the plate at  **1700 rpm, 00:02:00** . 2m
- 17 Insert the plate into the instrument and initiate the run. 1h 45m

### Interpretation

- 18 Analysis of the PCR results begins by checking the outcome of the negative control wells. If all negative controls show no amplification, you may proceed with the analysis. If one or more negative control wells show amplification, the results are invalid and the experiment must be repeated. An invalid outcome is typically due to contamination of the reaction mix or of a reagent.
- 19 Genotyping is based on the presence or absence of signal from the 2 probes of each gene. The example below shows representative amplification plots for the possible genotypes of  $\alpha 1$  (**A**) and  $\alpha 2$  (**B**). The graphs show results from a sample's  $\alpha 1^R$  (black) and  $\alpha 1^S$  (red) wells plotted together (**A**), and a sample's  $\alpha 2^S$  (black) and  $\alpha 2^R$  (blue) wells plotted together (**B**).



Representative amplification curves

Samples considered positive for a specific allele show  $C_T$  values of 26-31 for the corresponding probe, while negative samples do not reach threshold in 35 cycles. Heterozygous samples show amplification of both probes within 2  $C_T$  of each-other.