

May 09, 2024

## Native Nanobleach of native nanodiscs

DOI

**dx.doi.org/10.17504/protocols.io.ewov19nbklr2/v1**

Caroline Brown<sup>1,2,3</sup>, Gerard Walker<sup>1,2,4</sup>, Snehasish Ghosh<sup>1,2,3</sup>, Kallol Gupta<sup>1,2,3</sup>, Moitrayee Bhattacharyya<sup>4</sup>

<sup>1</sup>Nanobiology Institute, Yale University, West Haven, CT, USA;

<sup>2</sup>Department of Cell Biology, Yale University School of Medicine, New Haven, CT, USA;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815;

<sup>4</sup>Department of Pharmacology, Yale University, New Haven, CT, USA



Caroline Brown

Yale University

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.ewov19nbklr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov19nbklr2/v1)

**Protocol Citation:** Caroline Brown, Gerard Walker, Snehasish Ghosh, Kallol Gupta, Moitrayee Bhattacharyya 2024. Native Nanobleach of native nanodiscs. protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov19nbklr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 09, 2024

**Last Modified:** May 09, 2024

**Protocol Integer ID:** 99521

## Abstract

This is a protocol to conduct Native Nanobleach of native nanodiscs.



- 1 Note, this method is adapted from Walker, Brown, et al. (DOI: [10.1038/s41565-023-01547-4](https://doi.org/10.1038/s41565-023-01547-4)).
- 2 Prepare flow chamber by treating with poly-L-lysine (PLL) PEG and PEG-biotin.
- 3 Add streptavidin to flow chamber.
- 4 Add biotinylated GFP-nanobody to flow chamber.
- 5 Add native nanodiscs to the flow chamber for immobilization.
- 6 Conduct step-photobleaching analysis using TIRF microscopy.