




Aug 26, 2021

GUV preparation and assay

Forked from [GUV assay](#)Chunmei Chang¹¹Team Hurley, University of California, Berkeley

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bxrgpm3w Chunmei Chang

ABSTRACT

LC3 lipidation on GUVs

DOI

dx.doi.org/10.17504/protocols.io.bxrgpm3w

PROTOCOL CITATION

Chunmei Chang 2021. GUV preparation and assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bxrgpm3w>

FORK NOTE

FORK FROM

[Forked from GUV assay, Chunmei Chang](#)

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2021

LAST MODIFIED

Aug 26, 2021

PROTOCOL INTEGER ID

52744

1 GUV Preparation

1.1 Clean the coverslips of 25 mm diameter.

1.2 Coat cleaned coverslips with 60 μ L 5% (w/w) polyvinyl alcohol (PVA) with a molecular weight of 145,000 (Millipore).

Place the coated coverslip in a heating incubator at 60 °C to dry the PVA film for 30 min.

1.3

1.4 Spread a lipid mixture with a molar composition of 64.8% DOPC, 20% DOPE, 10% POPI, 5% DOPS and 0.2% Atto647N DOPE at 1 mg/ml uniformly onto the PVA film.

1.5 Put the lipid-coated coverslip under vacuum overnight to evaporate the solvent.

1.6 Use 400 μ L 400 mOsm sucrose solution for swelling for 1 h at room temperature

1.7 Harvest the GUVs and use them with 12 h.

2 GUV Assay

2.1 Set up the reaction in an eight-well observation chamber (Lab Tek) at room temperature.

2.2 Coat the chamber with 5 mg/ml β casein for 30 min.

2.3 Wash the coated chamber three times with reaction buffer (20 mM HEPES at pH 8.0, 190 mM NaCl and 1 mM TCEP).

2.4 Make a 120 μ L reaction mixtures with the proteins and 50 μ M ATP. The final concentration the each protein is 5 μ M GST-4xUb, 500 nM cargo receptors, 25 nM ULK1 complex, 25 nM PI3KC3-C1 complex, 100 nM WIPI2d, 50 nM ATG12-ATG5-ATG16 complex, 100 nM ATG7, 100 nM ATG3, and 500 nM mCherry-LC3B.

2.5 Add 10 μ L GUVs to initiate the reaction.

2.6 Pick random views for imaging within 5 min.

2.7 Acquire time-lapse images in multitracking mode on a Nikon A1 confocal microscope.

