



Jul 29, 2024

# NEB #T4010 Monarch Mag Viral DNA/ RNA Extraction Kit Protocol for KingFisher Flex Automated Isolation of Viral DNA/RNA from Wastewater Samples Following Ceres Nanotrap Enrichment

DOI

[dx.doi.org/10.17504/protocols.io.n92ld87o9v5b/v1](https://doi.org/10.17504/protocols.io.n92ld87o9v5b/v1)

Juliet Bonnevie<sup>1</sup>

<sup>1</sup>New England Biolabs

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)



Juliet Bonnevie

New England Biolabs

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.n92ld87o9v5b/v1](https://doi.org/10.17504/protocols.io.n92ld87o9v5b/v1)

**Protocol Citation:** Juliet Bonnevie 2024. NEB #T4010 Monarch Mag Viral DNA/ RNA Extraction Kit Protocol for KingFisher Flex Automated Isolation of Viral DNA/RNA from Wastewater Samples Following Ceres Nanotrap Enrichment. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ld87o9v5b/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 18, 2024

**Last Modified:** July 29, 2024

**Protocol Integer ID:** 103659



## Abstract

The Monarch Mag Viral DNA/RNA Extraction Kit provides a rapid and reliable magnetic bead-based process for extracting viral nucleic acids from saliva and respiratory swab samples. The kit combines the efficiency of silica-based nucleic acid purification with the ease of use of magnetic beads. Manual and automated workflows allow samples to be processed in microfuge tubes or 96-well plates. Kit sizes align to 96-well formats (100 preps, 600 preps, and 1800 preps), and the protocol is compatible with high throughput automation on a variety of platforms, including the KingFisher Flex magnetic particle processor and Agilent Bravo and MGLSP liquid handler platforms.

## Materials

### Reagents and Materials Supplied by User:

- 100% ethanol
- 100 % isopropanol
- Nuclease-free water
- RNase-free tips, tubes, and plastics.
- Adhesive seals for 96-well plates (KingFisher Flex automation protocol)

### Required Equipment For the Automation Protocol

- Vortex mixer
- KingFisher Flex or liquid handler (configured to align with the protocol)
- Automation platform-compatible plastics (e.g., 96-well deep well plates, 96-well microplates)  
Thermal mixer containing block for 96-well plates or plate shaker (may be required, depending on the automation platform)

### Required plastics

- KingFisher 96-deep well plates, v-bottom, (2.0 ml), Catalog # 95040450
- KingFisher 96 microplate (200 µl), Catalog # 9700254
- KingFisher 96 deep-well tip comb and plate, Catalog # 97002820



## Before start

### Important Notes Before You Begin

- Review Reagent Preparation section in the manual on NEB.com.
- Store Proteinase K at  $-20^{\circ}\text{C}$  upon receipt.
- Prepare Monarch Carrier RNA based on kit size used: Add 125  $\mu\text{l}$  (NEB #T4010S) or 750  $\mu\text{l}$  (NEB #T4010L/X) nuclease-free water, invert or pipette to mix, and transfer to an RNase-free microfuge tube. Keep on ice. Prepare single-use aliquots and store at  $-20^{\circ}\text{C}$ . Avoid multiple freeze-thaw cycles.
- Prepare 80% ethanol: 80% ethanol should be prepared fresh using 100% ethanol (user supplied) and nuclease-free water (user supplied). Prepare 1 ml of 80% ethanol per reaction and add overage.
- Perform all steps at room temperature unless directed otherwise.

### Required plastics

- KingFisher 96-deep well plates, v-bottom, (2.0 ml), Catalog # 95040450
- KingFisher 96 microplate (200  $\mu\text{l}$ ), Catalog # 97002540
- KingFisher 96 deep-well tip comb and plate, Catalog # 97002820

### Starting Material Notes

This protocol has been optimized for use with wastewater samples that have been pre-processed with Ceres Nanotrap Particles.



## Part I. Prepare the KingFisher Flex Instrument

- 1 Ensure the instrument is equipped with the KingFisher Flex 96 Deep Well head and the KingFisher Flex 96 heating block.  
IMPORTANT: The heat block must be compatible with the KingFisher 96 microplate (200  $\mu$ l).
- 2 Ensure the MagMAX Pathogen RNA/DNA (High Volume) program is loaded onto the instrument's connected computer and that the program has been modified to perform three 500  $\mu$ l wash steps, a 2-minute bead drying step, and a 33–100  $\mu$ l elution.
- 3 Enter sample, wash, and elution volumes into the program.
- 4 Select plate sizes for the run: KingFisher 96-deep well plates (2.0 ml) for sample and wash plates; KingFisher 96 microplate (200  $\mu$ l) for elution.

## Part II. Buffer Preparation

- 5 Prepare fresh Viral DNA/RNA Wash Buffer in a user-supplied tube or bottle (free of nucleases) according to the table.  
Add components in order, as listed. Prepare up to 15% excess to ensure a sufficient volume is available for each reaction.
- 6 Prepare Lysis Buffer Bead Mix immediately before use, according to the table.
  - a. Vortex magnetic beads to form a homogeneous solution before use.
  - b. Add components in order, as listed.
  - c. For a master mix, prepare up to 15% excess to ensure a sufficient volume of buffer/bead mix is available for each reaction.
  - d. Store Lysis Buffer Bead Mix at room temperature. Periodically invert or vortex to keep beads in suspension.

Viral DNA/RNA Wash Buffer:

| A                              | B                   |
|--------------------------------|---------------------|
|                                | Volume per Reaction |
| a. Combine the following:      |                     |
| Monarch Buffer BX              | 167 $\mu$ l         |
| Nuclease-free Water            | 83 $\mu$ l          |
| b. Vortex to mix and then add: |                     |
| Isopropanol                    | 250 $\mu$ l         |
| c. Vortex to mix               |                     |
| Total Volume                   | 500 $\mu$ l         |

### Lysis Buffer Bead Mix

| A                                | B                   |
|----------------------------------|---------------------|
|                                  | Volume per Reaction |
| a. Combine the following         |                     |
| Monarch StabiLyse DNA/RNA Buffer | 200 µl              |
| Monarch Carrier RNA              | 1 µl                |
| b. Vortex to mix and then add:   |                     |
| Isopropanol                      | 200 µl              |
| c. Vortex to mix and then add:   |                     |
| Monarch Mag Beads M1             | 20 µl               |
| d. Gently vortex to mix          |                     |
| Total Volume                     | 421 µl              |

## Part III. Prepare Wash and Elution Plates

- 7 Aliquot 500 µl Viral DNA/RNA Wash buffer to wells in a 96-well deep well plate.
- 8 Aliquot 500 µl 80% ethanol to wells in each of two 96-well deep well plates.
- 9 Aliquot 33–100 µl nuclease-free Water to wells in a 96-well microplate.
- 10 Seal plates with an adhesive film until ready to use.

| A                    | B                                             | C                                           | D                             | E                             | F                                        | G                              |
|----------------------|-----------------------------------------------|---------------------------------------------|-------------------------------|-------------------------------|------------------------------------------|--------------------------------|
| Plate Position       | 1                                             | 2                                           | 3                             | 4                             | 5                                        | 6                              |
| Plate type           | 96 deep well                                  | 96 deep well                                | 96 deep well                  | 96 deep well                  | 96-well microplate                       | Tip comb in 96-well microplate |
| Plate Identification | Sample plate                                  | Wash plate 1                                | Wash plate 2                  | Wash plate 3                  | Elution plate                            | N/A                            |
| Plate Contents       | Sample/Lysis Buffer Bead Mix (approx. 621 µl) | Viral DNA/RNA Wash Buffer (500 µl per well) | 80% ethanol (500 µl per well) | 80% ethanol (500 µl per well) | Nuclease-free water (33-100 µl per well) | N/A                            |



## Part IV. Elute target microbes from the Nanotrap pellet

- 11 *These steps will begin at the end of enrichment, specifically where Ceres protocol adds a Lysis Buffer. For example, for Protocol APP-042, you would follow the Nanotrap protocol till Step 11 and perform the steps below instead of Step 12.*

Add 100 µl Monarch StabiLyse DNA/RNA Buffer to the Nanotrap Particle pellet, pipette to resuspend the pellet.

- 12 Incubate for 10 mins at room temperature.

- 13 Use a magnetic rack that is compatible with sample tubes to separate Nanotrap Particles from the sample.

- 14 Transfer the supernatant to a deepwell plate compatible with processing on KingFisher.

15 **Sample Lysis (Sample Plate)**

Add 100 µl nuclease-free water to the Nanotrap enriched sample from Step 14 above.

- 16 111 Add 5 ul Proteinase K.

**Note**

Instead of individually performing Step 15 and 16, a master mix of nuclease water and Proteinase K can be made if desired.

- 17 Pipette mix or use a thermomixer for 30 seconds after sealing the plate.

- 18 Incubate for 15 min at room temperature.

- 19 Gently vortex the prepared Lysis Buffer Bead Mix and add 421 µl to each sample well.

- 20 Seal plate with adhesive film until ready to load onto the KingFisher Flex instrument.



## Part V. Viral Nucleic Acid Purification (Bind, Wash, Elute)

- 21 Carefully remove adhesive film from sample, wash, and elution plates.
- 22 Load sample, wash, elution plates, tip comb and plate, into the appropriate positions on the KingFisher Flex worktable.
- 23 Run the modified MagMAX program.
- 24 Upon completion of the run, seal the elution plate with adhesive film and place on ice for immediate use or freeze for storage.