



NOV 07, 2022

SHARE

WORKS FOR ME 1

## Cyanobacteria Total Lipid Extraction from Cell Pellets

 Forked from [Cyanobacteria Total Lipid Extraction from Polycarbonate Filters](#)

DOI

[dx.doi.org/10.17504/protocols.io.4r3l2oe3jv1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l2oe3jv1y/v1)[Robbie Martin](#)<sup>1</sup>, [Steven W Wilhelm](#)<sup>1</sup>, Katarina A. Jones<sup>1</sup>, [Shawn Campagna](#)<sup>1</sup><sup>1</sup>University of Tennessee, Knoxville[The Aquatic Microbial Ecology Research Group - AMERG \(The Buchan, Zinser and Wilhelm labs\)](#)[CyanoHABS](#)[1 more workspa](#)

Robbie Martin

COMMENTS 0

## ABSTRACT

This protocol is designed/used for extraction of total cellular lipids from cyanobacteria samples (either lab cultures or field samples) collected via centrifugation for use in lipid analysis and quantification *via* mass spectrometry.

Please contact Dr. Steven Wilhelm ([wilhelm@utk.edu](mailto:wilhelm@utk.edu)) or Robbie M. Martin ([rmarti49@vols.utk.edu](mailto:rmarti49@vols.utk.edu)) for additional information regarding this protocol.

Modified from Guan, X. L., Riezman, I., Wenk, M. R., & Riezman, H. (2010). Yeast lipid analysis and quantification by mass spectrometry. *Methods in Enzymology*, 470, 369-391.

DOI

[dx.doi.org/10.17504/protocols.io.4r3l2oe3jv1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l2oe3jv1y/v1)

## PROTOCOL CITATION

Robbie Martin, Steven W Wilhelm, Katarina A. Jones, Shawn Campagna 2022. Cyanobacteria Total Lipid Extraction from Cell Pellets. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.4r3l2oe3jv1y/v1>



MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Guan, X. L., Riezman, I., Wenk, M. R., & Riezman, H. (2010). Yeast lipid analysis and quantification by mass spectrometry. *Methods in Enzymology*, 470, 369-391.

## FORK NOTE

FORK FROM

[Forked from Cyanobacteria Total Lipid Extraction from Polycarbonate Filters, Robbie Martin](#)

KEYWORDS

cyanobacteria, lipids, extraction

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 12, 2022

LAST MODIFIED

Nov 07, 2022

PROTOCOL INTEGER ID

62491

- 1 Prepare the three separate solutions needed for this extraction protocol as follows:
  - lipid extraction solvent: a 15 : 15 : 5 : 1 : 0.18 ratio by volume of 95% ethanol, water, diethyl ether, pyridine, and 4.2 N ammonium hydroxide, respectively.
  - water-saturated butanol: a 1:1 ratio of butanol and Milli-Q water
  - purified lab water: (Milli-Q water)
- 2 Concentrate cells via centrifugation using methods appropriate for your available lab equipment and for your species of interest. Decant supernatant. Transfer pelleted cells to a 2-mL centrifuge tube. Cells may need to be re-pelleted and residual supernatant removed after transfer to final 2-mL centrifuge tube.  
  
Note: Appropriate volume of lab culture or field samples to concentrate and extract depends on cell concentration. As a guideline, we have been successful concentrating 10-25 mL of lab cultures of *Microcystis aeruginosa* and ~50 mL of either raw lake water or mesocosm samples.
- 3 Add 1 mL of extraction solvent, ~100 µL of glass beads, and vortex ~5 s.
- 4 Incubate sample in 60 °C water bath for 20 min.
- 5 Centrifuge sample at 10,000 x g for 10 min.
- 6 Remove supernatant and place into a 1-dram glass vial (dram vial #1). The first two extractions from a sample will be placed in this vial (#1).
- 7 Repeat steps 3-6, except DO NOT ADD more glass beads.

- 8 Dry the collected supernatant in dram vial #1 under a stream of nitrogen.
- 9 Re-suspend dried sample in 300  $\mu$ L of water-saturated butanol and 150  $\mu$ L of Milli-Q water.
- 10 Vortex and transfer to a 2-mL centrifuge tube.
- 11 Centrifuge at 10,000 x g for 2 min.
- 12 Remove top butanol phase and place into a NEW 1-dram glass vial (dram vial #2).
- 13 Wash original dram vial (#1) with 300  $\mu$ L saturated butanol and transfer to residual aqueous phase in 2-mL centrifuge tube from step 10. Vortex.
- 14 Centrifuge at 10,000 x g for 2 min. Remove top butanol phase and place into dram vial #2.
- 15 Dry the collected butanol phase in dram vial #2 under a stream of nitrogen.
- 16 Re-suspend dried sample in 300  $\mu$ L of 9:1 methanol:chloroform.
- 17 The samples are now ready for analysis *via* LC/MS.