




APR 15, 2024

In vivo Electrophysiology Protocol

 In 2 collections

Sasha Burwell¹

¹Duke University

ASAP Collaborative Research Network



Sasha Burwell
Duke University

OPEN  ACCESS



ABSTRACT

This protocol details the procedure used to collect extracellular single unit dopamine recordings for this paper.

MATERIALS

- Intan RHD 16-channel headstage (C3335)
- Intan ultra-thin SPI interface cable (C3211)

DOI:

dx.doi.org/10.17504/protocols.io.14egn68eqI5d/v1

Protocol Citation: Sasha Burwell 2024. In vivo Electrophysiology Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn68eqI5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Apr 02, 2024

Last Modified: Apr 15, 2024


In vivo Electrophysiology Protocol

2h

- 1 Headfix mouse on a treadmill.
- 2 Attach prepared internal cannula (see DART Infusion protocol for details on how to set this up) to the implanted guide cannula.
- 3 Attach implanted electrode bundle to the Open Ephys Acquisition board using an Intan RHD 16-channel headstage (C3335) and an Intan ultra-thin SPI interface cable (C3211).
- 4 Open the Open Ephys GUI.
 - 4.1 Open Ephys GUI must have
 - (1) an Acquisition Board node followed immediately by
 - (2) a Recording Node, saving the raw data in a binary format.



- 4.2 Optional:** downstream in the processors graph, you can include a bandpass filter, common average reference node, LFP viewer, spike sorter, and spike viewer to view channels with spikes online during the recording. 1h 30m



1. Begin recording (press "Play" then "Record").
2. After 15 minutes of baseline recording, begin the infusion ($1.5 \mu\text{L}$ at $0.1 \mu\text{L} / \text{min}$, taking 15 minutes).
3. After the infusion is complete, continue recording for $01:30:00$, then stop the recording.
4. Advance the electrodes by 26 undetermined .
5. Return the mouse to its home cage to recover, and wait 2 weeks before recording again.