



Nov 05, 2020

Quantification of the molar concentration of the NGS libraries by qPCR

PLOS One

Yoko Kato-Unoki¹¹Center for Advanced Instrumental and Educational Supports, Faculty of Agriculture, Kyushu University, Fukuoka, Japan**1** Works for me dx.doi.org/10.17504/protocols.io.bmgck3sw

Yoko Kato-Unoki

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0240823>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Kato-Unoki Y, Umemura K, Tashiro K (2020) Fingerprinting of hatchery haplotypes and acquisition of genetic information by whole-mitogenome sequencing of masu salmon, *Oncorhynchus masou masou*, in the Kase River system, Japan. PLoS ONE 15(11): e0240823. doi: [10.1371/journal.pone.0240823](https://doi.org/10.1371/journal.pone.0240823)

DOI

dx.doi.org/10.17504/protocols.io.bmgck3sw

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0240823>

PROTOCOL CITATION

Yoko Kato-Unoki 2020. Quantification of the molar concentration of the NGS libraries by qPCR.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bmgck3sw>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Kato-Unoki Y, Umemura K, Tashiro K (2020) Fingerprinting of hatchery haplotypes and acquisition of genetic information by whole-mitogenome sequencing of masu salmon, *Oncorhynchus masou masou*, in the Kase River system, Japan. PLoS ONE 15(11): e0240823. doi: [10.1371/journal.pone.0240823](https://doi.org/10.1371/journal.pone.0240823)

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0240823>

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 17, 2020

LAST MODIFIED

Nov 05, 2020

PROTOCOL INTEGER ID

42212

MATERIALS TEXT

MATERIALS TEST

MATERIALS

⊗ Agencourt AMPure XP Beckman

Coulter Catalog #A63880

⊗ KAPA Library Quantification Kit for Illumina® Platforms Kapa

Biosystems Catalog #KK4835

⊗ Agarose S Nippon

Gene Catalog #312-01193

BEFORE STARTING

Note: use low-binding tubes and pipet chips.

- 1 Using the NGS library concentration obtained from the Agilent BioAnalyzer, dilute the libraries to 100 pg/μl with 10 mM Tris-HCl, pH 8.0. Then prepare a 1:10000 dilution of each diluted library (final concentration, 0.01 pg/μl) by conducting 10-fold dilution series using 10 mM Tris-HCl, pH 8.0.

*If the library concentration is unknown, conduct a 10-fold dilution series to achieve an appropriate final dilution (Note: for a 1:3 diluted library, we conducted a 1:500,000 dilution with 10 mM Tris-HCl, pH 8.0).

- 2 Prepare a reaction mix and performs qPCR in accordance with the manufacturer's manual of KAPA Library Quantification kit (KAPA Biosystems).
- 3 Check the product size and presence of adaptor dimers in a 1-μl aliquot of the PCR products by conducting electrophoresis through a 3% agarose-TBE gel. For example, in the Fig. 1, the lanes marked "x" show contamination with adaptor dimer.




Fig. 1. Example agarose electrophoresis.

- 4 Further clean-up the library with Agencourt AMPure XP beads (Beckman Coulter), if the library contains adaptor dimer. Refer to the Agencourt AMPure XP beads manual for reagent amounts and clean-up procedures. Elute to the same final volume as the original library and re-quantify and check again (i.e, repeat steps 1–3).
- 5 Calculate the molar concentration of each libraries from the standard curve.