



Jun 28, 2022

Chloramphenicol 25 mg/mL Stock Solution

Bailey Clark¹¹Florida State University iGEM

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.eq2lyn7ewvx9/v1

FSU iGEM

Bec18

ABSTRACT

The abstract will be added later.

DOI

dx.doi.org/10.17504/protocols.io.eq2lyn7ewvx9/v1

PROTOCOL CITATION

Bailey Clark 2022. Chloramphenicol 25 mg/mL Stock Solution. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2lyn7ewvx9/v1>

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 18, 2022

LAST MODIFIED

Jun 28, 2022

PROTOCOL INTEGER ID

58447

Preparing the Stock Solution

1 Weigh **25 mg** (**0.025 g**) of Chloramphenicol.

- 2 Add the chloramphenicol to a 1.5 mL tube.
- 3 Add 100% ethanol to the 1mL level on the tube.
- 4 Vortex the tube for 5 seconds, repeat 3 times.
- 5 Store the tube in the -20°C freezer.