



2 ▼

Feb 02, 2022

Digestion for BioBrick Assembly Kit (E0546) V.2

[New England Biolabs¹](#)¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.bddsi26e**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This protocol explains methods for assembling multi-component genetic systems using BioBrick[®] parts.

DOI

dx.doi.org/10.17504/protocols.io.bddsi26e<https://www.neb.com/protocols/0001/01/01/digestion-protocol-e0546>

New England Biolabs 2022. Digestion for BioBrick Assembly Kit (E0546).

protocols.io<https://dx.doi.org/10.17504/protocols.io.bddsi26e>[New England Biolabs](#)

Biobrick, Assembly, E0546

_____ protocol ,


Mar 08, 2020

Feb 02, 2022

33938

The BioBrick® Assembly Kit provides a streamlined method for assembly of BioBrick parts into multi-component genetic systems. BioBrick parts are DNA sequences that encode a defined biological function and can be readily assembled with any other BioBrick part. The process for assembling any two BioBrick parts is identical and results in a new composite BioBrick part.

The BioBrick Assembly Kit contains EcoRI-HF®, XbaI, SpeI, PstI, T4 DNA Ligase and NEBuffer 2.1.



BioBrick Assembly Overview

MATERIALS

[BioBrick Assembly Kit - 50 rxns New England](#)

Biolabs Catalog #E0546S

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

The BioBrick® Assembly Kit was developed in partnership with Ginkgo BioWorks. What follows is an abbreviated set of protocols for the use of the BioBrick® Assembly Kit (to assemble an Upstream Part with a Downstream Part into a Destination Plasmid). For more details and for technical questions, please see [here](#).

Digest Upstream Part with EcoRI-HF[®] and SpeI:

A	B
Reagent	Volume
Upstream Part Plasmid	500 ng
EcoRI-HF	1 µl
SpeI	1 µl
10X NEBuffer 2.1	5 µl
H2O	to 50 µl

2 

Digest Downstream Part with XbaI and PstI:

A	B
Reagent	Volume
Upstream Part Plasmid	500 ng
XbaI	1 µl
PstI	1 µl
10X NEBuffer 2.1	5 µl
H2O	to 50 µl

3 

Digest the Destination Plasmid with EcoRI-HF[®] and PstI:

A	B
Reagent	Volume
Destination Plasmid DNA	500 ng
EcoRI-HF	1 µl
PstI	1 µl
10X NEBuffer 2.1	5 µl
H2O	to 50 µl

The Destination Plasmid DNA should either be prepared with PCR or contain a toxic gene (e.g. ccdB, sacB) in the cloning site to avoid the need for gel purification. The Destination Plasmid should also have a different antibiotic resistance marker from both the plasmid containing the Upstream Part and the plasmid containing the Downstream Part to avoid the need to purify the Upstream and Downstream Parts.

4



Incubate all three restriction digest reactions at **37 °C** for **00:10:00**.

5



Heat inactivate at **80 °C** for **00:20:00**.