



Jun 22, 2022

# Addition of Antibiotic Supplement to Media

Mariam Alkattan<sup>1</sup><sup>1</sup>Fulbright Program*In Development*

Share

This protocol is published without a DOI.



Mariam Alkattan

## DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## ABSTRACT

The following protocol is for the addition of antibiotic stock solution to molten media.

## PROTOCOL CITATION

Mariam Alkattan 2022. Addition of Antibiotic Supplement to Media .  
**protocols.io**  
<https://protocols.io/view/addition-of-antibiotic-supplement-to-media-cbsksncw>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 22, 2022

LAST MODIFIED

Jun 22, 2022

PROTOCOL INTEGER ID

65068

MATERIALS TEXT

### Equipment

- thermometer; laser infrared (preferred) or stem
- adjustable pipette
- pipette tips

### Consumables

- antibiotic stock solution

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 Retrieve molten media from autoclave using heat resistant gloves. Place media in laminar flow hood if one is available. If laminar flow hood is not available, place media next to open flame on lab bench.



To further protect from burn injuries, keep heat resistance gloves on and use box or bin to transport hot media if you plan to move media more than a few feet from autoclave.

- 1 Wait for media to cool to **50 °C**. Measure temperature using laser infrared thermometer or a stem thermometer. If using stem thermometer, sterile using 70% ethanol and allow to dry right before use.

- 2 When media reaches 🔥 50 °C , use an adjustable pipette to add appropriate volume of antibiotic stock solution to molten media.
- 3 Gently swirl media in bottle to mix antibiotic.

If media is being used to pour plates, pour plates immediately after last step.