



May 23, 2022

Labeled microtubules for single-molecule imaging

John Salogiannis¹¹Department of Cellular and Molecular Medicine, University of California, San Diego, La Jolla, CA 92093

1

dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1 Mariusz Matyszewski

This protocol is for making labeled taxol-stabilized microtubules to be used for single-molecule imaging assays by adhering to biotin slides.

DOI

dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1

John Salogiannis 2022. Labeled microtubules for single-molecule imaging.

protocols.io<https://dx.doi.org/10.17504/protocols.io.bp2l6bdedgqe/v1>

single-molecule, microtubule, imaging, ASAPCRN

 protocol ,

Jan 19, 2021

May 23, 2022

46496

Required Buffers:

5x BRB80:

- [M]400 millimolar (mM) PIPES, pH 6.8 with KOH
- [M]10 millimolar (mM) MgCl₂
- [M]5 millimolar (mM) EGTA

2x Polymerization Mix:

- [M]2 x BRB80
- [M]2 millimolar (mM) DTT
- [M]2 millimolar (mM) GTP (Add last)
- [M]2 millimolar (mM) MgCl₂ (Add second to last)
- [M]20 % DMSO
- Mix well between adding each ingredient


For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Please take notice of the buffer preparation in the Materials section

Create taxol-stabilized labeled microtubules (can be reused for multiple weeks):

50m

1 

In a prechilled 1.5 mL Eppendorf tube, make a 10 µL mixture of [M]10 mg/mL tubulin . The mixture should be 80% unlabeled, 10% biotin-tubulin, and 10% 405-tubulin. Mix by gently flicking.

2 Let it sit  On ice for 00:10:00 .

10m

3 

Add equal volume of [M]2 x polymerization buffer (10 µL). Mix by gently flicking.

4 

30m

Incubate at  37 °C for 00:30:00 . Make a [M]1 x BRB80 + [M]1 millimolar (mM) DTT

+ [M] **20 micromolar (μ M) Taxol** stock and incubate at **37 °C** at the same time.

5 

Add equal volume of prewarmed [M] **1 x BRB80** + DTT + Taxol (**20 μ L**). Mix by gently flicking.

6 

10m

Incubate at **37 °C** for at least **00:10:00** (solution will be stable for hours at this point).

7 Store in the dark at **Room temperature** . Should be usable for several weeks, but more aggregates will appear over time.