



Sep 07, 2020

XPRIZE SANATA Protocol for Swab LFIA Test

Mario Thomas¹, Jasmine Sollen¹, Natalia Ivanova¹, Heidi Abdilla¹, Reda Fayek¹, Michelle Feng¹, Stephanie Lim¹, Amanda Naaum¹

¹Precision Biomonitoring Inc.

1

Works for me

dx.doi.org/10.17504/protocols.io.bkzjx4n

XPRIZE Rapid Covid Testing

SANATA

heidi.abdilla

ABSTRACT

This procedure outlines the protocol using any swab suitable for a SARS-CoV-2 sample collection. The purpose of this test is to detect low levels of SARS-CoV-2 antigen at a higher sensitivity. Precision Biomonitoring Inc. developed an ultra-rapid digital, disposable, highly-sensitive and inexpensive testing device used for screening purposes. The mobile app complementary to this medical device is connected through Bluetooth. Using this innovation, the patients can be tested at the point-of-care (POC), and obtain qualitative results. It is recommended that positive results should be confirmed by RT-PCR.

DOI

dx.doi.org/10.17504/protocols.io.bkzjx4n

PROTOCOL CITATION

Mario Thomas, Jasmine Sollen, Natalia Ivanova, Heidi Abdilla, Reda Fayek, Michelle Feng, Stephanie Lim, Amanda Naaum 2020. XPRIZE SANATA Protocol for Swab LFIA Test. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bkzjx4n>

KEYWORDS

Lateral Flow, Covid-19, SARS-CoV-2, Antigen Testing, Pandemic, Global Pandemic, Virus

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 07, 2020

LAST MODIFIED

Sep 07, 2020

PROTOCOL INTEGER ID

41739

GUIDELINES

- The digital device kits should be stored at room temperature and should never be exposed to extreme temperatures.
- The digital device kits are for single use. Do not reuse the kits. Dispose of all used materials in a biohazardous waste container.
- It is recommended that positive results should be confirmed by RT-PCR.

MATERIALS

NAME	CATALOG #	VENDOR
SWAB Lysis Tube	SWALT202001	Precision Biomonitoring Inc.
Disposable Graduated Transfer Pipette	13-711-9AM	Fisher Scientific
Flocked Swab	22-025-195	Fisher Scientific
LFIA Testing Device	N/A	Precision Biomonitoring Inc.

MATERIALS TEXT

Bluetooth Smartphone

EQUIPMENT

NAME	CATALOG #	VENDOR
Flocked Swab	22-025-195	
Disposable Graduated Transfer Pipette	13-711-9AM	
LFIA device	N/A	
SWAB Lysis Tube	SWALT202001	

SAFETY WARNINGS

When working with human bodily fluids, pathogens may be present. Wear the correct personal protection equipment (ie. disposable gloves) and wash your hands immediately after removing the gloves. Dispose of materials in biohazardous waste containers.

Setting up the test

10s

1

10s



Ensure the smartphone is fully charged and Bluetooth on the smartphone is turned on. The mobile app should be downloaded and ready to run.



When working with human bodily fluids, pathogens may be present. Wear the correct personal protection equipment (ie.gloves) and wash hands immediately after removing the gloves.

Collecting and Preparing Sample

1m

2

Collect sample using any swab suitable for SARS-CoV-2 sample collection. For example:

30s



Flocked Swab


Fisher Scientific

22-025-195

3

Rinse the used swab in

30s



SWAB Lysis Tube

Precision Biomonitoring Inc.


SWALT202 001

by submerging the swab in the lysis buffer and rotating 5-7 times against the tube walls. Discard the swab in a biohazardous waste container.

Using Testing Device 13m 10s

4 Using

1m




Disposable Graduated Transfer Pipette

Fisherbrand

13-711-9AM

apply  **80 µl** (3 drops) of the contents of the SWAB Lysis Tube into the sample port of the




LFIA device

Precision Biomonitoring Inc.

N/A

5 Place device on a flat surface. Let sample mixture run undisturbed for  **00:12:00** at  **Room temperature** .^{12m}

6 Ensure the smartphone is connected to the testing device through Bluetooth and read the results using the mobile phone app.^{10s}



The result will appear on screen as positive, negative, or inconclusive. If the result is inconclusive, conduct another test with a new device.