

May 03, 2024

Patient PBMC collection and cryopreservation and cryorecovery

DOI

dx.doi.org/10.17504/protocols.io.n92ld895ov5b/v1

Rebecca Wallings¹

¹University of Florida



Rebecca Wallings

University of Florida

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.n92ld895ov5b/v1**

Protocol Citation: Rebecca Wallings 2024. Patient PBMC collection and cryopreservation and cryorecovery. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ld895ov5b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 24, 2024

Last Modified: May 03, 2024

Protocol Integer ID: 98735

Abstract

Isolation of PBMC from human blood using CPT tubes, cryopreservation and cryorecovery of cells

PBMC collection and cryopreservation

- 1 Collect 8 mL of blood per BD Vacutainer CPT Cell Preparation Tube with Sodium Citrate (BD Biosciences, 362761).

Invert tubes 8-10 times and centrifuged at room temperature at $1500\times g$ for 20 minutes with no brake.

Using a 10mL stripette , transfer the PBMC enriched layer to a new 50-mL conical tube, being sure not to touch the ficoll layer, and add MACS buffer (PBS, 0.5% bovine serum albumin, 20 mM EDTA, pH 7.2) to a final volume of 50 mL, followed by centrifugation at $1800\times g$ for 10 minutes at room temperature.

Following removal of the supernatant, resuspend PBMCs in 10mL MACS buffer and count using a hemocytometer with Trypan blue at a 1:20 dilution (10uL of cells mixed with 10uL of 1:20 Trypan blue solution).

Centrifuge PBMCs for 10 min at $300\times g$ at 4°C .

Aspirate supernatant and gently resuspend cell pellet in cryopreservation media (RPMI 1640, with 20% FBS and 10% DMSO) at a final concentration of 1×10^7 cells/mL in cryovials (Simport, T311-2).

Place cryovials in a room-temperature Mr. Frosty freezing container with isopropanol as per manufacturer's instructions and store at -80°C overnight.

Remove cryovials from freezing containers and immediately placed into liquid nitrogen for long-term storage.

Cryorecovery

- 2 For cryorecovery, retrieve PBMCs from liquid nitrogen, and thaw at 37°C in water bath.

Once thawed (roughly 1 minute with agitation) slowly add contents of one cryovial to 25mL 37°C filter sterilized complete culture media (RPMI 1640 media, 10% heat-inactivated FBS, 1mM Penicillin-Streptomycin) in 50mL falcon and pellet via centrifugation at $300\times g$ for 10 min at room temperature.

Remove supernatant and resuspend cells in 37°C complete culture media for cell counting using Trypan blue.