


Apr 15, 2021

MIBI and IHC solutions

Forked from [IHC staining](#)Marc MB Bosse¹, Christine Camacho¹, Sean Bendall¹, Mike Angelo¹¹Department of Pathology, Stanford University1 Works for me dx.doi.org/10.17504/protocols.io.bhmej43eHuman BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu Marc Bosse
Stanford University

ABSTRACT

These protocols are designed to prepare solutions to be used for the Multiplex Ion Beam Imaging Time of Flight instrument (MIBI_TOF) and IHC related protocols in the Bendall and Angelo labs. Particularly, the IHC protocol uses the same reagents and steps used for MIBI staining. Therefore, It can be used to validate the compatibility of a new target antibody before use for MIBI staining. These protocols have been successfully use for MIBI and are the result of extensive optimization experiments.

DOI

dx.doi.org/10.17504/protocols.io.bhmej43e

PROTOCOL CITATION

Marc MB Bosse, Christine Camacho, Sean Bendall, Mike Angelo 2021. MIBI and IHC solutions. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bhmej43e>

FORK NOTE

FORK FROM

Forked from [IHC staining](#), Marc Bosse

KEYWORDS

IHC, Immunohistochemistry, antibody, MIBI

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 17, 2020

LAST MODIFIED

Apr 15, 2021

PROTOCOL INTEGER ID

38278

GUIDELINES

Use sterile single use 50 mL tube for measurement

Use the measurement marks on the Nalgene bottle to complete with water

Do not use any other vessels then single use graduated plastic containers to prevent cross contamination

MATERIALS TEXT

A	B	C
Products	Provider	Catalogue No.
TBS IHC Wash Buffer with Tween 20	Cell Marque	935B-09
PBS IHC Wash Buffer with Tween 20	Cell Marque	934B-09
Target Retrieval Solution, pH 9, (3:1)	Agilent (Dako)	S2375
Tris Base, BioReagents	Fisher Scientific	BP154-1
UltraPure water	Invitrogen	10977-015
Ultrapure (type 1, >18 MOhms) water	In house	NA
Hydrochloric acid	Sigma-Aldrich	H1758-100ML
Gelatin (cold water fish skin)	Sigma-Aldrich	G7765-250
Bovine Albumin (BSA), heat shock treated	Fisher Scientific	BP1600-100
Horse serum	Vector Labs	S-2000
VectaMount Permanent Mounting Medium	Vector Labs	H-5000
Nalgene™ Rapid-Flow™ Sterile Disposable Filter Units, 150 mL	Fisher Scientific	09-740-28E
Nalgene™ Rapid-Flow™ Sterile Disposable Filter Units, 500 mL	Fisher Scientific	09-740-28C
Nalgene, bottle 1L	Fisher Scientific	03-312-5
Millex-HV (0.45µm)	Millipore	SLHV033RS
BD Syringes without Needle, 50 mL	Fisher Scientific	13-689-8
Brady label White polyester/polypropylene	Brady	M21-500-7425
BMP21-LAB Printer Kit	Brady	BMP21-LAB

SAFETY WARNINGS

All organic solvents should be manipulated under a chemical hood

BEFORE STARTING

Verify the stocks of all reagents prior starting an experiment

If some reagents are running low, ask to place an order or prepare solutions

Target retrieval solution (Antigen retrieval)

- 1 Make a fresh 1:10 dilution
▢ **2.5 mL** of target retrieval solution 10x (3-in-1) pH 9 (Dako)
in ▢ **22.5 mL** of ultrapure (type 1, >18 MOhms) water

Total volume (mL)	Volume (mL) Target retrieval	Volume (mL) H2O
25	2.5	22.5
50	5	45
100	10	90

MIBI or IHC 1X PBS wash buffer

- 2 Verify stock of **MIBI or IHC 1x PBS wash buffer** and prepare accordingly if running low

- 2.1 Use the same 1 L Nalgene Bottle, if before expiring date. Replace by a new 1 L Nalgene bottle on monthly basis

Use single use 50 mL tube for measurement

Use the measurement marks on the bottle to complete with water

Do not use any other vessels than single use graduated plastic containers to prevent cross contamination

- 2.2 **Label name** (Brady label 0.5", M21-500-7425, Font 14 pts bold, format: Banner landscape)

First line: **MIBI or IHC 1X PBS wash buffer**

Second line: **Exp.: MM/DD/YY** (Expiring date, 1 month from the present date)

- 2.3 Prepare 1x PBS wash buffer as follow:

Add reagent in row 2 and row 3 first, mix gently

than complete with row 4, ultrapure (type 1, >18 MOhms) water

A	B
Reagents	Qty for 1000 mL
PBS IHC Wash Buffer with Tween 20 (mL)	50
Bovine Albumin (BSA), heat shock treated (g)	1
Ultrapure (type 1, >18 MOhms) water (mL)	949

- 2.4 Store at 4°C

Reuse the same Nalgene bottle to make new 1x PBS wash buffer until the initial expiring date (<1 month)

Antibody Diluent

- 3 Prepare Antibody Diluent (NHS 3%)

- 3.1 Prepare antibody diluent preparation buffer (1xTBS-Tween)

Add **100 mL** UltraPure Water (Invitrogen, 10977-015) in a Nalgene filter unit (Fisher Scientific, 09-740-28E)

3.2 Add the stock solution **7.5 mL** of TBS IHC Wash Buffer with Tween 20 (Cell Marque, 20x)

3.3 Complete to **150 mL** with UltraPure Water (Invitrogen, 10977-015)

3.4 Filter under vacuum

3.5 Label as Antibody diluent preparation buffer (1xTBS-Tween)

3.6 Make a 3% (volume/volume) of normal Horse serum

Add **1.5 mL Normal Horse Serum (Vector labs, S-2000)** in **48.5 mL** of Antibody diluent preparation buffer (1xTBS-Tween) in 50 mL tube

3.7 Filter with a 50 mL syringe (Fisher Scientific, 13-689-8) and filter (Millipore, Millex-HV, SLHV033RS) into a new 50 tube

3.8 Make some small aliquots by distributing 1 mL with a 10 mL serological pipette in 1.5 or 2 mL tubes and stock aliquots of 10 mL in 15 mL tubes

3.9 Label cap tubes using dot printed labels **Antibody Diluent (NHS 3%)**

3.10 Store at -20°C small aliquots in the staining box reagents and stock aliquots in the specified location

Gelatin 10% stock, 100x

4 Prepare stock 10% weight/ volume cold fish water gelatin

4.1 Warm up the gelatin stock bottle at **37 °C** prior use to reduce viscosity

Put a 15 mL tube on a precision balance and tare it to zero

4.2

4.3 Take gelatin by aspirating ~ **1 mL** using a large opening pipette tip >4 mm OD

4.4 Add slowly the viscous solution to the tube until it weigh **1 g**

4.5 Warm up **20 mL** Ultrapure Water (Invitrogen, 10977-015) in a 50 mL

4.6 Take **8 mL** of warm Ultrapure Water (Invitrogen, 10977-015) and add into the gelatin 15 mL

4.7 Make sure that all the gelatin is dissolved and complete the volume to **10 mL**

4.8 Make small aliquots by distributing ~1 mL with a 10 mL serological pipette in 1.5 or 2 mL tubes

4.9 Store at -20°C small aliquots in the staining box reagents

Triton 10% stock, 100x

5 Prepare Triton stock 10% volume/ volume

5.1 Add **9 mL** of Ultrapure Water (Invitrogen, 10977-015) to a **50 mL** tube

5.2 Take Triton by aspirating slowly **1 mL** using a large opening pipette tip >4 mm OD

5.3 Transfer into the **50 mL** tube

5.4 Make sure that all the Triton is dissolved

5.5 Make small aliquots by distributing ~1 mL with a 10 mL serological pipette in 1.5 or 2 mL tubes

5.6 Store at -20°C small aliquots in the staining box reagents

Sodium azide (20 mg/mL), 100x

6

Prepare 100x (2%) of sodium azide

6.1 Add 48 mL of Ultrapure Water (Invitrogen, 10977-015) in a 50 mL

6.2 Weigh 1 g of sodium azide

6.3 Add sodium azide in the 50 mL tube

6.4 Make sure that all the sodium azide is dissolved and complete the volume to 50 mL

6.5 Make small aliquots by distributing ~1 mL with a 10 mL serological pipette in 1.5 or 2 mL tubes

6.6 Store at 4°C small aliquots in the reagents box

Blocking Buffer

7

Prepare the blocking buffer following the table:

A	B	C	D
Reagents	Qty for 10 mL	Qty for 30 mL	Qty for 50 mL
1xTBS IHC Wash Buffer with Tween 20 (mL)	9.5	28.5	47.5
Triton 10% (mL)	0.1	0.3	0.5
Gelatin 10% (mL)	0.1	0.3	0.5
Horse serum (mL)	0.2	0.6	1
Sodium azide 20 mg/mL (mL)	0.1	0.3	0.5

- 7.1 Filter the blocking buffer with a filter (Millipore, Millex-HV, SLHV033RS) into a new tube
- 7.2 Make small aliquots by distributing ~1 mL with a 10 mL serological pipette in 1.5 or 2 mL tubes
- 7.3 Store at -20°C small aliquots in the staining box reagents

PBS 10x low barium

- 8 Use a new 1 L Nalgene bottle
- 9 Add 900 mL of ultrapure (type 1, >18 MΩcm) water in 1 L Nalgene bottle

Use the measurement marks on the bottle

Do not use any other vessel to measure the water

- 10 Weight the following chemicals:

For 1 L add:

☐ **80 g** NaCl

☐ **2 g** KCl

☐ **26.8 g** Na₂HPO₄·7H₂O

☐ **2.4 g** KH₂PO₄

Make sure that the chemical used correspond to the formula, particularly, for the sodium phosphate dibasic heptahydrate

- 11 Mix thoroughly. Do not adjust pH
- 12 Filter using 2 x 500 mL vacuum filter unit (Fisher Scientific, 09-740-28C)
- 13 Label each bottles with Brady lab printer (font 14, bold, banner landscape, tape 1/2") as:

First line: **PBS low barium 10x stock**

Second line: **Date (MM/DD/YY)**

Tris Buffer pH 8.5 100mM

- 14 Use a new 1 L Nalgene bottle, label with Brady lab printer (font 14pt, bold, banner landscape, 1/2" tape) as follow:

First line: **100 mM TRIS BUFFER pH8.5**

Second line: **Exp.: MM/DD/YY** (expiration date: add 1 month to the preparation date in a new bottle)
- 15 Add 900 mL of ultrapure (type 1, >18 MOhms) water in 1 L Nalgene bottle

Use the measurement marks on the bottle

Do not use any other vessel to measure the water
- 16 Add 12.114 g of Tris Base
- 17 Mix thoroughly
- 18 Adjust to 1000 mL with ultrapure (type 1, >18 MOhms) water
- 19 Add 3.0 mL of Hydrochloric acid (Sigma-Aldrich, H1758-100ML) and mix well
- 20 Calibrate pH meter
 - 20.1 Press Calibrate on the pH meter
 - 20.2 Wash the pH probe with ultrapure (type 1, >18 MOhms) water, wipe, submerge in buffer 1 (pH 10) and Press Read
 - 20.3 Wait until it stabilizes and says second buffer
 - 20.4 Wash the pH probe with ultrapure (type 1, >18 MOhms) water, wipe, submerge in buffer 2 (pH 7) and Press Read
 - 20.5 If the calibration is successful, the screen should indicate CALIBRATED

21 Adjust to pH 8.5

21.1 Wash the pH probe with ultrapure (type 1, >18 MOhms), wipe, submerge in the Tris solution

21.2 Press read and wait for the reading to stabilize

21.3 Acceptable pH values range between 8.4 to 8.6, if pH value is in that range, the solution is ready

21.4 If the pH is above 8.6 add 100 µL of HCl

21.5 If HCl has been added measure the pH again → **go to step #21.1**

Reuse the same Nalgene bottle to make new Tris Buffer pH 8.5 100mM until the initial expiring date (<1 month)

PBS low barium 1x

22 Use a new 1 L Nalgene bottle, label with Brady lab printer (font 14pt, bold, banner landscape, 1/2" tape) as follow:

First line: **PBS low barium 1x**

Second line: **Exp.: MM/DD/YY** (expiration date: add 1 month to the preparation date in a new bottle)

23 Add 900 mL of ultrapure (type 1, >18 MOhms) water in 1 L Nalgene bottle

Use the measurement marks on the bottle

Do not use any other vessel to measure the water

24 Use 50 mL disposable Falcon tube to measure  **100 mL** (2x50 mL) of PBS 10x Low Barium stock solution

25 Mix thoroughly

26 Reuse the same Nalgene bottle to make new PBS low Barium until the initial expiring date (<1 month)