




Jun 12, 2021

Expression and purification protocol of linear GST-4xUbiquitin

Chunmei Chang¹¹James Hurley Lab, UC, Berkeley

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvjbn4in

Chunmei Chang

ABSTRACT

This protocol details the expression and purification of linear GST-4x Ubiquitin.

ATTACHMENTS

[246-481.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.bvjbn4in

PROTOCOL CITATION

Chunmei Chang 2021. Expression and purification protocol of linear GST-4xUbiquitin. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bvjbn4in>

KEYWORDS

linear GST-4xUbiquitin, Protein purification, Protein expression

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 05, 2021

LAST MODIFIED

Jun 12, 2021

OWNERSHIP HISTORY

Jun 05, 2021 Urmilas

Jun 11, 2021 Chunmei Chang

PROTOCOL INTEGER ID

50499

General information:
[Sequence Analyzer: ATP1A1_plasmid_donor_RD Sequencing](#)
[Result addgene Catalog #171421](#)

A	B
Expression system	E.Coli BL21DE3
Medium	Luria Bertani
Plasmid origin	Addgene 171421
Backbone	pGEX5
Resistance	Amp
Insert	Mus Ubiquitin (NM_001313984.1)
Tags & cleavage sites	N-term GST no TEV cleavage site
Ext coeff	46060 M-1cm-1, MW 58.9 kDa

Lysis Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	300 mM
TCEP	1 mM
Protease Inhibitors (Roche)	

Wash Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	300 mM
TCEP	1 mM

Elution Buffer:

A	B
Hepes pH=7.5	50 mM
NaCl	300 mM
TCEP	1 mM

SEC Buffer:

A	B
Hepes pH=8	20 mM
NaCl	200 mM
TCEP	1 mM

Columns/Resin: - Glutathione Sepharose 4B
S200 (GE Helathcare)

Protein expression

4h 20m

1 Transform the E.Coli BL21DE3 cells with a plasmid encoding for GST-tetraubiquitin and plate them on Amp plate.

2  4h

Carry out the protein expression in **1.5 L** LB medium, induce with **100 Micromolar (μM)** IPTG (isopropyl-β-D-thiogalactopyranoside) to an OD₆₀₀ of 0.8 and grow at **18 °C** **Overnight**.

3  20m

Harvest the cells by spinning at **4500 x g** for **00:20:00** at **4 °C** and stock at **-80 °C** until purification.

Protein purification 5h

4 Follow the GST batch purification by Size Exclusion Chromatography.

5 Resuspend the pellets in Lysis Buffer, sonicate for cell lysis and clear at **16000 rpm** at **4 °C** for **01:00:00**^{1h}

6  4h

Incubate the supernatant with Glutathione Sepharose 4B (GE Healthcare) at **4 °C** with gentle shaking for **04:00:00**, apply to a gravity flow column, and wash extensively with Wash Buffer.

7 Elute the protein of interest with Elution Buffer and then apply onto a Superdex 6 column (10/300 Increase) pre-equilibrated in SEC Buffer operating at **4 °C**.

8 Pool the peak fractions containing pure protein, store the snap-frozen in liquid nitrogen at **-80 °C**.



