



Aug 26, 2020

COOMASSIE COLLOIDAL STAIN - CHEM 584

Forked from [VISUALIZING PROTEINS ON PAGE GEL WITH COOMASSIE COLLOIDAL STAIN - CHEM 584](#)Ken Christensen¹¹Brigham Young University

In Development

dx.doi.org/10.17504/protocols.io.bkcqkswvKen Christensen
Brigham Young University

ABSTRACT

Staining polyacrylamide gels using a colloidal Coomassie Blue stain. Adapted from a JoVE article.

DOI

dx.doi.org/10.17504/protocols.io.bkcqkswv

PROTOCOL CITATION

Ken Christensen 2020. COOMASSIE COLLOIDAL STAIN - CHEM 584. [protocols.io](https://dx.doi.org/10.17504/protocols.io.bkcqkswv)
<https://dx.doi.org/10.17504/protocols.io.bkcqkswv>

FORK FROM

Forked from [VISUALIZING PROTEINS ON PAGE GEL WITH COOMASSIE COLLOIDAL STAIN - CHEM 584](#), Ken Christensen

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2020

LAST MODIFIED

Aug 26, 2020

PROTOCOL INTEGER ID

41072

- 1 Wash gel with ~50 mL ddH₂O 3x for ~ 00:10:00 for each wash on a shaker
NOTE: insufficient washing causes poor sensitivity because the remaining SDS on the gels disturbs the bounding of the dye to the protein
- 2 Shake the colloidal Coomassie solution before use to disperse particles evenly.
- 3 Incubate the gels with the Coomassie solution by agitation on a shaker up to 12:00:00 . This staining generates marginal background so you can observe the staining progress.

Note: after 10 minutes you should see some protein bands appearing, within 2 hours of incubation about 80% maximum staining is completed. For nearly maximum sensitivity, an overnight incubation is recommended. In some cases, when the amount of protein to be stained is large, the solution turns a bright blue color. If this happens, you

should refresh the staining solution as necessary.

- 4 After staining, remove the Coomassie solution and rinse the gels twice with ddH₂O.

Note: you can reuse the staining solution as long as particles still remain, otherwise discard it appropriately.

- 5 Remove any dye particles from the staining dish with a lint-free paper towel and destain for up to 🕒01:00:00 in ddH₂O or a Coomassie destain solution (10% (v/v) 96% ethanol, 2% (v/v) orthophosphoric acid 85% in ddH₂O.)
- 6 Finally rinse the gels twice with ddH₂O: the gels will reach their original thickness (the alcohol-acid media makes the gels shrink). Neutralization in water also enhances the Coomassie stain's color intensity.
- 7 You can record your gel using the gel documentation system or scan it on the Licor using the 680/700 nm channel.