



Version 2

Dec 07, 2020

SDS-PAGE gel electrophoresis V.2

Steven J Burgess¹¹University of Illinois at Urbana-Champaign

In Development

This protocol is published without a DOI.

Ag SynBio Lab UIUC



Steven Burgess

University of Illinois at Urbana-Champaign

ABSTRACT

SDS-PAGE gel electrophoresis protocol for analyzing samples from field-grown tissue via immunofluorescence. In this protocol no Coomassie blue is added to samples, the reason is that this interferes with the fluorescent signal during immunoblot. Instead, samples have already been prepared in Laemmli buffer (minus coomassie, see [protein extraction procedure](#)), the leading edge of samples can be visualized due to the presence of chlorophyll.

Note

- When using 15 well, 0.75 mm comb, try to limit the volume loaded to 10 µL to minimize the risk of spillover of protein between wells.
- Ensure to wipe the tip on the rim of the sample tube to remove the sample stuck to outside of the tip.

Literature:

http://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6040.pdf<https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10026447.pdf>

PROTOCOL CITATION

Steven J Burgess 2020. SDS-PAGE gel electrophoresis. **protocols.io**
<https://protocols.io/view/sds-page-gel-electrophoresis-bqhsmt6e>
Version created by Steven Burgess

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2020

LAST MODIFIED

Dec 08, 2020

PROTOCOL INTEGER ID

45330

MATERIALS TEXT

- 4–20% Mini-PROTEAN[®] TGX[™] Precast Protein Gels, 15-well, 15 µL (Bio-Rad Laboratories; [4561096](#))
- Opening lever (Bio-Rad Laboratories; 456-0000)
- Chameleon[™] Duo Pre-stained Protein Ladder (LI-COR Biosciences; [NC0738562](#))
- Mini-PROTEAN[®] Tetra Vertical Electrophoresis Cell, 4-gel (Bio-Rad Laboratories; [1658004](#))
- 10x Tris/Glycine/SDS (Bio-Rad Laboratories; [1610732](#))
- Fisherbrand[™] Gel-Loading Tips, 1-200 µL (Fisher Scientific; [02-707-181](#))
- Imperial[™] Protein Stain (Thermo Scientific; [24615](#))s
- Bio-Rad Gel-Doc Imager (optional)

ABSTRACT

SDS-PAGE gel electrophoresis protocol for analyzing samples from field-grown tissue via immunofluorescence. In this protocol no Coomassie blue is added to samples, the reason is that this interferes with the fluorescent signal during immunoblot. Instead, samples have already been prepared in Laemmli buffer (minus coomassie, see [protein extraction procedure](#)), the leading edge of samples can be visualized due to the presence of chlorophyll.


Note

- When using 15 well, 0.75 mm comb, try to limit the volume loaded to 10 µL to minimize the risk of spillover of protein between wells.
- Ensure to wipe the tip on the rim of the sample tube to remove the sample stuck to outside of the tip.

Literature:

http://www.bio-rad.com/webroot/web/pdf/lsr/literature/Bulletin_6040.pdf<https://www.bio-rad.com/webroot/web/pdf/lsr/literature/10026447.pdf>

Prepare gel tank and buffers

- 1 Create a working dilution of Tris-Glycine running buffer (~  1 L is required per gel tank) by diluting 1:10 with d H₂O.
- 2 Carefully remove the comb from the precast gel and the tape across the bottom.
- 3 Assemble the Mini-PROTEAN electrophoresis cell and fill the inner chamber with buffer and the outer chamber up to the recommended mark

the volume varies depending on whether running 2 or 4 gels, the level is marked on the tank

- 4 Wash the wells with running buffer by pipetting up and down





This is done to remove residual acrylamide that may have collected in wells

Prepare Samples

10m

- 5 In fresh centrifuge tubes, create a dilution of each sample using 1x PEB, such that each sample is set at a concentration of 3 µg /mL of total soluble protein.

Recommended final volume ~  100 µl µL (this will allow for 10 samples) but will depend on the application

- 6 Heat samples at  95 °C for  00:05:00 . 5m
- 7 Spin down samples at  10000 x g for  00:05:00 to pellet debris. 5m

- 8 Load  3 µl of Chameleon™ Duo Pre-stained Protein Ladder to the first well

- 9 Load 10 μ L of each sample (30 μ g of total soluble protein) per lane.

Running Gel 10m

- 10 Run precast gels at 200 V for ~ ⌚ 00:30:00 .

30m

Or until the samples have reached the end of the gel. For self-made gels, run at 80-120 V

- 11 Carefully open precast gel case using an opening lever, by inserting where the black arrows indicate on the gel case.

- 12 Remove stacking gel with a blade

- 13 Proceed either directly to immunoblot, or if the gel is as a loading control, place in a container and cover with ImperialTM Protein Stain.

Visualizing Gel (optional) 10m

- 14 Gently agitate on a rocking platform for ⌚ 00:30:00 .

30m

- 15 Pour off the InstantBlue stain (collect as hazardous waste)

- 16 Rinse the gel with dH₂O.

- 17 Cover the stained gel with dH₂O and gently agitate in a rocking platform for ⌚ 00:30:00 .

30m

- 18 Repeat steps 15-17 until the background of the gel is clear and blue protein bands can be clearly visualized

- 19 Image the gel on the Bio-Rad Gel Doc system trans-white illumination