

Sep 10, 2021

Variant (E484K) ALERT - Ligation-Dependent Loop-Mediated Isothermal Amplification

Ali Bektas¹, Anitha Jayaprakash¹¹Oakland Genomics Center

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bxzjpp4n

Ali Bektas

ABSTRACT

A 2-hour, 2-temperature protocol, using RNA templated DNA ligation, for the visual detection of the E484K mutation of concern pertaining to SARS-CoV-2.

N.B. This is the first version of this protocol, stay tuned for increased sensitivity and multiplexing (A.B.)

DOI

dx.doi.org/10.17504/protocols.io.bxzjpp4n

PROTOCOL CITATION

Ali Bektas, Anitha Jayaprakash 2021. Variant (E484K) ALERT - Ligation-Dependent Loop-Mediated Isothermal Amplification. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bxzjpp4n>

FUNDERS ACKNOWLEDGEMENT

Just One Giant Lab
Grant ID: Covid19 Microgrant

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 05, 2021

LAST MODIFIED

Sep 10, 2021

PROTOCOL INTEGER ID

53003

Ligation Primers

- 1 Oligonucleotides used for RNA templated DNA ligation using
[SplintR Ligase - 1,250 units New England Biolabs Catalog #M0375S](#)

All sequences shown 5' to 3. Underlined segments are complementary to the target and bold nucleotide is specific to the G>A substitution at position 23012 resulting in the E484K amino acid change.

E484K-Donor
/5Phos/AACACCATTACAAGCTACGTACATGACATAATCCAAACTCATAAATTCTCGCATTTTAGATCCGTCTCCTTACAGGA
CACATCATCC

E484-Acceptor
CATCTCTAACTCTACTAAGACTTCCATTATCAACAATAGCTGACATGTTCTAGCAGCGACAGGACACACGACAGGAGGAAA
GTAACAATTAACCTTT

Ligation

- 2
 - Prepare ligation mix to use with previously isolated RNA. A standard Trizol/Chloroform method or column purification method for RNA isolation works well.

X1 Ligation mix (multiply for number of samples processed)

▢ 2 µl E484K-Donor (100nM)

▢ 2 µl E484K-Acceptor (100nM)

▢ 1 µl SplintR Buffer (NEB)

⊗ SplintR Ligase - 1,250 units **New England**

▢ 0.2 µl **Biolabs Catalog #M0375S**

(NEB)

▢ 1.8 µl molecular biology grade water

- add ▢ 3 µl RNA per ▢ 7 µl of reaction mix for a ▢ 10 µl reaction volume

- 3 Incubate at ⚡ 37 °C for ⌚ 00:20:00 followed by a ⌚ 00:20:00 at ⚡ 65 °C
inactivation step.

40m

- 4 Place ligation reactions on ice.

LAMP Primers

- 5 Oligonucleotides used for Loop-Mediated Isothermal Amplification
all sequences shown 5' to 3'

FIP

CTTTCCTCCTGTCGTGTGCTCCTCCATTATCAACAATAGCTGAC

BIP

ACCTTTAACACCATTACAAGCTACGACGGATCTAAAATGCGAGAA

F3

CATCTCTAACTCTACTAAGACT

B3

GGATGATGTGCTCCTGTAAGG

LAMP 50m

- 6
 - Prepare LAMP mix.

X1 LAMP mix (multiply for number of samples processed)

2.5 µl Isothermal Buffer 1 (NEB)

Magnesium Sulfate (MgSO₄) Solution - 6.0 ml New England

0.5 µl Biolabs Catalog #B1003S

(100mM)

3.5 µl dNTPs (10µM)

0.8 µl FIP (50µM)

0.8 µl BIP (50µM)

0.5 µl F3 (10µM)

0.5 µl B3 (10µM)

SYTO 9 fluorescent nucleic acid stain Life

0.1 µl Technologies Catalog #S-34854

(100µM)

Bst 2.0 WarmStart DNA Polymerase - 8,000 units New England

1 µl Biolabs Catalog #M0538L


11.8 µl Molecular Biology Grade water

- add 3 µl of the ligation product from the Ligation step to 22 µl of LAMP mix for a 25 µl reaction.

7 Incubate at 65 °C for 00:50:00

50m

8 Visualize reaction tubes over a UV transilluminator. Samples containing a G>A mutation at position 23012 will exhibit fluorescence.



RNA samples are synthetic controls from Twist Biosciences, at various (approximate) copy numbers