



MAR 07, 2024

Tissue culture of HeLa cells

In 1 collection

Dan Tudorica¹

¹University of California, Berkeley



Dan Tudorica
Hurley Lab, QB3, UC Berkeley

ABSTRACT

Growth in DMEM supplemented with FBS and antibiotics

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.kxygx31oog8j/v1

Protocol Citation: Dan Tudorica 2024. Tissue culture of HeLa cells. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygx31oog8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 09, 2023

Last Modified: Mar 07, 2024

- 1 Medium: To maintain HeLa cells, grow in Dulbecco's Modified Eagle Medium supplemented with 10% FBS and 1x Pen-strep. Grow in a sterile incubator at 37 C and 5% CO₂. Cells double roughly every 24 H. Grow on standard tissue culture-treated plasticware
- 2 Splitting instructions: Split 1:4 on Mondays and Wednesday, and 1:8 on Fridays. Split no more than 1:8.