



Apr 03, 2021

# RNA Extraction Protocol for Shorea

mpfsum<sup>1</sup><sup>1</sup>University of Malaya

1

Works for me

This protocol is published without a DOI.

MP2Lab



mpfsum

## ABSTRACT

RNA extraction protocol using CTAB method optimized for leaf and bud samples from *Shorea* sp.

## PROTOCOL CITATION

mpfsum 2021. RNA Extraction Protocol for Shorea. **protocols.io**  
<https://protocols.io/view/rna-extraction-protocol-for-shorea-btw9nph6>



## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Apr 03, 2021

## LAST MODIFIED

Apr 03, 2021

## PROTOCOL INTEGER ID

48833

### Lysis 1h 10m

- 1 Weigh BUD **40 mg (25-40 mg)** or LEAF **60 mg (50-70 mg)** and transfer into a 1.5 mL tube. 15m
- 2 Ground in CTAB buffer (3% CTAB, 1.4 M NaCl, 20 mM EDTA, 100 mM Tris-HCl pH 8.0, 0.2% b-mercaptoethanol) using <sup>10m</sup>tissuelyzer.
- 3 Incubate at **60 °C** for **00:45:00**. 45m

### Precipitation 25m

- 4 Add 1 volume of chloroform.  
Mix by vortex but not too vigorous.  
Ensure that the two phases of the mixture are homogenized well. 5m
  - 5 Centrifuge at **5000 rpm, Room temperature, 00:01:00**.  
Transfer the supernatant into a new 1.5 mL tube. 5m
- Add 2/3 volume of isopropanol. 5m

6 Mix gently (turning the tubes upside down ⌚ 00:02:00 ).

10m

7 Centrifuge ⌚ 5000 rpm, Room temperature , 00:05:00  
Discard solution carefully.



Formation of white to translucent pellet.

#### Washing 25m

5m

8 Add 🧴 700 µl ethanol.  
Mix gently (flicking the tube) until the pellet no longer stick to the tube.

10m

9 Centrifuge ⌚ 5000 rpm, Room temperature , 00:05:00 .  
Discard solution carefully.

5m

10 Air dry the pellet ⌚ 00:01:00 .  
Repeat Washing if necessary.

#### Resuspension 25m

11 Resuspend the pellet in 🧴 50 µl RNase-free water.  
Flick the tube until the pellet is no longer visible.  
  
If the pellet is difficult to dissolve, leave it for ~30 min before flicking again.  
If it still does not dissolve, add more RNase-free water AND adjust the subsequent steps accordingly.

#### Genomic DNA Digestion 55m

12 [🔗 Turbo DNA-free Kit Invitrogen - Thermo](#)  
Removed gDNA using [Fisher Catalog #AM1907](#)  
Follow the manufacturer's protocol.

#### Purification 40m

5m

13 [🔗 RNeasy Plant Mini](#)  
Purification is done using [Kit Qiagen Catalog #74904](#) with minor adjustment to the protocol.  
Add 🧴 350 µl RLT buffer.  
Mix gently (flicking).














5m

14 Add 🧴 250 µl absolute ethanol.  
Mix gently (flicking).

3m

15 Transfer the solution into RNeasy Mini Spin column.  
Centrifuge ⌚ 10000 rpm, Room temperature , 00:00:15 .

Discard the flow through.

- 16 Add  **700 µl** RW1 buffer to the column. 3m  
Centrifuge  **10000 rpm, Room temperature , 00:00:15** .  
Discard the flow through.
- 17 Add  **500 µl** RPE buffer to the column. 3m  
Centrifuge  **10000 rpm, Room temperature , 00:00:15** .  
Discard the flow through.
- 18 Add  **500 µl** RPE buffer to the column. 5m  
Centrifuge  **10000 rpm, Room temperature , 00:02:00** .  
Discard the flow through.
- 19 Replace the collection tube. 3m  
Centrifuge  **14000 rpm, Room temperature , 00:01:00** to further dry the membrane.
- 20 Replace the collection tube with a new 1.5 mL tube. 8m  
Add  **50 µl** RNase-free water directly onto the membrane.  
Incubate at  **Room temperature** for  **00:05:00** .
- 21 Centrifuge  **10000 rpm, Room temperature , 00:01:00** . 5m  
Store the extracted RNA in  **-20 °C** or  **-80 °C** for longer storage.

#### Quality Assessment

- 22 Evaluate the RNA quality using Agilent 2100 Bioanalyzer.