



Dec 03, 2020

Large Volume Immunostaining for Cleared Samples

Seth Currlin¹, Marda Jorgensen¹, Jerelyn Nick¹¹University of Florida

1

Works for me

dx.doi.org/10.17504/protocols.io.bprtm6n

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu

Seth Currlin

ABSTRACT

This is a guide for immunostaining CLARITY-processed samples. dx.doi.org/10.17504/protocols.io.8jihuke

These steps are meant to be a guide for immunostaining large samples and should be optimized to suit your particular tissues and reagents. Large tissue volumes and dense tissue types will require longer incubation and wash times. The parameters suggested below are for a piece of human tissue (thymus, spleen, or lymph node) approximately 5 mm³ in size.

A useful link: <http://wiki.claritytechniques.org/index.php/Immunostaining>

DOI

dx.doi.org/10.17504/protocols.io.bprtm6n

PROTOCOL CITATION

Seth Currlin, Marda Jorgensen, Jerelyn Nick 2020. Large Volume Immunostaining for Cleared Samples.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bprtm6n>

KEYWORDS

null, immunofluorescence, IHC

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 16, 2020

LAST MODIFIED

Dec 03, 2020

PROTOCOL INTEGER ID

44563

MATERIALS TEXT

Solutions:

Wash buffer (PBST): 1x phosphate buffered saline (PBS), 0.5% Triton X-100, pH 7.4; at least 20x volume of tissue.

Blocking buffer: 1x PBS, 1% horse serum, 0.5% Triton-X 100, 0.05% Tween-20, pH 7.4; at least 10x volume of tissue.
oAnother serum type should be used based on the secondary antibody.

Staining buffer: 1x PBS, 0.5% Triton X-100, 0.05% Tween-20, pH 7.4; at least 10x volume of tissue.

ABSTRACT

This is a guide for immunostaining CLARITY-processed samples. [dx.doi.org/10.17504/protocols.io.8jihuke](https://doi.org/10.17504/protocols.io.8jihuke)

These steps are meant to be a guide for immunostaining large samples and should be optimized to suit your particular tissues and reagents. Large tissue volumes and dense tissue types will require longer incubation and wash times. The parameters suggested below are for a piece of human tissue (thymus, spleen, or lymph node) approximately 5 mm³ in size.

A useful link: <http://wiki.claritytechniques.org/index.php/Immunostaining>

BEFORE STARTING

Ensure the tissue is sufficiently transparent for your imaging requirements. If the sample has not been successfully cleared then subsequent imaging will be difficult due to regions retaining light-obscuring lipids. Samples should be imaged prior to immunostaining to confirm adequate clearing, particularly among excitation/emission settings intended for later use.

Immunohistochemistry for large (5 mm³) tissue volumes

3w 5d

3d

1 Wash out residual clearing solution

Three washes in PBST, 24 hours each, 37°C with gentle agitation. Large wash volumes are best, a 50 ml conical tube with 50 ml fresh PBST is recommended for each wash step.

3d

2 Nonspecific Epitope Blocking

Blocking buffer incubation: three days, 37°C with gentle agitation. Blocking buffer volume should be at least 5 volumes greater than tissue volume.

Note: Serum type should be used based on the secondary antibody host species.

1w

3 Primary Antibody Incubation

1:100 dilution (~10 ug/mL) in Staining Buffer.

Primary antibody incubation: five days, 37°C with gentle agitation; then two days, 4°C with gentle agitation.

For best results a large volume for primary antibody incubation is recommended however. A 5 ml total volume for antibody staining is recommended for a 5 mm³ piece of tissue.

3d

4 Wash out Primary Antibody

Three washes, 24 hours each, at 37°C with gentle agitation.

1w

5 Secondary Antibody Incubation

1:200 dilution (~5 ug/mL) in Staining Buffer.

Secondary antibody incubation: five days at 37°C with gentle agitation; two days at 4°C with gentle agitation.

Note: Protect samples from light.

3d

6 Wash out Secondary Antibody

Three washes, 24 hours each, at 37°C with gentle agitation.

Note: Protect samples from light.

7 Prepare for 3D Imaging or Storage

Prepare for imaging or store samples in washing buffer (PBST) at 4°C.

For 3D imaging guide, see: <https://dx.doi.org/10.17504/protocols.io.begajbse>