

Jun 07, 2024

Production of lentiviruses in Lenti-X cells

DOI

dx.doi.org/10.17504/protocols.io.n92ldmdyol5b/v1

Cole S Sitron¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry

ASAP Collaborative Rese...



Cole Sitron

Max Planck Institute of Biochemistry

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.n92ldmdyol5b/v1

Protocol Citation: Cole S Sitron, F Ulrich Hartl 2024. Production of lentiviruses in Lenti-X cells. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldmdyol5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 16, 2024

Last Modified: June 07, 2024

Protocol Integer ID: 94535

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: ASAP-000282

Abstract

This is a protocol to produce lentivirus in a HEK cell line selected to produce high-titer virus.





Attachments



958-2489.docx

15KB

Materials

-  Opti-MEM™ Reduced Serum Medium **Thermo Fisher Scientific Catalog #31985062**
-  X-tremeGENE™ HP DNA Transfection Reagent **Merck MilliporeSigma (Sigma-Aldrich) Catalog #6366244001**
-  Lenti-X™ Concentrator **Takara Bio Inc. Catalog #631231**
-  Lenti-X™ 293T Cell Line **Takara Bio Inc. Catalog #632180**


Equipment

Milllex®-GP Filter Unit (Sterile)	NAME
Filter	TYPE
Milllex®	BRAND
SLGP033RB	SKU
https://www.merckmillipore.com/IN/en/product/Milllex-GP-Filter-Unit-Sterile,MM_NF-SLGP033RB ^{LINK}	


-  psPAX2 **addgene Catalog #12260**
-  pMD2.G **addgene Catalog #12259**
-  DMEM, high glucose, pyruvate **Thermo Fisher Catalog #11995073** supplemented with 10% FBS
( Fetal Bovine Serum **Gibco - Thermo Fischer Catalog #10270106**) and 1X Pen/Strep
( Penicillin-Streptomycin (10,000 U/mL) **Gibco - Thermo Fisher Catalog #15140122**)
- 1X PBS  PBS, pH 7.2 **Thermo Fisher Scientific Catalog #20012068**
- Appropriate transfer vector



Transfection

1 Plate Lenti-X cells 2.5×10^6 /10 cm  24:00:00 before transfection.

1d

2 Warm OptiMEM, X-tremeGENE HP DNA Transfection Reagent, and DNA to warm to  Room temperature .



3 Add  1 mL OptiMEM to a tube.





4 Add plasmid DNA to the OptiMEM, gently pipetting to mix.



A	B	C
psPAX2	0.65 pmol	4600 ng
pMD2.G	0.36 pmol	1400 ng
Transfer vector	0.82 pmol	

5 Add X-tremeGENE HP DNA Transfection Reagent to OptiMEM/DNA ( 4 μ L per  1 μ g of DNA).



6 Incubate for  00:15:00 at  Room temperature .

15m




7 Add OptiMEM/DNA/X-tremeGENE to the plate of cells dropwise, then gently swirl the plate.



8 The day after transfection, replace the medium with



















Collection

3d 1h 35m

9  72:00:00 after transfection, collect the medium in a 15 ml conical.

3d



- 10 Centrifuge  500 x g x  00:05:00 , then pass the supernatant through a  .22 μm filter.  5m
- 11 Add  3 mL of LentiX concentrator to the supernatant, invert 6 times to mix, then incubate  Overnight at  4 °C .    45m
- 12 Centrifuge at  1500 x g for  00:45:00 at  4 °C .  45m
- 13 Remove the supernatant and resuspend the pellet in  120 μL ice cold PBS. 
- 14 Make aliquots in multiples of  14 μL and freeze at  -70 °C . 