

Jun 13, 2024



## **Organoid Drug Treatment**

DOI

### dx.doi.org/10.17504/protocols.io.x54v92rb1l3e/v1

yzhu<sup>1</sup>

<sup>1</sup>Stanford University

**NCIHTAN** 



## Joanna Bi

Stanford University

# OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x54v92rb1l3e/v1

Protocol Citation: yzhu 2024. Organoid Drug Treatment. protocols.io <a href="https://dx.doi.org/10.17504/protocols.io.x54v92rb113e/v1">https://dx.doi.org/10.17504/protocols.io.x54v92rb113e/v1</a>

**License:** This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's

working

Created: June 13, 2024

Last Modified: June 13, 2024

Protocol Integer ID: 101794

### **Abstract**

Protocol for drug experiments of organoid samples



- 1 Trypsinize organoids and plate at 30,000 cells/well in 24-well plates
- 2 Incubate 5-7 days
- 3 Incubate organoids in a complete medium containing 500 nM JQ1 for 24 hours (JQ cells) or 0.1% v/v DMSO (control cells)
- 4 Harvest organoids in Cell Recovery Solution (Corning) on ice for 1 hour, wash with PBS, and centrifuge to retrieve cell pellets
- 5 Proceed to process cell pellets for RNA extraction