



# Artificial sea water V.2

Simon Blanchoud<sup>1</sup>

<sup>1</sup>University of Fribourg

Version 2 ▾

Jan 29, 2021

1

Works for me

[dx.doi.org/10.17504/protocols.io.brxbm7in](https://dx.doi.org/10.17504/protocols.io.brxbm7in)

Blanchoud lab, UNIFR

Marta Wawrzyniak

University of Fribourg

## ABSTRACT

Three alternative solutions for artificial seawater (ASW) have been tested successfully on our *Botrylloides* colonies. For routine work, we use commercial sea salts (CSS), for most analyses we use the Cold Spring Harbor Protocols (CSPH) and for very clean work the K-depleted Phosphate-buffered saline (K-PBS). While, CSPH and K-PBS can be prepared at 10X, CSS should be prepared at 2X directly. Mix to dissolve and adjust pH as well as salinity to your local conditions.

## DOI

[dx.doi.org/10.17504/protocols.io.brxbm7in](https://dx.doi.org/10.17504/protocols.io.brxbm7in)

## PROTOCOL CITATION

Simon Blanchoud 2021. Artificial sea water . protocols.io

<https://dx.doi.org/10.17504/protocols.io.brxbm7in>

Version created by Marta Wawrzyniak



## KEYWORDS

artificial sea water, K-PBS, CSS

## LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jan 29, 2021

## LAST MODIFIED

Jan 29, 2021

## PROTOCOL INTEGER ID

46787

## PARENT PROTOCOLS

In steps of



[Whole colony fixation](#)

[Mitomycin C stem cell ablation](#)

## MATERIALS TEXT

- Commercial sea salts
- NaCl (CAS 7647-14-5)
- KCl (CAS 7447-40-7)
- $\text{CaCl}_2$  (CAS 10043-52-4)
- $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$  (CAS 7791-18-6)
- $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  (CAS 10034-99-8)
- $\text{NaH}_2\text{PO}_4 \cdot \text{H}_2\text{O}$  (CAS 10049-21-5)
- $\text{Na}_2\text{HPO}_4$  Sigma (CAS 7558-79-4)

### 2X CSS

1  **70 g** sea salt in  **1 L**  $\text{H}_2\text{O}$ .

2 Mix to dissolve.

3 Filter at 10  $\mu\text{m}$

### 10x CSHP

4 To prepare 10X CSHP take

A	B	C
NaCl	262.9g	4.5M
KCl	7.4g	100mM
$\text{CaCl}_2$	9.9g	90mM
$\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$	60.9g	300mM
$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$	39.4g	160mM
$\text{H}_2\text{O}$	1000ml	

### 10X K-PBS

5 To prepare 10x K-PBS take

A	B	C
NaH <sub>2</sub> PO <sub>4</sub> · H <sub>2</sub> O	8.4 g	0.06 M
Na <sub>2</sub> HPO <sub>4</sub>	34.2 g	0.24 M
NaCl	262.9 g	4.5 M
H <sub>2</sub> O	1000ml	