



Jan 15, 2021

VALAP protocol

Elizabeth Fozo¹¹In-house protocol

1

Works for me

This protocol is published without a DOI.

Eadewunm

ABSTRACT

VALAP is a biologically inert method for sealing coverslips, making it useful for live cell imaging.

PROTOCOL CITATION

Elizabeth Fozo 2021. VALAP protocol. **protocols.io**
<https://protocols.io/view/valap-protocol-brhkm34w>



KEYWORDS

VALAP, VALAP prtotoool, protocol

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 15, 2021

LAST MODIFIED

Jan 15, 2021

PROTOCOL INTEGER ID

46348

GUIDELINES

- Supplies
- Mix together
- Use:
 1. Method 1: Painting it on
 2. Method 2: Careful placement

MATERIALS TEXT

- Vaseline (also called petrolatum)
- Lanolin
- Paraffin wax (flakes work best)
- Disposable container that can be heated on a hot plate
- A disposable stirring device that can withstand heat (or glass rod that can be cleaned)

The Vaseline, lanolin, and paraffin can be purchased from Fisher/VWR or, oftentimes, at your local drug store.

SAFETY WARNINGS

- **NEVER leave unattended on the hot plate!** Valap left unattended on a hot plate will vaporize into the room and will re-solidify on other surfaces. This includes optics if there's a microscope in the room!
- Do not overheat!

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

VALAP is a biologically inert method for sealing coverslips, making it useful for live cell imaging.

BEFORE STARTING

Mix together:

- In a disposable but heat resistant container, combine equal parts by weight of the three components.
- In a fume hood, heat the container at a low temperature to slowly melt ingredients together, mixing occasionally.
- When it has combined into a homogenous liquid, carefully pour the liquid into small jars or bottles for storage. Allow it to cool before placing the cap on.

Use: Method 1 - "Painting" it on

- 1 Put the jar or bottle on a hotplate on the lowest temperature setting. **Do not leave unattended!!!**
- 2 Wait for the Valap to melt fully.
- 3 Using a cotton-tipped applicator, make sure it is well loaded with liquid Valap.
- 4 Quickly wipe the Valap along the open edge of the coverslip, half on the coverslip and half on the slide.
- 5 Repeat for the rest of the coverslip edges that need to be sealed.
- 6 Be careful not to get Valap on microscope optics! It is very difficult to remove and its presence will significantly degrade image quality.

Use: Method 2 - Careful placement

- 7 This method uses a small metal spatula and an alcohol burner.
- 8 Use the spatula to take a small amount of valap from the container.

9 Hold it over the flame to melt it.

10 Put a small bead of valap along the open edge of the coverslip.

11 Repeat for the other sides of the coverslip

Both methods take a bit of practice to do quickly. Method 2 creates a smaller profile and does not encroach on the coverslip surface as much as Method 1.