



1 ▼

Oct 21, 2021

# Agarose gel electrophoresis V.1

Yuichiroh Ikagawa<sup>1</sup><sup>1</sup>iGEM Gifu

protocol .

iGEM Gifu



Yuichiroh Ikagawa

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Electrophoresis is the most basic technique for identifying DNA amplified by PCR or restriction enzyme treatment. Agarose gels are used to separate the DNA by molecular weight. Here we describe the technique using a 1.5 % agarose gel.

Yuichiroh Ikagawa 2021. Agarose gel electrophoresis. **protocols.io**  
<https://protocols.io/view/agarose-gel-electrophoresis-bzccp2sw>



protocol ,

Oct 21, 2021

Oct 21, 2021

54372

## Reagents

Agarose S agarose powder (Nippon Gene)  
x50 TAE buffer  
DNA sample  
x6 DNA loading dye (New England Biolabs)

## Equipment

Mupid<sup>®</sup> exU electrophoresis machine (Mupid)

:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 Weigh out 0.75 g of agarose powder and add it to the triangular flask.
- 2 Add 49 ml of sterile water to the first flask and shake gently.
- 3 Add 1 ml x50 TAE buffer to the flask and heat in the microwave.
- 4 Allow the melted agarose to cool for 5 minutes before pouring it into a template set with a comb.
- 5 Cover with plastic wrap and leave for 15 minutes until the gel has solidified.

- 6 Mix 1 ul of loading dye with 5 ul of sample DNA.  
If necessary, the sample DNA solution can be diluted and added.
- 7 Carefully apply 6 ul of the DNA-dye mixture to the wells.
- 8 Run the electrophoresis 100 volts for 30 minutes.
- 9 Remove the gel from the electrophoresis bath and stain with ethidium bromide (EtBr) for 20 minutes.
- 10 Check the stained DNA with a transilluminator.