

Oct 04, 2024

## X-region PolyU/UC mutagenesis

DOI

**[dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1)**

Carolina Lopez<sup>1</sup>

<sup>1</sup>Washington University



Cecilia Escudero

wustl

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1)**

**Protocol Citation:** Carolina Lopez 2024. X-region PolyU/UC mutagenesis. **protocols.io**

**<https://dx.doi.org/10.17504/protocols.io.8epv5rxy6g1b/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 28, 2024

**Last Modified:** October 04, 2024

**Protocol Integer ID:** 100753

## Abstract

Protocol to PCR and in vitro transcribe the HCV X region w/out mutagenesis

## Materials

### **STAR METHOD**

| REAGENT                                 | SOURCE                               | IDENTIFIER    |
|---|--------------------------------------|---------------|
| Platinum® <i>Taq</i> DNA Polymerase kit | Invitrogen                           | Cat#10966-018 |
| Set of dATP, dCTP, dGTP, dTTP           | Promega                              | Cat# U1420    |
| MEGAscript® T7 Transcription Kit        | Ambion<br>Thermo Fisher Technologies | Cat#1334      |
| HCV polyU/UC plasmid                    | Gift from Michael Gale lab           | NA            |
| HCV X-region plasmid                    | Gift from Michael Gale lab           | NA            |
| UltraPure Water                         | Invitrogen                           | Cat#10977-015 |

### **Material:**

Primer: X-region F/ X-region R ( 10 micromolar ( $\mu$ M) )

PolyU/UC F and R ( 10 micromolar ( $\mu$ M) )

dNTPs ( 10 millimolar (mM) ) Fermentas –For PCR

HCV PolyU/UC X-region plasmid



## HCV X-region PolyU/UC PCR - T7 in vitro transcription

### 1 PCR:

- 1) Thaw component and cDNA, spin down before use
- 2) Assemble reaction: Master Mix: 10X (each reaction assemble 3 tubes)

| Reagent               | Volume (ul) |
|-----------------------|-------------|
| H <sub>2</sub> O      | 26          |
| 10X Buffer            | 5           |
| MgCl <sub>2</sub>     | 5           |
| dNTP (10uM)           | 2           |
| Forward Primer (10uM) | 5           |
| Reverse Primer (10uM) | 5           |
| Taq pol. (5unit/ul)   | 0.25        |
| Plasmid (10ng/ul)     | 1.5         |

Total:  50 µL

### 3) Run PCR

| Program         | Temperature (C) | Time   |           |
|-----------------|-----------------|--------|-----------|
| Denaturation    | 95              | 3 min  | Hold      |
| Denaturation    | 95              | 30 sec | 35 cycles |
| Annealing       | 55              | 30 sec |           |
| Extension       | 72              | 90 sec |           |
| Final extension | 72              | 5 min  | Hold      |
|                 | 4               | ∞      | Hold      |

## II. In vitro transcription:

**15m**

- 2
  - > Spin-down all reagents
  - > Keep the 10X Reaction Buffer at room temperature while assembling the reaction.
  - > Vortex all solutions until they are completely in solution.
  - > Assemble dNTP mix (2µL ATP+ 2µL CTP+ 2µL GTP+ 2µL TTP/reaction).
  - > Store the ribonucleotides on ice!

**15m**



Assemble the reaction in an appropriate RNase free PCR tube.

- 1) 1  $\mu\text{g}$  DNA plasmid (digested)
- 2) 8  $\mu\text{L}$  NTP Mix
- 3) 2  $\mu\text{L}$  10X Reaction Buffer
- 4) 0.5-1  $\mu\text{L}$  RNaseOUT (RNase Inhibitor)
- 5) 2  $\mu\text{L}$  T7 enzyme mix
- 6) RNase Free water up to 20  $\mu\text{L}$

- Pipette the mixture up and down gently and microfuge
- Incubate 3-6hrs at 37 °C using PCR machine. The shorter the ivtRNA, the longer the reaction time needed.
- (Optional but do most times) Template plasmid Digestion: 1  $\mu\text{L}$  Turbo DNase and mix well,

00:15:00 at 37 °C

#### IV. RNA precipitation: Transfer to an 1.5ml RNase free vial

35m

- 3 If you expect very low yields of RNA, do not dilute the transcription reaction with water prior to adding the LiCl.

35m

- Add 30  $\mu\text{L}$  H<sub>2</sub>O  
 30  $\mu\text{L}$  LiCl
- Mix well, incubate 00:30:00 at -20 °C
- Centrifuge at 4°C/20 min/max
- Remove supernatant
- Add 1 mL 70% COLD ethanol
- Centrifuge at 4°C/5 min/max
- Remove supernatant and allow the pellet to dry
- Resuspend in 15-20  $\mu\text{L}$  RNase-free H<sub>2</sub>O
- Place the tube 00:05:00 at 65 °C water bath
- Dilute 1  $\mu\text{L}$  in 19  $\mu\text{L}$  H<sub>2</sub>O to measure concentration

#### Protocol references

Jie Xu (Lopez Lab)