



Aug 14, 2020

Universal Sandwich ELISA for investigating the binding of Protein-LAG (SpLAG) to avian immunoglobulins using anti-IgY-peroxidase as conjugate.

Angel A Justiz-Vaillant¹¹University of the West Indies St. Augustine**1** Works for me dx.doi.org/10.17504/protocols.io.bjq5kmy6[University of the West Indies](#) angel.vaillant@sta.uwi.eduAngel Justiz-Vaillant
University of the West Indies St. Augustine

DOI

dx.doi.org/10.17504/protocols.io.bjq5kmy6

PROTOCOL CITATION

Angel A Justiz-Vaillant 2020. Universal Sandwich ELISA for investigating the binding of Protein-LAG (SpLAG) to avian immunoglobulins using anti-IgY-peroxidase as conjugate.. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bjq5kmy6>

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 14, 2020

LAST MODIFIED

Aug 14, 2020

PROTOCOL INTEGER ID

40445

MATERIALS

| NAME | CATALOG # | VENDOR |
|--|-----------|---------------|
| Anti-Chicken IgY, HRP Conjugate, 300ul | G1351 | Promega |
| Nunc™ 96-Well Polystyrene Round Bottom Microwell Plates, V 96 well plate, Non-Treated, clear, without lid, Sterile | 260210 | Thermo Fisher |
| Staphylococcal Protein-A | | Sigma Aldrich |
| Protein-L from P. Magnus | | |
| Streptococcal protein G by Sigma Aldrich | | |

- 1 This ELISA is used to study the interaction of protein-LAG (SpLAG) with diverse avian immunoglobulins.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of a mixture of SpL, SpA and SpG in carbonate-bicarbonate buffer pH 9.6.

- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of avian egg yolk or egg white (1 mg/ml) is added and incubated for 1.30h at RT and the microplate is then rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled-anti-IgY conjugate diluted 1:15000 in PBS-non-fat milk is added to each well and incubated for 1.30h at RT. After that the plate is washed 4X with PBS-Tween.
- 6 Pipette 50 µl of 3,3',5,5' - tetramethylbenzidine (TMB; Sigma-Aldrich) to each well.
- 7 The reaction is stopped with 50 µl of 3M H₂SO₄ solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 450 nm.
- 9 A cut-off point can be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher will be the binding affinity of SpLAG to avian immunoglobulins.