



Version 5

Sep 29, 2021

# Genomic DNA extraction from the diatom *Pseudo-nitzschia multistriata* for Illumina sequencing V.5

Francesco Manfellotto<sup>1</sup>, Monia Teresa Russo<sup>2</sup>, Pina Marotta<sup>1</sup>, Rossella Annunziata<sup>1</sup>, Anna Santin<sup>1</sup>, Antonella Ruggiero<sup>1</sup>, Mariella Ferrante<sup>1</sup>

<sup>1</sup>Stazione Zoologica Anton Dohrn; <sup>2</sup>SZN

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.byk7puzn](https://dx.doi.org/10.17504/protocols.io.byk7puzn)

Francesco Manfellotto  
Stazione Zoologica Anton Dohrn

## ABSTRACT

Genomic DNA extraction from the diatom *Pseudo-nitzschia multistriata* for Illumina sequencing

## DOI

[dx.doi.org/10.17504/protocols.io.byk7puzn](https://dx.doi.org/10.17504/protocols.io.byk7puzn)

## PROTOCOL CITATION

Francesco Manfellotto, Monia Teresa Russo, Pina Marotta, Rossella Annunziata, Anna Santin, Antonella Ruggiero, Mariella Ferrante 2021. Genomic DNA extraction from the diatom *Pseudo-nitzschia multistriata* for Illumina sequencing. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.byk7puzn>  
Version created by Francesco Manfellotto

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 28, 2021

## LAST MODIFIED

Sep 29, 2021

## PROTOCOL INTEGER ID

53631

- 1 Grow cells as described in: [dx.doi.org/10.17504/protocols.io.bgudjws6](https://dx.doi.org/10.17504/protocols.io.bgudjws6)
- 2 Detect the presence/absence of bacteria and collect the cultures as described in: [dx.doi.org/10.17504/protocols.io.btt5nng6](https://dx.doi.org/10.17504/protocols.io.btt5nng6)
- 3 Resuspend cells with 500 µL of TE buffer (10 mM TrisHCl pH 7.6 and 1 mM EDTA pH 8.0)

- 4 Add:
  - 400 mg of 0.2-0.3 mm diameter zirconia/silica beads
  - 500 µL phenol (pH 7.8).
- 5 Mix with vortex 3 times at 30 Hz for 85 seconds (Each time put sample in ice for 60 seconds before vortex)
- 6 Centrifuge at 11000 g for 5 minutes at 4°C. 5m
- 7 Recover aqueous phase in new 1.5 mL Eppendorf tubes (about 600 µl)
- 8 Add 500 µL of P.C.I (Phenol:Chloroform:Isoamyl alcohol 25:24:1 v/v) and mix by inversion
- 9 Centrifuge at 11000 g for 5 minutes at 4°C 5m
- 10 Move the aqueous phase in a new Eppendorf tube and add 5 µL of RNase-A 10 mg/mL
- 11 Incubate at 37 °C for 30 minutes. 30m
- 12 Add 500 µl di P.C.I (Phenol:Chloroform:Isoamyl alcohol 25:24:1 v/v) and mix by inversion.
- 13 Centrifuge at 11000 g for 5 minutes at 4°C. 5m
- 14 Move the aqueous phase in a new 2 mL Eppendorf tube and add:
  - 50 µL of 3 M NaAc (pH ± 5)
  - 1 mL of ethanol 96% (- 20 °C)
  - 2 µL glycogen (- 20 °C)
- 15 Incubate over night at -20°C. 12h
- 16 Centrifuge the samples at 13000 g for 30 minutes at 4°C 30m

- 17 Discard aqueous phase
- 18 Add 1 mL ethanol 70% and mix gently by inversion
- 19 Centrifuge at 13000 g for 10 minutes at 4°C 10m
- 20 Discard aqueous phase
- 21 Add 1 mL ethanol 70% and mix gently by inversion
- 22 Centrifuge at 13000 g for 10 minutes at 4°C 10m
- 23 Remove aqueous phase and dry pellet at R.T. for at least 20 minutes 20m
- 24 Add 50 µL of preheated (55°C) TE 1X (pH 8) or sterile MilliQ water
- 25 Incubate at 55 °C for 20 minutes 20m
- 26 Quantify DNA concentration by Nanodrop or Qubit
- 27 Check DNA integrity.  
Run a small amount of DNA with 1% agarose gel
- 28 Store DNA at +4 °C, or -20°C for longer storage times