

Apr 04, 2025

🌐 Protocol for 6mA labeling and HMW DNA extraction from fresh frozen human brain samples

DOI

dx.doi.org/10.17504/protocols.io.j8nlkw54wl5r/v1

Jonas Demeulemeester¹, Koen Theunis²

¹Laboratory of Integrative Cancer Genomics, VIB-KULeuven;

²Laboratory of Computational Biology, VIB-KULeuven

Lab of Computational Bi...

ASAP Collaborative Res...

1 more workspace



Koen Theunis

Laboratory of Computational Biology, VIB Center for AI & Com...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.j8nlkw54wl5r/v1

Protocol Citation: Jonas Demeulemeester , Koen Theunis 2025. Protocol for 6mA labeling and HMW DNA extraction from fresh frozen human brain samples. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkw54wl5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 09, 2022

Last Modified: April 04, 2025

Protocol Integer ID: 69773

Keywords: HMW DNA extraction, Nanopore sequencing, PromethION, ASAPCRN



Abstract

This protocol details the procedure of 6mA labeling and HMW DNA extraction of fresh frozen brain tissue. The protocol is inspired by Fiber-seq.

CITATION

Stergachis AB, Debo BM, Haugen E, Churchman LS, Stamatoyannopoulos JA (2020). Single-molecule regulatory architectures captured by chromatin fiber sequencing.. Science (New York, N.Y.).

LINK

<https://doi.org/10.1126/science.aaz1646>

Attachments




[irgkbewa7.pdf](#)



192KB

Materials

Prepare buffers (Volumes above are indicated per sample):

10x stocks of Wash buffer base (no spermidine/Tween) and Labeling buffer base (no spermidine/SAM/Hia5).

2 mL Wash buffer → take  488 µL and add Digitonin + 50XProtInh to make **Lysis buffer**.

 50 µL **labeling buffer** / 250.000 cells (typically ~ 1E6 cells →  200 µL).

Make  25 millimolar (mM) **spermidine stocks fresh monthly** and store at  -20 °C .

Add  2 µL spermidine (liquid, ~  6.38 Molarity (M)) to  350 µL  0.1 Molarity (M) HCl and

 160 µL H2O – **check pH with strips!**

| A | B | C | D |
|---|---------|---------------|--------|
| 10X Wash/Lysis base (10X-WLB) (200uL / sample) | | | |
| | Stock | Final for 10X | V (uL) |
| Tris-HCl pH 7.4 | 1000 mM | 100 mM | 100 |
| NaCl | 5000 mM | 100 mM | 20 |
| Water | | | 880 |
| Total | | | 1000 |

| A | B | C | D |
|---|---------|---------------|--------|
| 10X Hia5 labeling base (10X-H5B) (20uL / sample) | | | |
| | Stock | Final for 10X | V (uL) |
| Tris-HCl pH 8 | 1000 mM | 150 mM | 150 |
| NaCl | 5000 mM | 150 mM | 30 |
| KCl | 1000 mM | 600 mM | 600 |
| EDTA pH 8 | 500 mM | 10 mM | 20 |
| EDTA pH 8 | 250 mM | 5 mM | 20 |
| Water | | | 180 |
| Total | | | 1000 |

| A | B | C | D | E |
|-----------------------|---|---|---|---|
| 1X Wash buffer | | | | |

| A | B | C | D | E |
|-------------------|-------|--------|--------|-------|
| | Stock | Final | V (uL) | X5 |
| 10X WLB | 10 X | 1 X | 230 | 1150 |
| BSA | 10% | 0.10% | 23 | 115 |
| Spermidine pH 7.4 | 25 mM | 0.5 mM | 46 | 230 |
| Tween-20 | 10% | 0.10% | 23 | 115 |
| Water | | | 1978 | 9890 |
| Total | | | 2300 | 11500 |

| A | B | C | D | E |
|------------------------|-------|-------|--------|------|
| 1X Lysis buffer | | | | |
| | Stock | Final | V (uL) | X5 |
| Wash buffer | | | 488 | 2440 |
| Digitonin | 5% | 0.02% | 2 | 10 |
| Proteaselnh | 50 X | 1 X | 10 | 50 |
| Total | | | 500 | 2500 |

| A | B | C | D | E |
|--------------------------------|-------|-------|--------|------|
| 1X Hia5 labeling buffer | | | | |
| | Stock | Final | V (uL) | X5 |
| 10X-H5B | 10 X | 1 X | 20 | 100 |
| BSA | 10 | 0.1 | 2 | 10 |
| Spermidine pH 7.4 | 25 | 0.5 | 4 | 20 |
| SAM | 32 | 0.8 | 5 | 25 |
| Hia5 enzyme | 250 | 5 | 4 | 20 |
| Water | | | 165 | 825 |
| Total | | | 200 | 1000 |



6mA labeling and HMW DNA extraction

- 1 Place the Dounce homogenizer and pestles On ice , chill the centrifuge to 4 °C and preheat the ThermoMixer to 37 °C .
- 2 **Carefully** transfer 3-4 (2 mm diameter) tissue punch biopsies (~ 25 mg) to the chilled Dounce homogenizer.

Note

Keep the Dounce homogenizer on ice during the entire disruption process.

- 3 Add 500 µL of the 1X Lysis buffer and let the tissue thaw for 00:01:00 .

1m



- 4 Gently homogenize the tissue 10X with pestle A and 10X with B.
 - 4.1 Push the tissue with the pestle firmly into the bottom of the Dounce chamber with each stroke (Down + Up = 1X).
 - 4.2 Keep the tissue between tip of pestle and the bottom of the Dounce chamber for thorough homogenization.
 - 4.3 Homogenate may become foamy, but this is not a cause for concern.

Note

In the next step, transfer any foam that forms.

- 5 Incubate On ice for 00:05:00 before adding 1000 µL of 1X Wash buffer.


5m



- 6 Transfer the lysate to a 2 mL Protein LoBind microcentrifuge tube.







7 Rinse the pestles and homogenizer with the remaining  500 μ L 1X Wash buffer and add to the sample.




8 Pellet homogenate by centrifuging at  700 x g and  4 °C for  00:05:00 . Discard supernatant.

5m



9 Resuspend the pellet in  200 μ L of 1X Hia5 labeling buffer – use a  1 mL or wide bore tip.

10 Incubate on a ThermoMixer at  37 °C and  900 rpm for  00:30:00 .

30m





- 11
- **Continue with Circulomics CBB Tissue protocol from step 8 onwards.**
 - **Continue with NEB Monarch HMW.**





Nanopore sequencing (LSK-110, PromethION)

2h 5m



12 According to ONT protocol Genomic DNA by Ligation (SQK-LSK110) with the following modifications:

12.1 Start with  3 μ g -  4 μ g of HMW DNA in 150 uL and shear 25x with a 26G needle or in Megaruptor to 35kb.



12.2 Adjust volumes end-prep and FFPE repair step accordingly (i.e. vol x 3), omit control strand (CS).

12.3 Extend end-prep and FFPE repair steps from  00:05:00 to  00:30:00 (i.e. 30min at  20 °C and 30 min at  65 °C).

35m


12.4 Extend ligation step to  01:00:00 at  Room temperature .

1h

12.5 Elute the AMPure cleanups for  00:10:00 and  00:20:00 after the end-prep and ligation steps.

30m



- 12.6 This should yield a ~3x library, aim to load near the high end of the 5-50fmol range, typically ~  8 μ L .

Note

Note 1: Library prep yield is typically 30-50%.

Note 2: The amount loaded can be reduced during subsequent flushes to balance seq yield with # flushes.

Citations

Stergachis AB, Debo BM, Haugen E, Churchman LS, Stamatoyannopoulos JA. Single-molecule regulatory architectures captured by chromatin fiber sequencing.

<https://doi.org/10.1126/science.aaz1646>