



Oct 11, 2022

# 🌐 Preparation of bacterial cell lysate for SDS-PAGE (2022 iGEM)

Forked from [Preparation of bacterial cell lysate for proteomics \(LC-MS\) by freeze and thaw cycles](#)

Team Fudan iGEM<sup>1</sup>

<sup>1</sup>Fudan University

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.14egn2zwyg5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn2zwyg5d/v1)

2022



Team Fudan iGEM  
Fudan University

## ABSTRACT

Preparation of bacterial cell lysate for SDS-PAGE

## DOI

[dx.doi.org/10.17504/protocols.io.14egn2zwyg5d/v1](https://dx.doi.org/10.17504/protocols.io.14egn2zwyg5d/v1)

## PROTOCOL CITATION

Team Fudan iGEM 2022. Preparation of bacterial cell lysate for SDS-PAGE (2022 iGEM). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.14egn2zwyg5d/v1>



## FORK NOTE

## FORK FROM

Forked from Preparation of bacterial cell lysate for proteomics (LC-MS) by freeze and thaw cycles, alexrrus


## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 10, 2022

**Preparation of bacterial cell lysate for SDS-PAGE**

- 1 Grow the cells overnight, to an OD<sub>600</sub> above 1. If IPTG induction needed, at it around OD<sub>600</sub> 0.2-0.3 (keep IPTG stock solution at -20 freezer).
- 2  16h  
Take one 1 ml of OD<sub>600</sub> = 1 from each sample. Correct the sample volume to obtain an equivalent size pellet. For us, as quantified using NanoCym950 nanoparticles (equivalent to the size of *E. coli*), 1 OD<sub>600</sub> equals to 10<sup>8</sup> nanoparticles per mL. We estimate each well on SDS-PAGE contains proteins from 2\*10<sup>6</sup> bacteria.
- 3 Spin down the bacterial cells for 1 min at 13,000 rpm in a microfuge. 10m
- 4 Aspirate any trace of supernatant with a vacuum line avoiding touching the pellet, but removing all liquid. 1m
- 5 To each cell pellet, add 250 µL of distilled water, resuspend the pellet. Add 250 µL 2x SDS sample buffer (see below). Heat with a 95 degree heating block for 10 minutes. Store the samples at -20 degree freezer until SDS-PAGE. 1m

**Receipt**

- 6 2x SDS sample buffer, works for Tris-Cl or MOPS based gels
  - Glycerol 20 ml
  - SDS 4 g
  - 0.5 M Tris-HCl pH 6.8 stock solution 25 ml
  - Bromo Phenol Blue 1 mg; dissolve and volume to 95 ml
  - add 5% β-Mercaptoethanol before usage