



Aug 19, 2021

Endoglycosidase H (Endo H) digestion assay

Vik Khurana¹¹Brigham Womens Hospital

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bu3qnymw

hendersa

ABSTRACT

This protocol details the endoglycosidase H (Endo H) digestion assay.

ATTACHMENTS

[dh4abiq7.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bu3qnymw

PROTOCOL CITATION

Vik Khurana 2021. Endoglycosidase H (Endo H) digestion assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bu3qnymw>

KEYWORDS

Endoglycosidase H, Endo H, Endoglycosidase H digestion assay

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 17, 2021

LAST MODIFIED

Aug 19, 2021

OWNERSHIP HISTORY

May 17, 2021 Urmilas

Jun 23, 2021 hendersa

PROTOCOL INTEGER ID

50000

MATERIALS TEXT

Chemicals and supply:

[Protease Inhibitor Cocktail Sigma](#)

- 100X [Aldrich Catalog #P8340](#) ; "PiM" aliquots at -20°C

[Phosphatase Inhibitor Cocktail 1 Sigma](#)

- 100X [Aldrich Catalog #P2850](#) ; "PP" aliquots at **4 °C** , "Antibodies,

[Phosphatase Inhibitor Cocktail 2](#) **Sigma**

large vials" box); 100X **Aldrich Catalog #P5726**

; at **4 °C**,

"Antibodies, large vials" box)

- **Thickwall Polycarbonate Tubes (Beckman Coulter, 343775)**
- **10X Glycoprotein Denaturing Buffer (NEB, re-order B1704S)**
- **10X G5 Reaction Buffer (NEB, re-order B1702S)**

[Endo H - 50,000 units](#) **New England**

- **Biolabs Catalog #P0702L**
- **Novex® 10% Tris-Glycine Mini Gels, 1.5 mm, 10 wells (Life Technologies, EC6078BOX)**
- **Novex® 10% Tris-Glycine Midi Gels, 1 mm, 12+2 wells (Life Technologies, WT0101BOX)**
- **Novex® 10% Tris-Glycine Midi Gels, 1 mm, 20 wells (Life Technologies, WT0102BOX)**

[Tris-Glycine-SDS Running Buffer \(10X\)](#) **Boston**

- **Bioproducts Catalog #BP-150**

[Bio Rad Precision Plus protein Standard](#) **Bio-rad**

- **Laboratories Catalog #161-0374**
- 10X TRIS-Glycine Transfer Buffer (Boston BioProducts)

[Nicastrin Antibody](#) **Cell Signaling**

- **Technology Catalog #3632** ; **-20 °C**

[Anti-Glucocerebrosidase \(C-terminal\) antibody produced in rabbit](#) **Sigma**

- **Aldrich Catalog #G4171** ; **4 °C**,
"Antibodies, large vials" box)

[Anti-Glyceraldehyde-3-Phosphate Dehydrogenase Antibody clone 6C5](#) **Emd**

- **Millipore Catalog #MAB374**

4 °C

- **Gt anti-rb IgG, DyLight™ 800 conjugated (Rockland, 611-145-003; Lauren's fridge)**
- **Gt anti-ms IgG, DyLight™ 680 conjugated (Rockland, 610-144-003; Lauren's fridge)**

Buffers:

Table 1: Mazzulli lysis buffer:

Preparation of the lysis buffer: first add glycerol, NaCl and HEPES into a measuring cylinder, top up with diH₂O to

195 mL in total. Add EGTA, MgCl₂ and finally Triton X-100. Be very careful when adding the viscous Triton X-100 to

make sure that exactly **2 mL** go into the solution. Stir for ~ **00:20:00**, adjust the **pH7.4** and top up with diH₂O to

200 mL in total using a measuring cylinder. Filter with a SteriCup to ensure long-term storage – in this way the lysis

buffer can be kept at **4 °C** for up to 6 months or aliquoted at **-20 °C** (**40 mL** aliquots).

A	B	C	D
Ingredient	200 ml	Stock concentration	Final concentration
HEPES	4 ml	1 M	20 mM
NaCl	6 ml	5 M	150 mM
Glycerol	40 ml	50% (v/v)	10% (v/v)
EGTA	200 µl	1 M	1 mM
MgCl ₂	300 µl	1 M	1.5 mM
Triton X-100	2 ml		1% (v/v)
diH ₂ O	147.5 ml		

Recipe according to Mazzulli et al. (2006), PMID: 17005870.

Table 2: 5X Sample loading buffer (SB). Aliquot and store at **-20 °C**.

A	B	C	D
Ingredient	10 ml	Stock concentration	Final concentration
Tris-Cl, pH 6.8	1 ml	1 M	100 mM
Glycerol	4 ml	50% (v/v)	20% (v/v)
β-ME	630 μl		900 mM
Bromophenol blue	1.5 mg		0.015% (w/v)
SDS*	1 g		10% (w/v)
diH ₂ O	~ 4 ml		

*Work under the fume hood when adding SDS! Powder is toxic! Recipe according to Mazzulli et al. (2006), PMID: 17005870.

Endoglycosidase H (Endo H) digestion assay

12h 28m

1 PREPARE:

1.1 Make sure the lysis buffer (Table 1) is chilled.

Table 1 | Mazzulli lysis buffer. Preparation of the lysis buffer: first add glycerol, NaCl and HEPES into a measuring cylinder, top up with diH₂O to 195 ml in total. Add EGTA, MgCl₂ and finally Triton X-100. Be very careful when adding the viscous Triton X-100 to make sure that exactly 2 ml go into the solution. Stir for ~ 20 min, adjust the pH to 7.4 and top up with diH₂O to 200 ml in total using a measuring cylinder. Filter with a SteriCup to ensure long-term storage – in this way the lysis buffer can be kept at 4°C for up to 6 months or aliquoted at -20°C (40 ml aliquots).

Ingredient	200 ml	Stock concentration	Final concentration
HEPES	4 ml	1 M	20 mM
NaCl	6 ml	5 M	150 mM
Glycerol	40 ml	50% (v/v)	10% (v/v)
EGTA	200 μl	1 M	1 mM
MgCl ₂	300 μl	1 M	1.5 mM
Triton X-100	2 ml		1% (v/v)
diH ₂ O	147.5 ml		

Recipe according to Mazzulli et al. (2006), PMID: 17005870.

1.2



30m

Cool down ultracentrifuge to **4 °C** . Set parameters: TLA-100 rotor, **100000 x g** , **00:30:00** , **4 °C** ; apply vacuum.

1.3

Prepare lysis buffer MM: Add 100X protease inhibitor cocktail ("PiM") as well as 100X phosphatase inhibitor cocktail I ("PP") and II to the lysis buffer in a 1:100 ratio. Keep the lysis buffer MM **On ice** from this point onwards. The protease inhibitor cocktail can be briefly thawed using the **37 °C** heat block. Keep at **Room temperature** before re-freezing.

1.4 Get a large box and fill to approximately 2/3 with wet ice for the ice/water slurry (see below).

2 Get samples and keep on dry ice.

3

When ready, thaw samples on wet ice and add an appropriate volume of the lysis buffer MM as soon as the cell pellet is almost thawed.

Can also be thawed briefly at **Room temperature** but make sure that the cell pellet is not completely thawed when lysis buffer MM is added.

Lysis buffer MM volume: approximately **50 µl** - **150 µl** ; around 3x the volume of the cell pellet; better too little than too much.

4 Resuspend the cell pellet to a homogenous (!) solution by trituration and keep on wet ice.

Make sure the top of the Eppis is carefully labeled since side labels will likely disappear during a later incubation in the dry ice/EtOH slurry!.

5 Prepare the ice/water slurry: add cold water to approximately half of the height of the ice layer.

6

20m

Incubate samples for **00:20:00** in the ice/water slurry.

7 In the meantime PREPARE the following:

7.1 Label ultracentrifuge tubes on the side. Use P200 tip holder as a rack for the ultracentrifuge tubes.

7.2 Prepare dry ice/EtOH slurry approximately ⌚ 00:10:00 before the ⌚ 00:20:00 incubation on the ice/water slurry is completed: 30m

- Fill small box to 2/3 with EtOH (190 proof, 95%).
- Add dry ice until saturation (dry ice doesn't resolve quickly anymore; temperature of ⚡ -80 °C has stabilized).

8 

1m

Put Eppis into a floating rack and immediately place them in the dry ice/EtOH slurry in order to snap freeze the samples; incubate for ⌚ 00:01:00 .

Make sure that the lid does not touch the dry ice/EtOH slurry to prevent EtOH leaking into the Eppis; only cell lysate needs to be exposed to the dry ice/EtOH slurry. Incubation time of ⌚ 00:01:00 is minimum; can incubate longer.

9 

1m

Directly transfer Eppis onto a heat block and incubate at ⚡ 37 °C for exactly (!) ⌚ 00:01:00 (no longer!).

10 Repeat the described freeze/thaw cycle once → two freeze/thaw cycles (⚡ -80 °C / ⚡ 37 °C) in total.





11 Keep samples on wet ice afterwards.

12 Get TLA-100 rotor (Beckman Coulter) and keep on wet ice, including the lid. TLA-100 rotor usually kept in the ⚡ 4 °C cold room.

13 Transfer the cell lysate from the Eppis into the labeled ultracentrifuge tubes (thickwall polycarbonate tubes) and place them in the TLA-100 rotor ⚡ On ice ; keep the empty Eppis for a later step in case the labeling is sufficient and hasn't been removed by the EtOH.



14

Place TLA-100 rotor into the ultracentrifuge; make sure it's correctly placed (!); re-apply vacuum, make sure that the parameters are set correctly ( **100000 x g** for  **00:30:00** at  **4 °C**) and start. The ultracentrifuge might take a while to start because it will have to re-establish the vacuum; make sure that the ultracentrifuge gets back to  **4 °C** .

15 During the ultracentrifugation PREPARE the following:

15.1 Prepare a plate for the BCA test; each condition in triplicate.



15.2 Label Eppis and chill them on wet ice.

15.3 Get dry ice.







16 Get TLA-100 rotor and keep on wet ice.

17 Remove ultracentrifuge tubes and transfer the supernatant to a new, correctly labeled and chilled Eppi; keep Eppis on wet ice.

18 A small cell pellet should be visible; transfer to the appropriate old, correctly labeled Eppi and keep on dry ice.

This cell pellet can be re-used, e.g., to check for α Syn levels by Western blotting; in this case store cell pellet at  **-80 °C** and later resuspend in approximately  **10 μ l** of the appropriate lysis buffer for Western blotting.

19

For the BCA test either add  **1 μ l** of the cell lysate to  **24 μ l** water or dilute  **2 μ l** of the cell lysate in a 1:1 ratio with  **2 μ l** water before adding  **1 μ l** of the 1:2 dilution to  **24 μ l** water; in the latter case don't forget to change the dilution factor on the Endo H spreadsheet to 50!.

I got more accurate results with a dilution of 1:25.

20 

If the Endo H is proceeded with on another day, transfer the Eppis with the cell lysate directly onto dry ice and immediately store them at δ **-80 °C** (to reuse, thaw cell lysates on wet ice); if the Endo H is to be proceeded with on the same day, place Eppis on wet ice and store intermediately in the δ **4 °C** cold room.

21 

26m

During the incubation for the BCA test at δ **37 °C** for \odot **00:26:00** roughly label 2 Eppis for each condition (\rightarrow +/- Endo H).

22 Thaw 10X glycoprotein denaturing buffer and 10X G5 buffer at δ **Room temperature** (don't vortex denaturing buffer as this will create a lot of small bubbles!).

Fine to briefly keep them in the heat block if they are not thawed when needed in the steps below; make sure that these buffers are very well mixed before adding them to the reaction volume!.

23 After BCA test use Chee-Yeun's Endo H spreadsheet to calculate the protein concentration.

24 PREPARE: thaw Endo H (stored at δ **-20 °C**) on wet ice.

25 From this point onwards all in duplicate (\rightarrow +/- Endo H):

25.1 

Refer to the spreadsheet for mixing an appropriate volume of the cell lysate (\square **30 μ g** - \square **40 μ g** protein) with diH₂O for a total volume of \square **15.3 μ l** ; from this point onwards samples can be kept at δ **Room temperature** .




25.2 

Add \square **1.7 μ l** of the glycoprotein denaturing buffer for a total reaction volume of \square **17 μ l** .



25.3 

Briefly vortex and spin down.

25.4 Boil at  **100 °C** for  **00:10:00** . In the meantime transfer the cell lysate stock to  **-80 °C** . ^{10m}

Make sure that the lids of the Eppis are tightly closed!

25.5 In the meantime PREPARE the +/- Endo MMs:

MM-1 (+ Endo H): ( **1 µl Endo H** +  **2 µl** 10X G5 buffer)/rxn; can be kept at

 **Room temperature** .

MM-2 (- Endo H): ( **1 µl diH₂O** +  **2 µl** 10X G5 buffer)/rxn.


25.6  

Spin down samples briefly; add  **3 µl** of MM-1 or -2, briefly vortex and spin down.

25.7 

2h

Incubate at  **37 °C** for  **02:00:00** (minimum incubation time).

25.8 

Spin down briefly.

25.9 

Store samples at  **-80 °C** or  **-20 °C** if the gel is to be run on another day.

26 Thaw sample buffer ("5X SB", kept at -20°C; Table 2) and protein ladder at  **Room temperature** .

27 

In case +/- Endo H samples have been frozen down, thaw at  **Room temperature** and briefly spin down.

28  

10m

Add **5 µl** 5X SB to each **20 µl** rxn; briefly vortex and spin down; denature proteins at **100 °C** for **00:10:00** ; briefly spin down.

- 29 PREPARE Novex® 10% Tris-glycine gel with adapter. For 1 gel prepare **500 mL** running buffer: **50 mL** 10X Tris-glycine-SDS running buffer (Boston BioProducts) + **450 mL** diH₂O.

No need to recycle running buffer.

- 30 Load **8 µl** of the protein ladder and the entire reaction volume of the +/- Endo H samples – first all -Endo H samples followed by all +Endo H samples.

- 31 Run the gel:

31.1 Run at 100 V first for **00:20:00** – until the running front is completely condensed by means of the ^{20m} stacking gel.

31.2 Then increase voltage to 120 V and set timer to **02:00:00** . ^{3h}
▪ Run the gel until the 37 kDa protein ladder band is approximately **1 cm** - **1.5 cm** above the lower end of the gel cassette; check gel regularly after **01:00:00** ; don't run too far!!.

- 32 Wet electrophoretic transfer onto a PVDF membrane (Bio-Rad) with 10X Tris-glycine transfer buffer (Boston BioProducts): **700 mL** diH₂O + **100 mL** 10X Tris-glycine transfer buffer + **200 mL** MeOH [20% (v/v)].

32.1 Prepare **2 L** transfer buffer and keep in **4 °C** cold room; can be recycled 3-5x.

32.2 Transfer for **02:00:00** at 60 V and **4 °C** or for **02:30:00** at 65 V in the case of two ^{4h 30m} simultaneous transfers.

- 33 Continue with immunoblotting using the Odyssey infrared imaging system (LI-COR Biosciences).

- 34 Antibodies:

A	B	C	D
Primary antibody		Secondary antibody	
Rb anti-GCase	1:500	Gt anti-rb IgG, DyL 800	1:5,000
Rb anti-nicastrin	1:500		
Ms anti-GAPDH	1:3,000	Gt anti-ms IgG, DyL 680	1:10,000