



Oct 05, 2022

🌐 Sample preparation and lysis of homogenized malaise trap samples

Dominik Buchner¹¹University of Duisburg-Essen, Aquatic Ecosystem Research

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1

Dominik Buchner

University of Duisburg-Essen, Aquatic Ecosystem Research

ABSTRACT

This protocol describes the steps of sample preparation and lysis before DNA extraction for the Malaise trap metabarcoding protocol of the LeeseLab.

DOI

dx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1

PROTOCOL CITATION

Dominik Buchner 2022. Sample preparation and lysis of homogenized malaise trap samples. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.dm6gpjrmjgzp/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 29, 2022

LAST MODIFIED

Oct 05, 2022

PROTOCOL INTEGER ID

70661

PARENT PROTOCOLS

In steps of

[Guanidine based DNA extraction with silica-coated beads or silica spin-columns](#)

GUIDELINES

Follow general lab etiquette. Wear gloves to prevent contaminating the samples.
Clean the workspace before starting with 80% EtOH.

MATERIALS TEXT

Materials required:

Below all materials needed for the protocol are listed. Vendors and part numbers are listed but interchangeable depending on the supply situation.

Chemicals:

[Tris ultrapure](#)

Tris ultrapure 99.9% [99.9% Diagonal Catalog #A1086.1000](#)

[Sodium chloride](#) Fisher

Sodium chloride [Scientific Catalog #10112640](#)

[EDTA disodium salt](#) Sigma

EDTA disodium salt [Aldrich Catalog #E5134-50G](#)

[Sodium dodecyl](#)

SDS ultrapure [sulfate Diagonal Catalog #A1112.0500](#)

[Proteinase](#)

Proteinase K [K 7BioScience Catalog #RP100B](#)

Hydrochloric acid fuming 37%

[Hydrochloric acid fuming 37%](#) Sigma

[Aldrich Catalog #1003171011](#)

[Sodium hydroxide - pellets](#) Fisher

Sodium hydroxide [Scientific Catalog #S/4920/60](#)

[Calcium chloride 94%](#) Carl

Calcium chloride [Roth Catalog #A119.1](#)

Glycerol 87%

[Glycerol 87 % for molecular biology](#) Panreac

[AppliChem Catalog #A3739,1000](#)

Labware:

2 mm dia zirconia beads

[Zirconia Beads 2 mm dia](#) BioSpec

[Products Catalog #11079124zx](#)

 2 mL screwcap

2 mL screwcap tubes  [tube Sarstedt Catalog #72.693](#)

Wide-bore tips

 [ART Wide Bore Tip 1000 uL Thermo Fisher](#)

Scientific Catalog #2079G






Disposable PES Bottle Top Filters

 [Fisherbrand Disposable PES Bottle Top Filters Fisher](#)




Scientific Catalog #15973307

Stock solutions:





1 L Tris stock solution 1 Molarity (M) pH 7.5

- Add  121.14 g Tris ultrapure 99.9% to a beaker
- Adjust volume to  800 mL with ddH₂O
- Adjust pH to  pH 7.5 with HCl
- Adjust volume to  1 L
- Sterilize by filtering and store at  Room temperature




1 L NaCl stock solution 5 Molarity (M)

- Add  292.2 g Sodium chloride to a beaker
- Adjust volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature

1 L EDTA stock solution 0.5 Molarity (M) pH 8

- Add  186.12 g EDTA disodium salt to a beaker
- Adjust volume to  1 L with ddH₂O
- Adjust pH to  pH 8 with sodium hydroxide
- Sterilize by filtering and store at  Room temperature

1 L SDS stock solution 10 Mass Percent

- Add  100 g SDS ultrapure to a beaker
- Adjust volume to  1 L with ddH₂O
- Sterilize by filtering and store at  Room temperature

1 L Proteinase K storage buffer (50 millimolar (mM) Tris ,

[M]3 millimolar (mM) CaCl₂ [M]50 % (v/v) glycerol) pH7.8

- Add 50 mL of [M]1 Molarity (M) Tris stock solution pH7.5
- Add 333 mg calcium chloride
- Add 500 mL of glycerol 87%
- Adjust volume to 900 mL with ddH₂O
- Adjust pH to pH7.8 with sodium hydroxide
- Adjust volume to 1 L with ddH₂O
- Sterilize by filtering and store at 8 Room temperature

Working solutions:

1 L TNES buffer ([M]50 millimolar (mM) Tris , [M]400 millimolar (mM) NaCl , [M]20 millimolar (mM) EDTA , [M]0.5 Mass / % volume SDS , pH7.5)

- Add 50 mL of [M]1 Molarity (M) Tris stock solution pH7.5
- Add 80 mL of [M]5 Molarity (M) NaCl stock solution
- Add 40 mL of [M]0.5 Molarity (M) EDTA stock solution pH8
- Add 50 mL of [M]10 Mass / % volume SDS stock solution
- Adjust volume to 800 mL with ddH₂O
- Adjust pH to pH7.5 with HCl
- Sterilize by filtering and store at 8 Room temperature

200 mL Proteinase K working solution ([M]10 mg/mL Proteinase K)











- Dissolve 2 g of Proteinase K in 200 mL Proteinase K storage buffer
- Store at 8 -20 °C

SAFETY WARNINGS

Reagents are potentially damaging to the environment. Dispose waste responsibly.

BEFORE STARTING

Make sure all buffers are prepared before starting.

- 1 For each sample prepare a screwcap tube pre-filled with a few 2 mm zirconia beads.
 - 2 Shake the sample well.
Transfer  **800 µL** of the small size fraction and  **200 µL** of the large size fraction to a 2 mL screwcap tube. It might be beneficial to use wide-bore tips or cut off the tip when using regular pipette tips.
 - 3  **11.000 x g, 00:03:00** 3m
 - 4 Remove as much ethanol as possible with a  **1000 µL** pipette.
 - 5 Add  **900 µL** of **TNES buffer** and  **100 µL** of **Proteinase K working solution**. Vortex shortly.
- Depending on the amount of samples this can be prepared as a mastermix. We usually prepare TNES + Proteinase K in batches for 24 samples. Proteinase K tends to self-digest if the time for samples preparation takes too long.
- 6 Bead-beat for  **00:02:00** at  **2400 rpm** 2m
 - 7 Incubate  **1400 rpm, 56°C, 00:20:00**
 - 8 Store at  **-20 °C** until DNA extraction.

