



Mar 30, 2021

Methanol-Chloroform-Water Precipitation

Bryon Drown¹, Kelleher KRG Research Group¹¹Northwestern University**1** Works for me dx.doi.org/10.17504/protocols.io.biknkcve

Kelleher Research Group

Tech. support email: kelleher-ofc@northwestern.edu

Kelleher Research Group

Northwestern University, National Resource for Translational...

SUBMIT TO PLOS ONE

ABSTRACT

Methanol-chloroform-water (MCW) precipitation is a rapid method for removing interfering species from protein samples prior to LC-MS analysis. MCW precipitation can effectively remove both ionic and non-ionic detergents and salts. Precipitating less than 1 µg protein is challenging yet possible with careful pipetting.

DOI

dx.doi.org/10.17504/protocols.io.biknkcve

PROTOCOL CITATION

Bryon Drown, Kelleher KRG Research Group 2021. Methanol-Chloroform-Water Precipitation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.biknkcve>

KEYWORDS

protein, mass spectrometry, top-down proteomics

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 15, 2020

LAST MODIFIED

Mar 30, 2021

OWNERSHIP HISTORY

Jul 15, 2020



Bryon Drown

Northwestern University

Mar 30, 2021



Kelleher Research Group

Northwestern University

PROTOCOL INTEGER ID

39278

MATERIALS TEXT

MATERIALS

☒ Protein LoBind 1.5mL microcentrifuge

tubes **Eppendorf Catalog #0030108116**

☒ Water Optima™ LC/MS Grade Fisher Scientific **Contributed by**

users Catalog #10728098

☒ Pierce™ Formic Acid **Thermo**

Fisher Catalog #TS-28905

☒ Methanol Optima™ LC/MS Grade Fisher Chemical **Fisher**

Scientific Catalog #A456-4

☒ Acetonitrile Optima™ LC/MS Grade Fisher Chemical **Fisher**

Scientific Catalog #A955-500

☒ Chloroform (HPLC grade) **Sigma Aldrich**

Resuspension buffer: 94.9% water, 4.9% acetonitrile, 0.2% formic acid

SAFETY WARNINGS

Halogenated solvent waste should be disposed properly as it poses environmental hazards.

- 1 Dilute protein sample to 150 µL with Optima-grade water.
- 2 Precipitate protein at interface
 - 2.1 Add 600 µL Optima-grade methanol and mix well.
 - 2.2 Add 150 µL HPLC-grade chloroform and mix well.
 - 2.3 Centrifuge 🌀 **20000 x g, 4°C, 00:10:00** and discard top layer while being careful to not disturb interface between layers. Use a gel-loading tip to remove last bit.
- 3 Pellet protein at interface
 - 3.1 Add 600 µL and homogenize the solvent, being careful not to break up the precipitated protein.
 - 3.2 Centrifuge 🌀 **20000 x g, 4°C** and discard supernatant, being careful not to disturb pellet.

4 Wash pellet with methanol. (This step can be repeated if a high level of contaminating species is anticipated)

4.1 Add 600 µL methanol and mix by inverting tube.

4.2 Centrifuge  **20000 x g, 4°C** and discard supernatant, being careful not to disturb pellet.

5 Redissolve protein in LC buffer

5.1 Evaporate excess methanol by laying tube on its side for 5 min at room temperature.

5.2 Redissolve precipitated protein in 30 µL Buffer A. Incubate 10 min on ice. If protein is not completely dissolved, sample can be briefly sonicated with a bath sonicator.

5.3 Centrifuge  **20000 x g, 4°C** and transfer supernatant to LC sample vial.