



Jun 26, 2024 Version 1

Immunofluorescence staining on larval and adult *Drosophila* gonads V.1

DOI

dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1

Samantha Goetting¹

¹Johns Hopkins University



Samantha Goetting

Johns Hopkins University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1

Protocol Citation: Samantha Goetting 2024. Immunofluorescence staining on larval and adult *Drosophila* gonads. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvjnd7ngk5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 26, 2024

Last Modified: June 26, 2024

Protocol Integer ID: 102463

Abstract

Immunofluorescence staining protocol for *Drosophila* gonads



Materials

1x PBS

0.3% PBTx (0.3% Triton-X in PBS)

1% PBTx (1% Triton-X in PBS)

Paraformaldehyde or formaldehyde

Normal serum (usually NGS)

Primary antibodies

Secondary antibodies

DAPI



Before start

All steps are done with gentle rotation.



Day 1

2h 40m

- 1 Dissect tissue in 1x PBS. Transfer to a 1.5 mL tube containing **1x PBS**. If it is a quick dissection (<20 mins), no ice needed. If more time is needed, keep samples on ice.
- 2 Remove PBS and add fixative. Fix in 4% paraformaldehyde in **0.3% PBTx** for 20 min RT with gentle rotation.
500 uL fixative = 125 uL of 16% paraformaldehyde + 375 uL 0.3% PBTx 20m 
- 3 Aspirate the fixative and wash twice for 10 min in **1% PBTx**. Not getting rid of fix will affect your immunostaining. 20m
- 4 Aspirate the supernatant and block/permeabilize for at least 2 hours in **1% PBTx** + 5% normal serum, or overnight at 4°C.
1 mL block solution = 50 uL NGS + 950 uL 1% PBTx 2h
- 5 Primary antibodies are diluted in **0.3% PBTx** + 5% normal serum and incubate for 1 hour at RT or overnight at 4°C. Overnight will give better staining.
Primary antibody mix in 0.3% PBTx + 15 uL NGS (300 uL total) 12h 

Day 2

5h 20m

- 6 Remove the primary antibody mix and wash in **0.3% PBTx** three times for 20 min at RT. 1h
- 7 Wash in **0.3% PBTx** + 5% normal serum twice for 30 min at RT.
1 mL wash solution = 50 uL NGS + 950 uL 0.3% PBTx 1h
- 8 Secondary antibodies are diluted in **0.3% PBTx** + 5% normal serum and incubated for ~2 hours at RT or overnight at 4°C. Keep tubes covered from light. 2h
- 9 Aspirate the supernatant and add 500 uL DAPI to each tube. Incubate for 10 min at RT. Keep tubes covered from light. 10m
- 10 Aspirate the supernatant and wash in **0.3% PBTx** three times for 20 min at RT. 1h
- 11 Aspirate the supernatant and wash in **PBS** for 10 min at RT. 10m



12 Store in **PBS** at 4°C or proceed with mounting. Keep tubes covered from light.

Protocol references

Slaidina et al. (2020)