



Sep 11, 2022

GUS staining

lilyli¹¹CAU

1 Works for me

 Share

This protocol is published without a DOI.

 lilyli

ABSTRACT

GUS staining

PROTOCOL CITATION

lilyli 2022. GUS staining. **protocols.io**
<https://protocols.io/view/gus-staining-cgemttc6>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 11, 2022

LAST MODIFIED


Sep 11, 2022

PROTOCOL INTEGER ID

69805

1 Fix with 90%  **90 % volume** acetone for  **00:20:00**

20m

2 Add  **1 mL** Gus washing solution to wash off the acetone and wash it twice

3 Add Gus staining solution and vacuum on ice for 15-20 min

- 4 37 °C overnight
- 5 Remove the dyeing solution, add 1ml of 70% ethanol
- 6 Change the ethanol several times
- 7