

Aug 29, 2024

Seizures – like behavior assay for D. melanogaster

DOI

dx.doi.org/10.17504/protocols.io.81wgbzqd1gpk/v1

Natalie Kaempf¹, Patrik Verstreken¹

¹VIB-KU Leuven Center for Brain & Disease Research, 3000 Leuven, Belgium.

ASAP Collaborative Rese...



Natalie Kaempf

VIB KU Leuven

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.81wgbzqd1gpk/v1**

Protocol Citation: Natalie Kaempf, Patrik Verstreken 2024. Seizures – like behavior assay for D. melanogaster. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgbzqd1gpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 29, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 104257

Keywords: ASAPCRN, seizure-like behavior



Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: ASAP-000430

EMBO long-term postdoctoral

fellowship

Grant ID: ALTF_299-2019

Research project, FWO

Vlaanderen

Grant ID: G0A5219N

Research project, FWO

Vlaanderen

Grant ID: G0B8119N

Methusalem project

Grant ID: METH/21/05

(3M210778)

Research project, KU Leuven

Parkinson Fonds

Grant ID: EQZ-PARFON-O2010

Opening the Future grant,

Leuven Universiteitsfonds

(LUF)

Grant ID: EQZ-OPTFUP-O2010

Research project, FWO

Vlaanderen

Grant ID: G031324N


Abstract

This protocol allows to assess the susceptibility of flies to seizure-like behavior based on mechanically hyper-stimulating sensory inputs and assessing seizure-like activity.



Seizures – like behavior assay

30m

- 1 prepare groups of 5 ± 1 flies and transfer them to transparent empty vials
- 2 let them rest for at least  00:30:00 after CO2 anesthesia 30m
- 3 mechanically stimulate the flies by vortexing the vial for 10 s at maximum intensity
- 4 move the vial immediately on an even surface before a white background and start recording the behavior
- 5 Numbers of resistant flies and flies with seizure-like behavior are counted within 10 s post stimulation

Protocol references

Kuebler, D., and M. A. Tanouye. 2000. "Modifications of Seizure Susceptibility in *Drosophila*." *Journal of Neurophysiology* 83 (2): 998–1009.
<https://doi.org/10.1152/jn.2000.83.2.998>.