

Version 2

Jul 28, 2021

Hyperspectral core-logger image acquisition V.2

Kevin Jacq¹, Jeanne Auboiron¹, Kévin Humbert¹, Ruth Martinez-Lamas¹, Antonin Van Exem¹,
Maxime Debret¹

¹University of Normandy Rouen

1 Works for me

Share

This protocol is published without a DOI.

Kevin Jacq

Université Rouen Normandie, Centre national de la recherche ...

ABSTRACT

The following protocol describes how to acquire hyperspectral images with Specim's Single Core Scanner (SCS) with two hyperspectral cameras (VNIR, SWIR). It was developed thanks to the expertise of many people within the M2C laboratory, as well as the Specim manual and a Butz publication.

This protocol also highlights important properties (exposure time and pixel overlap) on which it is necessary to take time to obtain the most informative hyperspectral image with an optimal signal-to-noise ratio.

Butz, C., Grosjean, M., Fischer, D., Wunderle, S., Tylmann, W., Rein, B., 2015. Hyperspectral imaging spectroscopy: a promising method for the biogeochemical analysis of lake sediments. J. Appl. Remote Sens. 9, 1–20.
<https://doi.org/10.1117/1.JRS.9.096031>

PROTOCOL CITATION

Kevin Jacq, Jeanne Auboiron, Kévin Humbert, Ruth Martinez-Lamas, Antonin Van Exem, Maxime Debret
2021. Hyperspectral core-logger image acquisition. **protocols.io**
<https://protocols.io/view/hyperspectral-core-logger-image-acquisition-bwxhpfj6>
Version created by Kevin Jacq

KEYWORDS

Hyperspectral imaging, Specim, Cores, Sediment core

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 27, 2021

LAST MODIFIED

Jul 28, 2021

PROTOCOL INTEGER ID


51913

SAFETY WARNINGS

- The white reference panel is fragile. Do not touch its surface.
- Do not touch the lens.


General description of the acquisition bench

1




General description of the hyperspectral core-logger (Specim)

2



Description of camera support: (a) front, (b) back

3




Description of the calibration table

Core logger set-up


- 4 Check that the **emergency stop buttons** placed along the core logger are unlocked. To do this, press and turn at the same time to unlock the buttons.

- 5 Pull the "**main power**" locker on the block at the base of the computer (1).



Power supply units constituting the hyperspectral core logger

- 6 Turn on the computer: **O/I "power"** button at the top left of the computer (2).
- 7 Check which camera is on the acquisition support. Depending on the problem related to the sample, the choice of the camera is important. Switch on the corresponding camera (3):
- For the **VNIR camera**: Turn on the small grey "**VNIR power supply**" unit located on the top of the computer.
 - For the **SWIR camera**: Pull the "**main power**" and "**HSI POWER**" lockers on the central unit of the computer.
- 8 Turn on the **illumination unit** at the base of the camera support, check that the position of the lamp is on "**VNIR**" or "**SWIR**" depending on the camera used.



Camera support with the illumination unit

- 9 Turn on the **air conditioning** to maintain a constant temperature in the room, as the temperature affects the signal.
- 10 Close the **blinds** on the room windows to keep the room dark. The **lights** should be turned off before acquiring the sample.
- 11 Take the **white reference** and place it (the smoothest part towards the camera) on its location on the calibration table.



The white reference is very delicate. Handle it by touching only the slices and dust its surface once placed in the calibration table with the **brush**.

- 12 Remove the **shutter** from the camera.



Hold the lens with one hand and point the shutter with the other hand do not touch the lens !


Adjustment of the height of the calibration table and the rulers

- 13 Place a **core** on the **sample tray**. Place the TOP of the core towards the camera, butting against the calibration plate.

When acquiring multiple cores, it is important to keep the same acquisition parameters and heights (samples, light and camera) to be able to compare data with each other.

Tip: calibrate the calibration table and ruler heights on the sample that has the greatest height, so it will just have to compensate the difference for the others.

- 14 Adjust the height of the **calibration table** with the worm screw located at the back under the calibration table. The height must correspond to the height of the sample (not necessarily the height of the container).





Adjustment of the calibration table height

- 15 Adjust the **ruler** against the core and to the height of the calibration table by turning the adjustment screw. The purpose of this point is to position the ruler as horizontally as possible using a spirit level which will allow later to adjust the horizontality of the sample.

Illumination unit adjustment


16

To check the correct height, place a white paper under the illumination unit at the calibration table height and check that the two light beams overlap.



Alignment of the light beams as seen through the aperture provided for the camera. (a) bad, (b) good alignment

if the two beams do not overlap, adjust the height of the **illumination unit** with the orange handle.



Illumination unit adjustment with the orange handle



The illumination unit must be held with one hand when handling the orange handle. The illumination unit may be hot!

Connecting the camera to the computer

- 17 Open the **LUMO** software to control the core logger.

- 18 In the "**Setup**" tab at the top right, select **Sensor 1** and **connect** the camera "**VNIR PFD with NI**" or "**SWIR3 with NI**", then **connect** "**Motor 1**" as well.

In the tab " **Sensor 1** " in " **Capture folder** ", you can choose the folder where the different acquisition files will be saved. Check that the disk where the data will be saved has enough space available.

Adjusting the camera height

19

The camera height adjustment has an influence on the resolution, but also on the depth of field (range of "sharpness" of the image). Therefore, depending on the sample and the interests of the study, a compromise will have to be made when adjusting the camera height.

We usually work at a height of 18 cm or 20 cm for VNIR and SWIR cameras.
The height is measured between the underside of the camera and the calibration table.

A	B	C	D	E
	VNIR	VNIR	SWIR	SWIR
Focal length	14 cm	35 cm	14 cm	35 cm
Pixel size	50 μm	120 μm	150 μm	370 μm
Volume of data	8 Go	3 Go	2,5 Go	0,8 Go

Choice of camera height and influence on pixel size, analysis time and data volume for a 1.5m sample length (spectral binning of 8 and 1 for the VNIR and SWIR cameras)

To help you choose this height, open the spreadsheet **Parameters.xlsx**. Only the grey boxes need to be filled in, the others are calculated automatically. Define the **working distance** (1) value for your problem, which will allow to calculate the spatial resolution and the working height to have the same angle of view with the other camera. The other calculated parameters (acquisition time and data size) will be refined with the definition of the other acquisition parameters.

 **Parameters.xlsx**

The **Parameters.xlsx** includes formulas for the calculation of several acquisition parameters thanks to a design of experiment. It is therefore specific to the SCS with these sensors and lens, and should be done for another device.

The experimental design used focuses on the effects of working distance, frame rate and scanning speed. A quadratic relationship was found between the FR/v ratio and the working distance. Thus, thanks to the equation that defines the overlap, it is possible to define the frame rate and the speed as a function of the working distance and the exposure time.

To build this experimental design, we vary the working distance, the frame rate or speed is fixed, and the other is determined by the optimizer based on the calculation of the radius of a round object (cf 22.2).

Step	VNIR OLE23			Step	SWIR OLE22,5		
1	Working distance	18	cm	1	Working distance	20	cm
1'	Working distance (dSWIR)	20,10	cm	1'	Working distance (dVNIR)	17,91	cm
	FOV	8,22	cm		FOV	8,18	cm
	Pixel size	62,63	µm		Pixel size	212,96	µm
	Ground Sample Distance	490,17	µm		Ground Sample Distance	1893,00	µm
	Depth of field	4,53	mm		Depth of field	5,56	mm
2	Spatial binning	1		2	Spatial binning	1	
2	Spectral binning	8		2	Spectral binning	1	
3	Exposition time	13	ms	3	Exposition time	3,5	ms
4a	Pixel overlapping	2,00	24,48 %	4a	Pixel overlapping	2,00	-11,54 %
	Frame rate	8,46	24,17 nb line / s		Frame rate	40,70	6,85 nb line / s
	Scan speed	0,58	1,67 mm/s		Scan speed	9,90	1,67 mm/s
4b	Sample length	1,5	m	4b	Sample length	1,5	m
4a'	Acquisition time	42,84	15,00 min	4a'	Acquisition time	2,52	15,00 min
	Data size	6,28	5,70 Go		Data size	1,59	1,39 Go
Fill in the grey cells.							
Apply in Lumo software the selected frame rate and scan speed.							
Reference the grey and green cells.							

Parameters spreadsheet with all the acquisition parameters used and calculated

- 19.1 The height of the camera is measured between the bottom of the lens and the calibration table. To place the sample tray on the calibration table, in **Lumo** in the "**Adjust**" tab, you must define its position in the "**Sensor 1**" - "**Focus grid**" box, it is at 150 (if this is not the case, change it, click on **Go**, wait for the sample tray to be placed, and validate with **set**).
- 19.2 Reduce or increase the height of the **camera** with the hand crank at the back above the camera support and lock the rack with the brake screw (small and black).




Camera height ajustement with the hand crank

- 19.3 Note the chosen height in the **spreadsheet** so that the calculated parameters are correct.

Adjusting the focus and determining the acquisition parameters

- 20 Adjusting the focus is an important step and one on which it is necessary to take time to achieve the optimal setting to obtain sharp images.


- 20.1 Click on the "**Adjust**" tab at the top right to see the interface:



LUMO user interface


20.2

To do this, place the sample tray on the focus grid. In the **"Adjust"** tab, you must define its position in the **"Sensor 1" - "Focus grid"** box, it is at 150 (if this is not the case, change it, click on **Go**, wait for the sample tray to be placed, and validate with **set**). And **set** the **exposure time** (top left) to **13 ms** for the **VNIR** camera and **3 ms** for the **SWIR** camera. At this point, the spatial profile show in the software interface:



Focus grid adjustment

20.3 Rotate the **lens** (below the camera) so that the ripples on the profile in **"spatial view"** look like **slots** as much as possible. To do this, you can use the dark grey grid in the graph to find your way around. It is also possible to zoom, with a right click, in autoscale uncheck autoscale x, then in option choose a range in the center allowing to see some variations.



Focusing the focus by rotating the lens

20.4 An optimization can be made in the "**Adjust**" tab, then "**Focusing**", press **Run**. This tool will help you find the optimal focus setting. Then, turn the **lens** to scan different focus settings (slots (optimal) or sinusoids (non-optimal) in "**spatial view**"). Finally, rotate the lens until you reach a value close to 100%. If the percentage remains at 0%, redo manually and keep this adjustment.

21 Data correction:
In order to view images with colors close to what is observed, it is necessary to calibrate them with **black** (automatic shutter closing) and **white** (with the white reference positioned on the calibration table). To do this, go to the "**Adjust**" tab, then to "**Image correction**":

- Dark subtraction": Select **Shutter** and press **Run Dark**
- White reference": Select **Scanner**, press **Run White**

This correction may not have the desired effect. There may also be a difference between viewing in "**Adjust**" and "**Capture**".

22 Acquisition parameters:
In the "**Adjust**" tab, in the top left square, you will find the main parameters to be adjusted to perform the acquisition.

- Shutter: Which must always be with an orange square (Opened)
- Trigger mode: Who should mark: "Internal".
- **Spectral binning**: Binning allows to cumulate intensities of neighbouring wavelengths. In the case of a binning of 2, we combine 2 neighboring wavelengths. This also allows to reduce the number of wavelengths and thus the size of the output data file, but this can hide important spectral information. (Usually set to 8 for the VNIR camera and 1 for the SWIR camera, VNIR binning of 1 should be used if organic pigment are present in the sample).
- **Spatial binning**: Here too, binning allows you to average the intensities of neighbouring pixels. In the case of a binning of 2, we combine 2 neighboring pixels. This also allows to reduce the number of pixels and thus reduce the resolution and size of the output data file. (Usually set to 1 for both cameras).

A	B	C	D	E
	VNIR	VNIR	SWIR	SWIR
Spectral Binning	8	1	2	1
Spatial Binning	1	1	1	1
Volume of data	7,5 Go	60 Go	1 Go	2 Go

Impact of the binnings for a 150cm sample

Fill in these chosen values in the acquisition software, validate with **Apply**.

Note the chosen binning in the **spreadsheet (2)** so that the calculated parameters are correct.


22.1 Spectral quality (exposure time):

Exposure time: Adjusts the recorded light intensity that affects the signal-to-noise ratio. If the time is too long, there is a risk of signal saturation, and if the time is too long, the signal may be masked by noise.

- VNIR: minimum: 0.1 ms, maximum: 100 ms
- SWIR: minimum: 0.1 ms, maximum: 20 ms

Define the exposure time with the white reference. To do this, go to its level in the tab "**Adjust**", "**Sensor 1**", "**White**" (at the bottom right), it is usually defined between 90 (start) and 100 (end). Click on **Go** for one of the two positions of the white, wait for the sample tray to move.

- Adjust the **Exposure time** (top left) to find out when the signal is saturated (plateau at maximum values) and choose a value corresponding to 75% of it. (Usually set to 13 ms for the VNIR camera and 3.5 ms for the SWIR camera).



Signal saturation due to excessive exposure time

Fill in this chosen value in the acquisition software, validate with **Apply**.

Note the chosen exposure time in the **spreadsheet (3)** so that the calculated parameters are correct.

22.2 Spatial quality (overlap):

Overlap is the amount of a pixel that is found in the next pixel. In the perfect case, it would be 0%. In the case of our images this is not the case, so we try to approach it with some of the acquisition parameters that will constrain the value of the overlap.

The Excel spreadsheet (**Parameters.xlsx**) has been created to estimate the **frame rate** and **scan speed** values from the previous acquisition parameters (**working distance**, **binning**, **exposure time**). The frame rate and scan speed can be estimated either from the overlap or from the analysis time.

- **Frame rate** (Hz): Represents the number of acquisitions per second. (maximum: 150 fps VNIR, 450 fps SWIR)
- **Scan speed** (m/s): Speed of the traveling bench (minimum: 0.50 mm/s)

The overlap is usually set to 2% (4) for both cameras if the frame rate and scan speed range allow it, otherwise the overlap value is set higher.

Fill them in in the acquisition software (top left for the frame rate, top right and manual for the scan speed), validate with **Apply**.

Optimization of the values or manual selection of the frame rate and scan speed:

1. Place the **sample tray** at the dark level with "**Sensor 1**" - "**Dark start**" then **Go**.

2. Place a **round coin** on the **indicator plastic**.
3. In "**Sensor 1 positions**", check that the "**Target start**" is at 180 and the "**Target stop**" at 220 (if not, change it, click **Go**, wait for the measuring bench to be placed, and validate with **set**).
4. In "**Adjust**", "**Scanning speed Calculation**", select **Aspect ratio**, then **Run**. The software will then offer you speed or frame rate values close to those previously estimated. Select the **speed** change. Then select Manual to see the new speed.

Sediment core preparation

23

It is very important to keep the core out of the freezer as little time as possible. It should be based on a maximum time of 1 hour and a minimum of 15 minutes. Beyond 1 hour, the sample may crack. Below 15 minutes, the sample may be too humid, which will induce important absorptions with the SWIR camera, and reflection for the VNIR camera. These times vary from one sample to another depending on the properties of the sedimentary matrix.

The surface of the core should be as smooth as possible, because the depth of field of the cameras is small, in the order of a few millimetres. And any variations in the surface have an impact on the signal. Biofilms may also be present on the surface of the sample, which will partially mask the sediment signal. Similarly, oxidation of the sample surface may mask underlying sedimentary structures. Therefore, it is important to take the time to "clean" (shave the surface with a **spatula** or **knife**).



In areas with shells or sand, "cleaning" can very easily leave marks on the surface. To avoid this, it is recommended to "clean" these shell areas with very little pressure and then smooth the surface with your finger by wearing a **latex glove** slightly moistened (the smoothest part on the sediment). Very light finger pressure should be applied to the sediment. If the sediment is very clayey and very wet, avoid using the finger with the glove but use only the spatula with very little pressure.



It is also important to remove elements that are not part of the sediment (mould, small pieces of plastic...).

- 24 Once cleaned, place the **core** on the **sample tray**. Place the TOP of the core towards the camera, butting against the calibration plate.

Using the ruler, adjust the horizontality and height of the **sample** by placing objects (pieces of cardboard, polystyrene) between the tray and the sample.

Acquisition of hyperspectral images

- 25 In "**Sensor 1 positions**", define "**Target start**" and the "**Target stop**" corresponding to positions on the sample tray that cover the entire length of the core to be passed. Also be sure to click "**Set**" each time you change these numbers. (Start at 180, end less than 2,000)
- 26 Turn off the **lights** in the room.
- 27 Check that the **illumination unit** is turned on according to the **camera** used ("VNIR" or "SWIR"). Check one last time that the **acquisition parameters** are correct and validate them with **Apply**.

Make a final correction to the image. To do this, go to the "**Adjust**" tab and then to "**Image Correction**" :

- 28
- Dark subtraction ": Select **Shutter** and press **Run Dark**
 - White reference": Select **Scanner**, press **Run White**
- 29 Write the name of the core on a piece of paper and place it on the **indicating plastic** between the core and the focus grid in order to have a mark of the name of the core on the acquired image.
- 30 In the "**Capture**" tab, give the name of the sample (dataset name) and put some comments if necessary. Launch the capture, by clicking on "**Record**".



It is not advisable to place heavy objects on the table so as not to alter the horizontality of the sample tray.

Data standardization

31

Normalize the image using Dark and White acquisitions to remove instrumental noise and obtain calibrated reflectance.

Open the software: **ENVI classic+IDL**.

- 32 Click on "**Specim**," "**Scan Normalization**". Return to the acquisition folder and select the file: **manifest**.
- 33 In the new windows, select the files that are requested, for example the "white" or "dark reference". These files are located in the **Capture folder**.
- 34 Once all the files are selected. A new window will appear with a blue bar: "Normalization".
- 35 As soon as normalization is complete, new windows will be opened that can be closed. A new file is also created during normalization (_refl.dat). This file will be used for all future processing.
- 36 Then close **ENVI classic+IDL**. Then repeat the previous steps for a new acquisition or continue this protocol.

For the acquisition of multiple cores, it is also possible to complete all acquisitions before doing all normalization.

When acquiring multiple cores, try to keep the same software settings and heights (samples, illumination and camera) to be able to compare data with each other. For this, it is advisable to always work at the same focal length, so the other acquisition parameters are the same.

37 In the **LUMO** acquisition software,


37.1 Go to "**Setup**" in the upper right-hand corner, on the left in the "**Motor 1**" tab, click on "**Disconnect**".

37.2 Then click on "**Disconnect**" in the tab "**Sensor 1**" (upper left-hand corner).

37.3 Close the **LUMO** software.

38 Turn off the **illumination unit** at the base of the **camera support**.

39 Switch off the **computer** in a conventional way, and with the **O/I "power"** button (top block) (1).



Switching off the hyperspectral core logger

40 Switching off the **camera** (2)

- For the VNIR: Turn off the small grey "**VNIR power supply**" box located on the top of the computer.
- For the SWIR: Pull the "**HSI POWER**" and "**main power**" lockers on the central block of the computer.

41 Pull the lockers "**Main power**" at the base (3).

42 Turn off the **air conditioning** in the room.

43 Store the **white reference**.

44 Clean the various components of the room (tray, benches, tools) and put them away.