



DEC 05, 2023

🌐 CMV Resistance testing (UL54 and UL97)

Fernando Lazaro¹

¹Servicio de Microbiología. Hospital Universitario La Paz



Fernando Lazaro

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://dx.doi.org/10.17504/protocols.io.bp2l6xy4zlqe/v1) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://dx.doi.org/10.17504/protocols.io.bp2l6xy4zlqe/v1), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

OPEN  ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.bp2l6xy4zlqe/v1

Protocol Citation: Fernando Lazaro 2023. CMV Resistance testing (UL54 and UL97). [protocols.io](https://dx.doi.org/10.17504/protocols.io.bp2l6xy4zlqe/v1)
<https://dx.doi.org/10.17504/protocols.io.bp2l6xy4zlqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 30, 2023

Last Modified: Dec 05, 2023

ABSTRACT

This protocol is a procedure for the study of antiviral resistance in Cytomegalovirus by NGS techniques.

The primers have been designed using <https://primalscheme.com/> with the intention of covering the most relevant regions of the UL54 and UL97 genes.

BEFORE START INSTRUCTIONS

Take into account that the quality of the results is greatly affected by the time from sample extraction to DNA amplification.

The protocol may fail if the protocol is performed from refrigerated samples or DNA.

Prepare Reagents

- 1
- Q5® High-Fidelity DNA Polymerase (New England)
 - Agarose gel 1%
 - Ethanol 70%
 - Mag-Bind® TotalPure NGS (omega)
 - Elution Buffer

Set up primer pools (UL97 and UL54)

2 UL97 Primers

UL97_4_LE FT	TGCGCGCGGAAAGT CAG	UL97_4_RI GHT	CGGCATAACAGATCT TGTGGC
UL97_5_LE FT	CTCTGCGAGCTCTC TATCTCCT	UL97_5_RI GHT	AGCAGACAGCAGCCC GT
UL97_6_LE FT	TGGCGAGCAACAGC AGC	UL97_6_RI GHT	GCGCGCATGATCTCG CT
UL97_7_LE FT	TGCCACTTTGACAT TACACCCA	UL97_7_RI GHT	TCCGACATGCAATAA CGCCG
UL97_8_LE FT	TTTCCGACCCATGC CGCT	UL97_8_RI GHT	ATGCTCGCCCAGGAG ACAG
UL97_9_LE FT	CATGGGTACGGAGG CGTTG	UL97_9_RI GHT	GGCCAACAGACGCTC CA



UL97 Primers

UL54 Primers

UL54_1_LE FT	TGCAAAAACCTTGTCC TTGCGC	UL54_1_RI GHT	ATTCTGTAACCAACCG GCGTG
UL54_2_LE FT	GTAGTTGCACACGGC CGAC	UL54_2_RI GHT	CGTCAATCTAACCTG CCGCA
UL54_3_LE FT	CGTAAAAGACCCGAT CCCCG	UL54_3_RI GHT	TCTCGCTGCTCTTTG AGGATC
UL54_4_LE FT	CTTCATCGAGTGAGA GGCGC	UL54_4_RI GHT	AGGCTTTGGTGGCGC GT
UL54_5_LE FT	GCTTGACGGGCTCCA CAAAA	UL54_5_RI GHT	GGCGCGGTTTCATCAA AGACA
UL54_6_LE FT	TCCCGCGTTCCCACT ACATA	UL54_6_RI GHT	CAACAAGTGGGTTTC GCAGC
UL54_7_LE FT	ATACGGCGCACAGG GTCTT	UL54_7_RI GHT	GTGTTTGAGCCCGAG GTGG
UL54_8_LE FT	AGTAGCAGAGGTTGT GAGCCA	UL54_8_RI GHT	GGTTCTGTGGCGGCT ATGTT
UL54_9_LE FT	AAACGCCGTCCTGAC TCGA	UL54_9_RI GHT (V2)	CTTGCAATCTGCGCC GTC
UL54_10_L EFT	GGATCTGCTGTCCGT CAAAGA	UL54_10_ RIGHT	ATATTGCGGGTTCGG TGTT
UL54_11_L EFT	TGTTGAGCTTATAGT TGGGCGA	UL54_11_ RIGHT	CGGCCTTTGTGACCG GTTAC
UL54_12_L EFT	CCTTATACAGGTACT CGAGGCG	UL54_12_ RIGHT	GTGCTACGAGACGGG AGGA
UL54_13_L EFT	AAGTGCAGCCCCGA CCAT	UL54_13_ RIGHT	GGATCACCACGTTCG GCTG
UL54_14_L EFT	CCTCGATATCACAAG TCGACGC	UL54_14_ RIGHT	GGCGAACTAGTGCCC GAAC

UL54_15_L EFT	CCGTACCCGTTAGATG GAGGT	UL54_15_ RIGHT	GGGACCTATTCGTTT TCACACCTA
UL54_16_L EFT	ACGATAGCGCGGCG ACA	UL54_16_ RIGHT	CGGCGTCAGCGTTTG CA

UL54 Primers

- 2.1** Pool 1 UL97 (odd primers LEFT and RIGHT)
Pool 2 UL97 (even primers LEFT and RIGHT)
Use  10 µL for each primer at 10 micromolar (µM)
- 2.2** Pool 1 UL54 (odd primers LEFT and RIGHT)
Pool 2 UL54 (even primers LEFT and RIGHT)
Use  10 µL for each primer 10 micromolar (µM) except:
- UL54_14_LEFT/RIGHT 14: 2,5 µl
 - UL54_5_LEFT/RIGHT: 5 µl
 - UL54_13_LEFT/RIGHT: 5 µl
 - UL54_15_LEFT/RIGHT: 20 µl

DNA extraction

- 3** Perform DNA extraction with your method of choice. Preferably from a **plasma sample** collected on the **same day**. Perform DNA extraction with your method of choice. Preferably from a plasma sample collected on the same day.

The quality of the results is greatly affected by the time from sample extraction to DNA amplification. The protocol may fail if the protocol is performed from frozen samples or DNA.

DNA amplification

- 4** For each pool and sample, mix the following reagents (two reactions per sample):

Reagent	Volume / sample (µl)
Primer Pool	2
Q5 [®] Polymerase	0,25

Q5 Buffer	5
H2O	5
dNTP	0.5
Sample DNA	12,5

PCR Mix for Pool 1 and 2

Set up the PCR with the following program

A	B	C
Cycles	Temperature	Time
1	98°C	30"
5	98°C	10"
	65°C	3'
30	98°C	10"
	65°C	30"
	72°C	2'
1	72°C	2'
1	4°C	∞

PCR program

Confirm amplification of approximately 300bp fragments by **1% agarose gel**.

Product cleaning with Mag-Bind TotalPure NGS (omega)

- 5 Add 40 uL of beads to 20 uL of amplicon. Mix.

- 6 Incubate the mixture for 5 minutes at room temperature.
- 7 Place the tube in the magnet until the solution becomes clear.
- 8 Gently remove the supernatant by pipette.
- 9 Add 180 uL of 70% ethanol. Mix without breaking the pellet.
- 10 Gently discard the ethanol by pipette.
- 11 Incubate for 2 minutes at room temperature.
- 12 Remove any remaining ethanol. Note: The pellet must not be allowed to dry excessively. If it does occur, the pellet will appear black and cracked.
- 13 Remove the tube from the magnet.
- 14 Add 25 uL of EB. Mix. Note: do not break the pellet, just peel it away from the wall of the tube.

- 15 Incubate for 2 minutes at room temperature.
- 16 Place the tube back in the magnet until the solution clears.
- 17 Transfer 20 uL of the supernatant to a new tube.

NGS sequencing

- 18 Sequence the amplicons using the sequencer of choice according to the manufacturer's instructions.

Antiviral resistance

- 19 Generate consensus sequence from the amplicons using the bioinformatics procedure of choice. Take into account that in cases of previous exposure to antivirals or prolonged treatment, minority variants may appear.
The website <http://cmv-resistance.ucl.ac.uk/herpesdrg/> is helpful for the study of resistance mutations.