



Oct 18, 2021

RCA of Gotcha

Chia-Hsien Shih¹¹Chung Shan Medical University

1


dx.doi.org/10.17504/protocols.io.by53py8n

CSMU_Taiwan



Chia-Hsien Shih
Chung Shan Medical University

This protocol is to test that GotCha is functional and can works as designed.

DOI

dx.doi.org/10.17504/protocols.io.by53py8n

Chia-Hsien Shih 2021. RCA of Gotcha. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.by53py8n>



_____ protocol ,

Oct 18, 2021

Oct 18, 2021









54171

Preparation

- 1 Add **5 μ L** of GotCha(functional beads) into eppendorf
- 2 Centrifuge for **15000 rpm, 00:05:00** and remove supernatant. Make sure that eppendorf should put on DynaMag when removing supernatant. ^{5m}


Protocol

- 3 Add **3 μ L** of 10X phi29 polymerase reaction buffer into eppendorf with GotCha

- 4 Add  **20.4 μ L** RNase-free water
- 5 Add  **3 μ L** of 100nM miRNA
- 6 Add  **3 μ L** of 2mM dNTPs
- 7 Add  **0.6 μ L** of 10U/ μ l phi29 polymerase
- 8 Pipetting to mix well
- 9 Incubate for  **02:00:00** at  **Room temperature**
- 10 Add  **1.5 μ L** of 20X evagreen dye
- 11 Add  **1.5 μ L** of 0.5M EDTA buffer to suspend the reaction

2h

Measuring

- 12 Load  **20 μ L** of reaction solution into 384-well plate
- 13 Measure the fluorescence excitation and emission intensity

