

Jul 15, 2024

# Senescence induction by DNA-damage in PCLS

DOI

[dx.doi.org/10.17504/protocols.io.e6nvwdwpwlmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwdwpwlmk/v1)



Delphine Beaulieu<sup>1</sup>, Nayra Cardenes<sup>1</sup>, Melanie Königshoff<sup>1</sup>

<sup>1</sup>University of Pittsburgh

TriState SenNet

Cellular Senescence Net...

1 more workspace



Nayra Cardenes

Division of Pulmonary, Allergy, Critical Care and Sleep Med...

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.e6nvwdwpwlmk/v1](https://dx.doi.org/10.17504/protocols.io.e6nvwdwpwlmk/v1)

**Protocol Citation:** Delphine Beaulieu, Nayra Cardenes, Melanie Königshoff 2024. Senescence induction by DNA-damage in PCLS. protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvwdwpwlmk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 22, 2024

**Last Modified:** July 15, 2024

**Protocol Integer ID:** 95638

**Keywords:** PCLS, Senescence, SenNet, Bleomycin, Doxorubicin, TriState SenNet, Lung

**Funders Acknowledgement:**  
**TriState SenNET (Lung and Heart) Tissue Map and Atlas consortium - NIA**  
**Grant ID: U54AG075931**

## Abstract

Precision-Cut Lung Slices (PCLS) are uniform, in size and thickness, tissue slices generated from human lungs that can be used for disease modeling, drug discovery and preclinical validation, while retaining the lung's cellular complexity and tissue architecture.

This protocol describes two different methods for the induction of cellular senescence by DNA damage in PCLS, and collection of samples for different analysis and measurements.

## Image Attribution

CLAR (Center for Lung Aging and Regeneration) Logo designed by Nayra Cardenes, PhD.

## Guidelines

All steps involving human tissue must be conducted under BSL2 safety guidelines.

## Protocol materials

	Penicillin-Streptomycin <b>Merck MilliporeSigma (Sigma-Aldrich) Catalog #P00781</b>	Step 1
	Amphotericin B solution <b>Merck MilliporeSigma (Sigma-Aldrich) Catalog #A2942</b>	Step 1
	WST-8 Assay Kit (Cell Proliferation) <b>Abcam Catalog #ab65475</b>	Step 5.2
	Bleomycin <b>Fresenius Kabi Catalog #10361</b>	Step 4.2
	Dimethyl sulfoxide (DMSO) <b>Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2438-50ML</b>	Step 4.2
	Doxorubicin hydrochloride <b>Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1515-10MG</b>	Step 4.2
	DMEM/F-12, powder, HEPES <b>Thermo Fisher Scientific Catalog #12400024</b>	Step 1
	16% Formaldehyde (w/v), Methanol free <b>Life Technologies Catalog #28908</b>	Step 5.1
	Senescence $\beta$ -Galactosidase Staining Kit <b>Cell Signaling Technology Catalog #9860</b>	Step 5.3



## Safety warnings




! The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## Before start




To start this procedure, you must have PCLS (Precision-Cut Lung Slices) ready (check protocol DOI [dx.doi.org/10.17504/protocols.io.36wgq3xr5lk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgq3xr5lk5/v1)).



## Day -1: PCLS acclimation to growth conditions

- 1 Add 1 mL of media per well in 24-well plates DMEM F-12 media (  
 DMEM/F-12, powder, HEPES **Thermo Fisher Scientific Catalog #12400024** ) + 1%  
FBS + 1% Penicillin-Streptomycin  
 Penicillin-Streptomycin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P00781** +  
0.3 µg/mL Amphotericin B  
 Amphotericin B solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A2942** .
- 2 Place one PCLS (1 cm diameter / 300 µm thick) per well with a brush.
- 3 Incubate o/n (overnight) at 37 °C, 5% CO<sub>2</sub>.

## Treatment at Days 0, 2 and 4

- 4 Treatment: change media with 1ml/PCLS, diluted in DMEM media + 0.1% FBS + 1% Pen/Strep + 1% Amphotericin B
- 4.1 Prior to media exchange, collect supernatant (s/n) from each condition in 1.5 mL tubes and, store at -80°C for future analysis.
- 4.2 Change media for each condition:  
  
Control - Change media (above).  
Bleomycin - 15 µg/ml Bleomycin in media above.  
 Bleomycin **Fresenius Kabi Catalog #10361**  
DMSO - 1:100,000 DMSO in media above.  
 Dimethyl sulfoxide (DMSO) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D2438-50ML**  
Doxorubicin - 0.1 µM Doxorubicin in media above.  
 Doxorubicin hydrochloride **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1515-10MG**


## PCLS Harvesting for measurements



5 Collect PCLS at Day 0, 2, 4 and Day 6 for different baseline and final timepoint measurements.

5.1 Days 0, 2, 4 6:

Fix PCLS in 4% formaldehyde/PBS for 30 minutes (2x PCLS per condition).


 16% Formaldehyde (w/v), Methanol free **Life Technologies Catalog #28908**

Then wash twice with PBS and store at 4°C in PBS.

Use 1-2 slices for paraffin embedding, and posterior histological analysis.

5.2 Days 0 and 6:

Cell proliferation measurement: Use 2-3x PCLS to make 4x4 mm punches to perform WST-8 cell proliferation assay in a 96-well plate, following manufacturer's instructions.

 WST-8 Assay Kit (Cell Proliferation) **Abcam Catalog #ab65475**

5.3 Days 0 and 6:

Fix 2x PCLS with 1x fixative solution from  $\beta$ -galactosidase kit and stain for  $\beta$ -galactosidase as per manufacturer instructions.

 Senescence  $\beta$ -Galactosidase Staining Kit **Cell Signaling Technology Catalog #9860**

5.4 Day 0 and 6:

Snap freeze in liquid nitrogen: 4x PCLS per cryotube.

Store at -80°C for future analysis.